

OPTIM 33TB, exceptional cleaning
& disinfection capabilities.

A breath of fresh air.



- One step disinfectant cleaner – fast & convenient
- AHP – patented technology designed to deliver better cleaning and faster disinfection results (active ingredient 0.5% hydrogen peroxide)
- Hydrogen peroxide does not leave behind active residue which may lead to the development of resistant strains of microorganisms



Your Infection Control Specialist™

SciCan
A Sanavis Group Company



FAST CONTACT TIME:

- 1 minute bactericidal & virucidal claim
- 30 second sanitizer
- 5 minute TB claim



- No strong chemical odor
- Active ingredient breaks down into oxygen and water
- Non-flammable
- EPA Cat. 4 for skin & eye irritation
- 0 HMIS health hazard rating*
- EPA Cat. 4: lowest toxicity level on EPA toxicity category

Active Ingredient	EPA Toxicity Category	Signal Word Required
AHP	Cat.4	NONE
Quats+Low % Alcohol	Cat.3	CAUTION
Quats+High % Alcohol	Cat.2	WARNING
Phenolics	Cat.1	DANGER

* Refer to HMIS rating on MSDS & compare with our competitors' wipes

SURFACE DISINFECTANT CONTACT TIME**

Product Name	Active Ingredient	Mycobacteria	Non-enveloped virus			Enveloped virus			Bacteria				
		TB	Poliovirus	Norovirus	Rotavirus	HIV	HBV (Hep.B)	HCV (Hep.C)	MRSA	VRE	Influenza A2	E.coli	Pseudomonas
OPTIM 33 TB	AHP	5	1	1	1	1	1	1	1	1	1	1	1
Birex SE	Phenolics	10	10			10	10		10		10	10	10
Prospray	Phenolics	10	10			1	10	10	10		10	10	10
Caviwipes	Quats-alcohol	3				2	2	2	2	2			3
Opticide 3	Quats-alcohol	3	3			3	3	3	3	3	3	3	3
Sanicloth Plus	Quats-alcohol	3				2	2	2	2	2		2	3
Lysol IC	Quats-alcohol	10	10			10							10

** Contact time in minutes. Contact time was obtained from publicly available material & is valid when table was constructed. Information will be updated if there are changes in claims. If blank, no published information was available.

ORDERING INFORMATION FOR OPTIM 33TB

Ready-To-Use Pre-saturated Wipes:

OPT33-W12 RTU Wipes 6" x 7" (12 x 160 sheets/case)

Dry Wipes & Dry Wipes Refill:

OPT-EMPCAN STAT-Wipes Canister & Roll of Dry Wipes (12 x 160 sheets /case)

SCI33R18 Dry Wipes Refill, 6" x 7" (18 x 160 sheets/case)

SCI33R18XL Large Dry Wipes Refill 10" x 10" (18 x 60 sheets/case)

Liquid:

OPT33-12x32 1 case of 12 x 32 oz.

OPT33-4x1 1 case of 4 x 1 gallon

FSA01 Sprayer for 32 oz. bottle

SciCan Inc., 701 Technology Drive, Canonsburg, PA 15317 • Phone (724) 820-1600; Fax (724) 820-1479; Toll Free 1-800-572-1211
www.scican.com