

DACRIOS

PUNCTAL PLUGS

OPTIONS TO MEET YOUR PATIENT NEEDS AND EXPECTATIONS


Permanent plugs are manufactured from soft silicone in a single mold design, our plugs feature a low-profile cap design that virtually eliminates foreign-body sensation:

- Flexibility of cap design facilitates easy removal with jeweler's forceps, if necessary
- Individually packaged and preloaded on a sterile, disposable inserter / dilator



Absorbable plugs are available in short term dissolving collagen from 0.2mm to 0.4mm diameter and extended wear PCL materials from 0.2mm to 0.5mm diameter to assist with short and midterm needs like allergies and post-op recovery.

- Indicated for use in the treatment of Dry Eye Syndrome and the dry eye components of various ocular surface diseases
- May be used after ocular surgery to prevent dry eye related complications
- Ideal for Post-Op dry eye management of LASIK patients
- Treatment option for seasonal allergy patients experiencing dry eye symptoms
- May be beneficial for contact lens patients experiencing dry eye related symptoms
- Manufactured from E-Caprolactone-L-Lactide copolymer (PCL)



Essentially absorbed within 2-6 months (60-180 days)

VISION-SURGICAL.COM

INFO@VISION-SURGICAL.COM

DACRIOS SILICONE PUNCTAL PLUGS

Inserted into the punctal opening to block tear drainage through the canaliculus, silicone plugs provide long-term relief of dry eye symptoms, and enable the eye to retain more of its own natural tears. 2 per box.

- Low-profile cap design virtually eliminates foreign body sensation
- Flexibility of cap design facilitates easy removal with jeweler’s forceps, if necessary
- Individually packaged and preloaded on a sterile, disposable inserter / dilator
- Manufactured from soft silicone utilizing a single piece mold design for greater patient comfort

ITEM NUMBER	DESCRIPTION
DPPS04	.4MM PLUG/INSERTER
DPPS05	.5MM PLUG/INSERTER
DPPS06	.6MM PLUG/INSERTER
DPPS07	.7MM PLUG/INSERTER
DPPS08	.8MM PLUG/INSERTER



DACRIOS EXTENDED WEAR SYNTHETIC ABSORBABLE PUNCTAL PLUGS

Inserted into the punctal opening to block tear drainage through the canaliculus, synthetic absorbable plugs provide extended temporary relief of dry eye symptoms. 2 plugs per package /10 packages per box.

- Indicated for use in the treatment of Dry Eye Syndrome and the dry eye components of various ocular surface diseases
- May be used after ocular surgery to prevent dry eye related complications
- Ideal for post-op dry eye management of LASIK patients
- Treatment option for seasonal allergy patients experiencing dry eye symptoms
- May be beneficial for contact lens patients experiencing dry eye related symptoms
- Manufactured from E-Caprolactone-L-Lactide copolymer (PCL)
- Essentially absorbed within 2-6 months (60-180 days)



ITEM NUMBER	DESCRIPTION
DPXT02	0.2MM DIAMETER/2.0MM LENGTH
DPXT03	0.3MM DIAMETER/2.0MM LENGTH
DPXT04	0.4MM DIAMETER/2.0MM LENGTH
DPXT05	0.5MM DIAMETER/2.0MM LENGTH

DACRIOS COLLAGEN ABSORBABLE PUNCTAL PLUGS*

Inserted into the punctal opening to block tear drainage through the canaliculus, collagen plugs provide temporary relief of dry eye symptoms. 6 plugs per package /10 packages per box.

- Manufactured from precisely measured premium collagen
- Used as a diagnostic assessment for permanent punctal occlusion
- Each individually packaged card contains 6 plugs
- Essentially absorbed within 7 to 10 days

ITEM NUMBER	DESCRIPTION
DCP02	0.2MM DIAMETER/2.0MM LENGTH
DCP03	0.3MM DIAMETER/2.0MM LENGTH
DCP04	0.4MM DIAMETER/2.0MM LENGTH



*DACRIOS COLLAGEN ABSORBABLE PUNCTAL PLUGS ARE NOT CE MARKED