

# STERILE MARKERS

**Vision Surgical**

**[www.vision-surgical.com](http://www.vision-surgical.com)**

P:1.888.518.4746

F:1.214.722.2020

**[info@vision-surgical.com](mailto:info@vision-surgical.com)**



# ABOUT OUR MARKERS



## TRADITIONAL INK

Traditional ink is Gentian Violet ink and has been used in the industry to mark patients for over 40 years. It can withstand betadine and iodine scrubs. AllSkin® (pg 9) features a brilliant blue option available for visibility on all skin tones. Ophthalmic surgeons have used Traditional Ink intraocularly for over 40 years.

*Ideal for procedures above the neck or where betadine and iodine are used.*

*Featured on pages 8-12*



## GV FREE INK

Viscot Medical's GV Free ink is a gentian violet free alternative to Viscot's Traditional ink. Eliminates concerns over Gentian Violet while maintaining the classic purple color. To be used on skin only.

*Ideal for procedures above the neck or where betadine and iodine are used.*

*Featured on page 13*



## XL PREP RESISTANT INK

Viscot XL Prep Resistant Ink is visible up to 10x longer vs. traditional gentian violet ink. Able to withstand a 60 second prep of Chloraprep, XL Ink remains visible against all common alcohol based prep solutions. Available in gentian violet free red and green ink for enhanced visibility on all skin tones. To be used on skin only.

*Ideal for any procedure where aggressive preps are used.*

*Featured on pages 15-17*



## EZ REMOVABLE INK™

Viscot Medical's EZ Removeable Ink is designed to come off easily with an alcohol swab or even warm soap and water. Available in red, green, and white.

*Ideal for non-surgical aesthetic procedures.*

*Featured on page 19*



## WATERPROOF UTILITY INK

Black Utility ink is designed for writing on medication labels. Smudge-proof, waterproof, and fluid resistant.

*Ideal for writing on labels, glass, vinyl and other surfaces.*

*Featured on page 25-26*

# ABOUT OUR MARKERS



## ULTRA FINE TIPS

Soft rounded Ultra Fine tips create fine lines while eliminating concerns over causing trauma to delicate tissues and reducing frustration from sharp marker points snagging making it difficult to draw.

*Ideal for facial procedures such as blepharoplasty procedures.*



## REGULAR TIPS

This is the most common tip style found on the market. Designed to provide versatility and a middle-ground between fine lines and thick lines.

*Ideal for most general procedures and for correct site marking.*



## BOLD TIP

Broad tip is unique to our 1404 and 1405 markers. This tip was specifically designed to replace gentian violet ink stamp pads. The paintbrush like tip coats instruments easily and completely.

*Ideal for coating surgical instruments with ink but especially cataract procedures.*

# ABOUT OUR MARKERS



## MINI MARKERS

Viscot Mini Markers are economical and eco-friendly. They're about half the size of a standard marker so they're about half the cost and produce half the waste

*Perfect for promoting single patient use.*



## MULTI-TIP BARRELS

Features two tips one on each end. Multi-Tip barrels differ depending on the exact product numbers but are an economical, convenient and eco-friendly alternative by providing two markers in one.

Usually featuring different tips, the unique Blephmarker® features two ultra fine tips, offering a fresh tip per eyelid.



## PRECISION BARREL

Our standard barrels designed with a comfortable grip in mind. Precision Barrels feel and perform like a fountain pen. Featured on our Precision markers and full size XL markers.



## VALUE BARREL

This barrel feels like a pencil when gripped. Value barrels are thinner and shorter than Precision and Wide barrels, but not as small as a Mini Marker. Less plastic means more cost savings.



## WIDE BARREL

Thicker barrels that feel more like a standard marker and less like a pen. All Wide Barrel Markers have a clip cap.

# CATARACT- TORIC LENS MARKING FOR THE INSTRUMENT

## NEW

### TOR-MARK-10-5

#### Disposable Cataract Marking Kit

Mini Max Regular Tip Traditional Ink Marker, Dial with Level and Indexed Cap 5/cs



### 1404/1405 PRECISION SURGICAL MARKER

Broad Tip | Precision Barrel

#### PART #

**1404-50**

**1404-100**

**1405-100**

#### DESCRIPTION

Sterile Broad Tip w/ Ruler

Sterile Broad Tip w/ Ruler

Non-Sterile Broad Tip

#### CASE QTY

50

100

100

# CATARACT- TORIC LENS MARKING DIRECTLY TO THE EYE



## 1424 TWIN TIP BARREL

Blephmarker: 2 Ultra Fine Tips | Traditional Gentian Violet Ink

PART #	DESCRIPTION	CASE QTY
1424-100	Sterile 2 Ultra Fine Tips	100
1424SR-10	Sterile 2 Ultra Fine Tips w/ Ruler	10
1424SR-100	Sterile 2 Ultra Fine Tips w/ Ruler	100



## 1436 VALUE SURGICAL MARKER

Ultra Fine Tip | Value Barrel

PART #	DESCRIPTION	CASE QTY
1436-100	Non Sterile Ultra Fine Tip	100
1436SR-10	Sterile Ultra Fine Tip w/ Ruler	10
1436SR-50	Sterile Ultra Fine Tip w/ Ruler	50
1436SR-100	Sterile Ultra Fine Tip w/ Ruler	100
1436SRL9-100	Sterile Ultra Fine Tip w/ Ruler & 9 Blank Labels	100



## 1402/1403 PRECISION SURGICAL MARKER

Ultra Fine Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1402-50	Sterile Ultra Fine Tip w/ Ruler	50
1402-100	Sterile Ultra Fine Tip w/ Ruler	100
1403-100	Non-Sterile Ultra Fine Tip	100

# CATARACT- TORIC LENS MARKING

BOTH



## 1425 TWIN TIP BARREL

Regular Tip & Ultra Fine Tip | Traditional Gentian Violet Ink

PART #	DESCRIPTION	CASE QTY
1425-100	Sterile 1 Regular & 1 Ultra Fine Tip	100
1425SR-100	Sterile 1 Regular & 1 Ultra Fine Tip w/ Ruler	100



**1451 MINI SURGICAL MARKER**  
Regular Tip | Mini Barrel

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1451-200</b>	Non-Sterile Regular Tip - 2 Plastic Containers of 100	200
<b>1451-1000</b>	Non-Sterile Regular Tip - 10 Plastic Containers of 100	1,000
<b>1451SR-100</b>	Non-Sterile Regular Tip w/ Ruler	100
<b>1451W-200</b>	Non-Sterile Regular Tip Individually Wrapped (Proof of Use) - 4 Plastic Containers of 50 Markers	200
<b>1451W-1000</b>	Non-Sterile Regular Tip Individually Wrapped (Proof of Use) - 20 Plastic Containers of 50 Markers	1,000



**1437 VALUE SURGICAL MARKER**  
Regular Tip | Value Barrel

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1437-100</b>	Sterile Regular Tip	100
<b>1437SR-50</b>	Sterile Regular Tip w/ Ruler	50
<b>1437SR-100</b>	Sterile Regular Tip w/ Ruler	100
<b>1437SRL9-50</b>	Sterile Regular Tip w/ Ruler & 9 Blank Labels	50
<b>1437SRL9-100</b>	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100



**1422 DUAL TIP SURGICAL INK MARKER & PERMANENT INK MARKER**  
Regular Tip + Traditional Ink | Fine Tip + Permanent Utility Ink

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1422-100</b>	Non-Sterile Regular Tip & Non-Sterile Fine Tip	100
<b>1422S-100</b>	Sterile Regular Tip & Non-Sterile Fine Tip	100
<b>1422SR-50</b>	Sterile Regular Tip & Non-Sterile Fine Tip w/ Ruler	50
<b>1422SR-100</b>	Sterile Regular Tip & Non-Sterile Fine Tip w/ Ruler	100
<b>1422SRL9-50</b>	Sterile Regular Tip & Non-Sterile Fine Tip w/ Ruler & 9 Blank Labels	50
<b>1422SRL9-100</b>	Sterile Regular Tip & Non-Sterile Fine Tip w/ Ruler & 9 Blank Labels	100





**1460 STANDARD TRADITIONAL INK SURGICAL MARKER**  
Regular Tip | Wide Barrel

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1460-100</b>	Non-Sterile Regular Tip	100
<b>1460SRL-100</b>	Sterile Regular Tip w/ Ruler & Labels	100
<b>1460SRL9-100</b>	Sterile Regular Tip w/ Ruler & 9 Labels	100



**1481 MINI GVFREE INK SKIN MARKER | PROP 65 COMPLIANT**  
Regular Tip | Mini Barrel | Gentian Violet Free

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1481GVF-200</b>	Non-Sterile Regular Tip - 2 Plastic Containers of 100	200
<b>1481GVF-1000</b>	Non-Sterile Regular Tip - 10 Plastic Containers of 100	1,000
<b>1481GVFSR-100</b>	Sterile Regular Tip w/ Ruler	100



**1487 VALUE GVFREE INK SKIN MARKER | PROP 65 COMPLIANT**  
Regular Tip | Value Barrel | Gentian Violet Free

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1487GVFSR-100</b>	Sterile Regular Tip w/ Ruler	100
<b>1487GVFSRL-100</b>	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100



## SMALL 1/2" YES SITE MARKERS

Cataract, Pediatric & Oncology Site Marking

PART #	DESCRIPTION	CASE QTY
<b>TYS-100W</b>	Non-Sterile Temporary 1/2" in Green "Yes" - 100 Markers w/ 100 Pre-Moistened Wipes	400
<b>TYS-400</b>	Non-Sterile Temporary 1/2" in Green "Yes" - 100 Per Box / 4Box / 1 Cs	400



## LARGE 1" YES SITE MARKERS

Cataract, Pediatric & Oncology Site Marking


PART #	DESCRIPTION	CASE QTY
<b>TYL-100W</b>	Non-Sterile Temporary 1" in Green "Yes" - 100 Markers w/ 100 Pre-Moistened Wipes	400
<b>TYL-400</b>	Non-Sterile Temporary 1" in Green "Yes" - 100 Per Box / 4Box / 1 Cs	400



**LASER TAT MARKERS**

For Ophthalmic Laser Site Marking

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>TEL-400</b>	Non-Sterile "Laser" in Red w/ White Background 1' Diameter - 100 Per Box / 4 Box / 1 Cs	400



**EYE IMPLANT TAT MARKERS**

For Cataract Implant Site Marking

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>TEY-400</b>	Non-Sterile Ophthalmic - 100 Per Box / 4 Box / 1 Cs	400



**RIGHT AND LEFT SITE MARKERS**

For Ophthalmic, Pediatric and Oncology Site Marking

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>TLR5050</b>	Non-Sterile 50 Green "Left" & 50 Blue "Right" -100 Per Box / 4 Box / 1 Cs	400



## 1424 TWIN TIP BARREL

Blephmarker: 2 Ultra Fine Tips | Traditional Gentian Violet Ink

PART #	DESCRIPTION	CASE QTY
1424-100	Sterile 2 Ultra Fine Tips	100
1424SR-10	Sterile 2 Ultra Fine Tips w/ Ruler	10
1424SR-100	Sterile 2 Ultra Fine Tips w/ Ruler	100



## 1425 TWIN TIP BARREL

Regular Tip & Ultra Fine Tip | Traditional Gentian Violet Ink

PART #	DESCRIPTION	CASE QTY
1425-100	Sterile 1 Regular & 1 Ultra Fine Tip	100
1425SR-100	Sterile 1 Regular & 1 Ultra Fine Tip w/ Ruler	100



## 1436 VALUE SURGICAL MARKER

Ultra Fine Tip | Value Barrel

PART #	DESCRIPTION	CASE QTY
1436-100	Non Sterile Ultra Fine Tip	100
1436SR-10	Sterile Ultra Fine Tip w/ Ruler	10
1436SR-50	Sterile Ultra Fine Tip w/ Ruler	50
1436SR-100	Sterile Ultra Fine Tip w/ Ruler	100
1436SRL9-100	Sterile Ultra Fine Tip w/ Ruler & 9 Blank Labels	100



## 1436 VALUE SURGICAL MARKER Ultra Fine Tip | Value Barrel

PART #	DESCRIPTION	CASE QTY
1436-100	Non Sterile Ultra Fine Tip	100
1436SR-10	Sterile Ultra Fine Tip w/ Ruler	10
1436SR-50	Sterile Ultra Fine Tip w/ Ruler	50
1436SR-100	Sterile Ultra Fine Tip w/ Ruler	100
1436SRL9-100	Sterile Ultra Fine Tip w/ Ruler & 9 Blank Labels	100



## 1402/1403 PRECISION SURGICAL MARKER Ultra Fine Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1402-50	Sterile Ultra Fine Tip w/ Ruler	50
1402-100	Sterile Ultra Fine Tip w/ Ruler	100
1403-100	Non-Sterile Ultra Fine Tip	100



## 1448 XL PREP RESISTANT INK SKIN MARKER Ultra Fine Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1448XL-50	Non-Sterile Ultra Fine Tip - 1 Plastic Container of 50	50
1448XLSR-100	Sterile Ultra Fine Tip w/ Ruler	100



## 1456 MINI MARKER XL PREP RESISTANT INK SKIN MARKER Ultra Fine Tip | Mini Barrel

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1456XL-200</b>	Non-Sterile Ultra Fine Tip - 2 Plastic Containers of 100	200
<b>1456XLSR-100</b>	Sterile Ultra Fine Tip w/ Ruler	100
<b>1456XLW-200</b>	Non-Sterile Ultra Fine Tip Individually Wrapped (Proof of Use) - 4 Plastic Containers of 50	200

# FEMTOSECOND LASER MARKING



## 1411 WATERPROOF PERMANENT MARKER

Fine Tip | Wide Barrel

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1411-100</b>	Non-Sterile Fine Tip	100
<b>1411S-100</b>	Sterile Fine Tip	100
<b>1411SRL9-50</b>	Sterile Fine Tip w/ Ruler & 9 Blank Labels	50
<b>1411SRL9-100</b>	Sterile Fine Tip w/ Ruler & 9 Blank Labels	100

# NON-SURGICAL AESTHETIC PROCEDURES



## 1440 EZ REMOVABLE WHITE INK SKIN MINI MARKER

Regular Tip | Mini Barrel | White Ink

PART #	DESCRIPTION	CASE QTY
1440-30	Non-Sterile Regular Tip - 1 Plastic Container	30
1440-240	Non-Sterile Regular Tip - 8 Plastic Containers of 30	240



## 1444 EZ REMOVABLE GREEN INK SKIN MINI MARKER

Regular Tip | Mini Barrel | Green Ink

PART #	DESCRIPTION	CASE QTY
1444-30	Non-Sterile Regular Tip - 1 Plastic Container	30
1444-240	Non-Sterile Regular Tip - 8 Plastic Containers of 30	240



## 1445 EZ REMOVABLE RED INK SKIN MINI MARKER

Regular Tip | Mini Barrel | Red Ink

PART #	DESCRIPTION	CASE QTY
1445-30	Non-Sterile Regular Tip - 1 Plastic Container	30
1445-240	Non-Sterile Regular Tip - 8 Plastic Containers of 30	240



# MEDICATION LABELS & OR ACCESSORIES

<b>AIR</b>	Mitomycin _____ mg/ mL
Antibiotic _____	Moxifloxacin 0.5% (Vigamox)
<b>BSS</b>	Perfluoron _____
<b>BSS</b>	SF <sub>6</sub> _____ %
C <sub>3</sub> F <sub>8</sub> _____ %	<b>Sterile Water</b>
Decadron _____ mg/mL	Triamcinolone PF _____ mg/mL
Dexamethasone _____ mg/mL	<b>Trypan Blue</b>
Lidocaine _____ % MPF	<b>Vitrase</b> _____ units/mL
<b>Miochol</b>	
<b>Miostat</b>	

## OPHTHALMIC COLOR MEDICATION LABEL KIT

PART #	DESCRIPTION	CASE QTY
<b>VEYE-20</b>	Sterile 20 Color Kit	100
<b>VEYE-22R-20</b>	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100



## 1453 WATERPROOF PERMANENT MARKING MINI MARKER

Regular Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
<b>1453-100</b>	Non-Sterile Regular Tip	100
<b>1453ARL9-100</b>	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100



**1422 DUAL TIP SURGICAL INK MARKER & PERMANENT INK MARKER**

Regular Tip + Traditional Ink | Fine Tip + Permanent Utility Ink

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1422-100</b>	Non-Sterile Regular Tip & Non-Sterile Fine Tip	100
<b>1422S-100</b>	Sterile Regular Tip & Non-Sterile Fine Tip	100
<b>1422SR-50</b>	Sterile Regular Tip & Non-Sterile Fine Tip w/ Ruler	50
<b>1422SR-100</b>	Sterile Regular Tip & Non-Sterile Fine Tip w/ Ruler	100
<b>1422SRL9-50</b>	Sterile Regular Tip & Non-Sterile Fine Tip w/ Ruler & 9 Blank Labels	50
<b>1422SRL9-100</b>	Sterile Regular Tip & Non-Sterile Fine Tip w/ Ruler & 9 Blank Labels	100



**1411 WATERPROOF PERMANENT MARKER**

Fine Tip | Wide Barrel

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1411-100</b>	Non-Sterile Fine Tip	100
<b>1411S-100</b>	Sterile Fine Tip	100
<b>1411SRL9-50</b>	Sterile Fine Tip w/ Ruler & 9 Blank Labels	50
<b>1411SRL9-100</b>	Sterile Fine Tip w/ Ruler & 9 Blank Labels	100



**1416 WATERPROOF PERMANENT MARKER**

Ultra Fine Tip | Wide Barrel

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>1416-100</b>	Non-Sterile Superfine Tip	100
<b>1416S-100</b>	Sterile Superfine Tip	100




**3" x 9" SOFFZONE NEUTRAL ZONE MAT \***

**PART #**  
**FNZ309RED**

**DESCRIPTION**  
3"x9" Red Tray w/ Raised Contours for Sharps

**CASE QTY**  
6



**SOFFZONE SYRINGE TRAY \***

**PART #**  
**FSY309RED**

**DESCRIPTION**  
3"x9" Red Tray w/ Tabs for Holding Syringes

**CASE QTY**  
6

[www.vision-surgical.com](http://www.vision-surgical.com)