



OPHTHALMOLOGY
www.Vision-Surgical.com

“Better Products, Better Value, Better Customer Service.”

Vision Surgical is an ophthalmic surgical (and some diagnostic) distribution company that specializes in Anterior Segment (Cornea, Glaucoma, Plastic & Cataract) and Posterior Segment instruments, knives, lasers, laser probes, cannulas and lenses. Better prices without lower quality. Most of our products are made in the U.S., Switzerland or Germany. Please visit our website at www.Vision-Surgical.com to see our complete line of products.



Owner's Great Grandmother
Elizabeth M. (Oberlender) Paden
WW1 Army Nurse Corps 1918

ABOUT US
www.Vision-Surgical.com

Terms of Sale

Ophthalmology



Terms of Sale

- This price list supersedes all price lists.
- All prices are subject to change without notice
- Our terms are Net 30 days to approved accounts
- List prices are subject to change
- All Vision Surgical instruments are guaranteed to be free of defects in materials and workmanship. Any instrument which proves defective in workmanship or material will be repaired or replaced at our discretion with out charge. This guarantee is void for instruments serviced by any company other than Vision Surgical instrument repair service.
- All Vision Surgical instruments are carefully inspected before shipment. Any claims must be made within (2) two days of receipt of shipment. Most orders are shipped same day and, at the latest, by UPS within 72 hours.

Return Policy

- With our prior written authorization, returns may be made within 30 days of purchase for full credit. Instruments must be unused and in their original package.
- To insure proper credit, you must request our "Return Authorization Form" and return the completed to us. A restocking and inspection charge of 20% will be charged on all authorized returns more than (30) days after invoice date. No returns will be accepted after 90 days. A restocking charge may be assess on any return that in our determination was not an error of Vision Surgical.
- A restocking and inspection charge of 35% will be charged for on all authorized returns more than (60) sixty days after our invoice date.
- Goods accepted for return will be on a merchandise credit basis only and we can only accept prepaid shipments.
- Any merchandise returned without approval will be subject to an automatic restocking fee with a minimum of 35%.

Procedure

1. Notify our Customer Service Department in writing or by phone of your request for return authorization.
2. We will mail a "Return Authorization Request" form. Please complete the form with the following information:
 - (A) Our stock number
 - (B) Invoice Number
 - (C) Date of invoice
 - (D) Reason for return
3. When your form has been approved, we will issue you an authorization to return and a label.
4. The following types of goods are not eligible for credit:
 - (A) Custom designed instruments
 - (B) Instruments which have been used, marked or altered in anyway
 - (C) Engraved instruments
 - (D) Discontinued products
 - (E) Sterile products
 - (F) Color coated instruments

Prices subject to change without notice. Vision Surgical is not accountable for any errors or omissions that may appear in the catalog.

forceps



OPHTHALMOLOGY

www.Vision-Surgical.com

Forcep Contents

Ophthalmology



Bipolar Forceps	1-9
Pencils Super Reusable	10
Titanium Bipolar Forceps	11-12
Bipolar Cords and Tray	13
Other Forcep Products	15-71



VS122-000

Jewelers Style Bipolar Forcep Straight

With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm. Overall Length 4" (101.6mm).



VS122-000-3

Jewelers Style Bipolar Forcep Straight

With A (High Polished 10mm Long Area) With A Tip Size Of 0.3mm. Overall Length 4" (101.6mm). Overall Weight 0.5 oz



VS122-200

Jewelers Style Bipolar Forcep Straight

Overall Length 4" (101.6mm) Insulated With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm. Overall Weight 0.6 oz



VS19-3161-1D

Jewelers Style Bipolar Forcep Straight

Insulated With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm. Overall Length 4.7" (119.38mm)

VS417827C

Jewelers Style Bipolar Forcep #5 Straight

Overall Length 4 3/4" (120.65mm) With A (High Polished 10mm Long Area) With A Tip Size 0.5mm. Overall Weight 0.6 oz



VS122-100

Jewelers Style Bipolar Forcep Curved In

Overall Length 4" (101.6mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm And Both Tips Angled In 15 Degrees Towards Each Other



Bipolar Forceps

Ophthalmology



VS122-100-3

Jewelers Style Bipolar Forcep Curved In

Overall Length 4" (101.6mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.3mm And Both Tips Angled In 15 Degrees Towards Each Other. Overall Weight 0.5 oz



VS122-150

Jewelers Style Bipolar Forcep Angled

Overall Length 4" (101.6mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm And Both Tips Angled 45 Degrees.



VS122-250

Jewelers Style Bipolar Forcep Angled

Overall Length 4" (101.6mm) Insulated With A (High Polished 10mm Long Area) With A Tip Size (0.4mm) And Both Tips Angled 45 Degrees



VS122-300

Tenzel Style Bipolar Forcep Straight

Overall Length 4" (101.6mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm And Both Tips Angled In Towards Each Other With A Spacer Stop To Create A 0.5mm Gap Between Tips. Overall Weight 0.5 oz



VS122-400

Jewelers Style Bipolar Forcep #7 Curved

Overall Length 4 3/4" (120.65mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm. Overall Weight 0.6 oz



Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS122-500

Adson Style Bipolar Forcep

Wide Serrated Handle Overall Length 4 3/4" (120.65mm) With A (High Polished 10mm Long Area) With A Tip Size Of 1.0mm. Overall Weight 0.8 oz



VS122-500-1.5

Adson Style Bipolar Forcep

Wide Serrated Handle Overall Length 4 3/4" (120.65mm) With A (High Polished 10mm Long Area) With A Tip Size Of 1.5mm. Overall Weight 0.8 oz



VS19-3161-11

Adson Style Bipolar Forcep

Wide Serrated Handle Overall Length 4 3/4" (120.65mm) With A (High Polished 10mm Long Area) With A Tip Size Of 1.0mm, Inside Jaws Is Serrated.



VS122-600

Adson Style Bipolar Forcep

Wide Serrated Handle Overall Length 4 3/4" (120.65mm) Insulated With A (High Polished 10mm Long Area) With A Tip Size Of 1.0mm



VS122-650

Adson Tissue Style Bipolar Forcep

Wide Serrated Handle Overall Length 4 3/4" (120.65mm) Insulated With A (High Polished 10mm Long Area) With 1 x 2 Teeth.

Bipolar Forceps

Ophthalmology



VS19-3161-11D

Adson Style Bipolar Forcep

Wide Serrated Handle Overall Length 4 3/4" (120.65mm) Insulated With A (High Polished 10mm Long Area) With A Tip Size Of 1.0mm, Inside Jaws Is Serrated



VS320-800

Iris Style Bipolar Forcep Straight

Overall Length 3 1/2" (88.9mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm. Overall Weight 0.4 oz



VS912317C

Iris Style Bipolar Forcep Straight

Overall Length 3 5/8" (92.08mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.6mm. Overall Weight 0.4 oz



VS320-850

Iris Style Bipolar Forcep Straight

Overall Length 3 1/2" (88.9mm) Insulated With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm. Overall Weight 0.4 oz



VS320-800

Iris Style Bipolar Forcep Angled

Overall Length 3 1/2" (88.9mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm Tip Is Angled 45 Degrees. Overall Weight 0.4 oz



VS912318C

Iris Style Bipolar Forcep Curved

Overall Length 3 5/8" (92.08mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.6mm. Overall Weight 0.4 oz

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS320-950

Iris Style Bipolar Forcep Angled

Overall Length 3 1/2" (88.9mm) Insulated With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm Tip Is Angled 45 Degrees.



VS912315C

Jewelers Style Bipolar Forcep Straight

Overall Length 3 5/8" (92.08mm) With A (High Polished 10mm Long Area) With A 0.5mm Angled Tips. Overall Weight 0.4 oz



VS19-3161-10A

Iris Style Bipolar Forcep Angled


Overall Length 3 1/2" (88.9mm) Insulated With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm Tip Is Angled 45 Degrees.



VS19-3161-10S

Mcpherson Style Bipolar Forcep Straight

Straight 3 1/2" (88.9mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.5mm And 4mm Tying Platforms. Overall Weight 0.4 oz



VS19-3161-14A

Shepard Style Bipolar Forcep Angled

Angled 4 1/2" (114.3mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.5mm And 4mm Tying Platforms, Tip Is Angled 45 Degrees. Overall Weight 0.6 oz



VS19-3161-14S

Shepard Style Bipolar Forcep Straight

Straight 4 1/2" (114.3mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.5mm And 4mm Tying Platforms. Overall Weight 0.5 oz



Bipolar Forceps

Ophthalmology



VS321-140

Nadler Style Coaptation Bipolar Forcep

With Removable Spacer Stop Overall Length 4 1/4" (107.95mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm. Overall Weight 0.5 oz

VS912310C

Nadler Style Coaptation Straight Bipolar Forcep

Overall Length 4" (101.60mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.25mm. Overall Weight 0.5 oz

VS321-120

Nadler Style Coaptation Bipolar Forcep

With Removable Spacer Stop Overall Length 4 1/4" (107.95mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.2mm. Overall Weight 0.5 oz

VS321-150

Nadler Style Coaptation Bipolar Forcep

With Removable Spacer Stop Overall Length 4 1/4" (107.95mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.5mm. Overall Weight 0.5 oz



VS321-020

Precision Style Coaptation Bipolar Forcep

Overall Length 4 1/4" (107.95mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.2mm, Tip Length From Bend To Tip Is 3mm Long. Overall Weight 0.5 oz

VS321-030

Precision Style Coaptation Bipolar Forcep

Overall Length 4 1/4" (107.95mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.3mm, Tip Length From Bend To Tip Is 3mm Long. Overall Weight 0.5 oz

VS321-040

Precision Style Coaptation Bipolar Forcep

Overall Length 4 1/4" (107.95mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm, Tip Length From Bend To Tip Is 3mm Long. Overall Weight 0.5 oz

VS321-050

Precision Style Coaptation Bipolar Forcep

Overall Length 4 1/4" (107.95mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.5mm, Tip Length From Bend To Tip Is 3mm Long. Overall Weight 0.5 oz

VS19-3161-9

Coaptation Bipolar Forcep Curved

With A Wide Serrated Handle Overall Length 4 3/5" (117mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm And 3.5mm Tying Platforms. Overall Weight 0.5 oz

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS321-240

Precision Style Coaptation Bipolar Forcep

Overall Length 4 1/4" (107.95mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm, Tip Length From Bend To Tip Is 2mm Long. Overall Weight 0.5 oz

VS321-220

Precision Style Coaptation Bipolar Forcep

Overall Length 4 1/4" (107.95mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.2mm, Tip Length From Bend To Tip Is 2mm Long. Overall Weight 0.5 oz

VS321-250

Precision Style Coaptation Bipolar Forcep

Overall Length 4 1/4" (107.95mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.5mm, Tip Length From Bend To Tip Is 2mm Long. Overall Weight 0.5 oz



VS19-3161-15

Osher Style Coaptation Bipolar Forcep

Overall Length 4 1/2" (114.3mm) With A (High Polished 10mm Long Area) With A Tip Size Of 0.4mm, Tip Length From Bend To Tip Is 2mm. Overall Weight 0.5 oz



VS19-3161-13

Cozean Style Colibri Bipolar Forcep

Overall Length 3 3/8" (85.73mm) With A (High Polished 10mm Long Area) With A 1.5mm Angled Tip With 0.12mm Teeth. Overall Weight 0.4 oz



VS19-3161-16

McCullough Style Jewelers Bipolar Forcep Straight

Overall Length 4 1/2" (114.3mm) With A (High Polished 10mm Long Area) With A Rounded Tip Size Of 1.0mm. Overall Weight 0.5 oz



VS19-3161-16I

McCullough Style Jewelers Bipolar Forcep Straight

Overall Length 4 1/2" (114.3mm) Insulated With A (High Polished 10mm Long Area) With A Rounded Tip Size Of 1.0mm.



Pencils Super Reusable

Ophthalmology



VS421-100

Reusable Bipolar Pencil
18 Gauge Straight Tip



421-120

Reusable Bipolar Pencil
18 Gauge With Wide Tip Curved 45 Degrees



VS421-140

Reusable Bipolar Pencil
20 Gauge With Straight Tip. Overall Weight 0.3 oz



421-160

Reusable Bipolar Pencil
23 gauge With Tapered Sharp Tip. Overall Weight 0.3 oz

421-170

Reusable Bipolar Pencil
23 Gauge With Tapered Blunt Tip.

421-180

Reusable Bipolar Pencil
18 Gauge With Wide Tip



VS520-100

Blue Titanium Jewelers Bipolar Forcep
3 1/2" (88.9mm) With A (High Polished Gold 10mm Long Area) With A 0.4mm Tip 30 Degree Angled. Overall Weight 0.3 oz



VS520-200

Blue Titanium Jewelers Bipolar Forcep Curved
3 1/2" (88.9mm) With A (High Polished Gold 10mm Long Area) With A 0.4mm Tip Curved In. Overall Weight 0.3 oz



Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS520-300

Blue Titanium Jewelers Bipolar Forcep Straight

3 1/2" (88.9mm) With A (High Polished Gold 10mm Long Area) With A 0.4mm Tip Straight. Overall Weight 0.3 oz



VS520-400

Blue Titanium Iris Angled Bipolar Forcep

3 3/4" (95.25mm) With A (High Polished Gold 10mm Long Area) With A 0.5mm Tip. Overall Weight 0.3 oz



VS520-500

Blue Titanium Jewelers Angled Bipolar Forcep

Angled 15 Degree 3 3/4" (95.25mm) With A (High Polished Gold 10mm Long Area) With A 0.4mm Tip.



VS520-600

Blue Titanium Coaptation Bipolar Forcep

4" (101.6mm) With A (High Polished Gold 10mm Long Area) With A 0.5mm Tip. Overall Weight 0.3 oz



VS520-700

Blue Titanium Jewelers Bipolar Forcep Straight

Wide Serrated Handle Straight Overall Length 4 3/4" (120.65mm) With A (High Polished Gold 10mm Long Area) With A 0.4mm Straight Tip. Overall Weight 0.5 oz



Titanium Forceps

Ophthalmology



VS520-800

Blue Titanium Jewelers Bipolar Forcep


4 1/4" (107.95mm) With A (High Polished Gold 10mm Long Area) With A 0.4mm Tip 90 Degree Angled. Overall Weight 0.3 oz



VS520-900

Blue Titanium Jewelers Bipolar Forcep

4 1/4" (107.95mm) With A (High Polished Gold 10mm Long Area) With A 0.4mm. Overall Weight 0.3 oz



VS521-000

Blue Titanium Jewelers Bipolar Forcep

With Fine Serrations 4 3/4" (120.65mm) With A (High Polished Gold 10mm Long Area) With A 0.4mm Tip.



VS521-100

Blue Titanium Jewelers Bipolar Forcep

4 1/2" (114.3mm) With A (High Polished Gold 10mm Long Area) With A 0.4mm Reversed Tip. Overall Weight 0.3 oz



VS521-200

Blue Titanium Jewelers Bipolar Forcep

4 1/2" (114.3mm) With A (High Polished Gold 10mm Long Area) With A 0.4mm. Overall Weight 0.3 oz



Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



Bipolar Cords
VS201102

Disposable Bipolar Cord Blue
10 Sterile Cords Per Box

VS201103

Reusable Bipolar Cord Blue Angled Instrument End
1 Non-Sterile Cords Per Box


VS201103-3

Reusable Blue Bipolar Cord
With Angled Instrument End Sold Non-Sterile 3 Per
Package

Bipolar Sterilization Tray

VST-4000

Bipolar Forcep Sterilization Tray With Lid
Holds 6 Forceps. Dimensions 9 5/8" x 5 1/8" x 1"
For Both Titan And MSi





Suture	15-22
Fixation	23-26
Capsulorhexis	27-31
Colibri	32-37
Tying	38-46
Pierson	47
Nucleus Cracker	48
Capsule	48
Dressing	49-50
Tissue	51-54
Jewelers	54-56
Chalazion	57-59
Muscle	60-61
Cilia	62-63
Muscle Clamps	64
Utility	65
Hemostat	66
Iris	67
IOL	67-69
Retinal	69
Foreign Body	70
Miscellaneous	71

VSTDK-3-221

Titanium Flat Handled Square Grip Suture Forcep
Straight 0.12mm 1 x 2 Teeth 6mm Platforms, Overall Length 85mm

VSTDK-3-221-1

Titanium Flat Handled Square Grip Suture Forcep
Straight 0.10mm 1 x 2 Teeth 6mm Platforms, Overall Length 85mm



VSTDK-3-222

Titanium Flat Handled Square Grip Suture Forcep
Curved 0.12mm 1 x 2 Teeth 6mm Platforms, Overall Length 83.5mm

VSTDK-3-224

Titanium Flat Handled Square Grip Suture Forcep
Curved 0.12mm 1 x 2 Teeth 6mm Platforms, Overall Length 113.5mm

VSTDK-3-225

Titanium Flat Handled Square Grip Suture Forcep
Straight 0.12mm 1 x 2 Teeth 6mm Platforms, Overall Length 115mm

VSTDK-3-227

Titanium Flat Handled Square Grip Suture Forcep
Straight 0.3mm 1 x 2 Teeth 6mm Platforms, Overall Length 115mm

VSTDK-3-228

Titanium Flat Handled Square Grip Suture Forcep
Straight 0.3mm 1 x 2 Teeth 6mm Platforms, Overall Length 115mm

VSTDK-3-253

Titanium Flat Handled Square Grip Suture Forcep
Angled 45 Degrees 0.12mm 1 x 2 Teeth 6mm Platforms 7mm Bend To Tip, Overall Length 82mm

VSTDK-3-255

Titanium Flat Handled Square Grip Suture Forcep
Angled 45 Degrees 0.12mm 1 x 2 Teeth 6mm Platforms 10mm Bend To Tip, Overall Length 82mm

VSTDK-3-259

Titanium Flat Handled Square Grip Suture Forcep
Angled 45 Degrees 0.12mm 1 x 2 Teeth 6mm Platforms 10mm Bend To Tip, Overall Length 111mm

VSTDK-3-268

Titanium Flat Handled Square Grip Suture Forcep
Angled 45 Degrees 0.12mm 1 x 2 Teeth 6mm Platforms 7mm Bend To Tip, Overall Length 112mm

VS5-2894

Castroviejo Suturing Forcep
Wide Serrated Handle With Polished Finish, Bayonet Shafts, 0.3mm Teeth, 6mm Tying Platform, And Overall Length Of 4 1/4" (106mm)

Suture Forceps

Ophthalmology



VS5-2896

Castroviejo Suturing Forcep

Castroviejo Suturing Forcep, Wide Serrated Handle With Polished Finish, Bayonet Shafts, 0.5mm Teeth, And 6mm Tying Platform, And Overall Length Of 4 1/4" (11 cm)



VS5-2899

Castroviejo Suturing Forcep

Castroviejo Suturing Forcep, Wide Serrated Handle With Polished Finish, 40 Degree Angled Shafts, 0.3mm Teeth, 6mm Tying Platform, And Overall Length Of 4 1/8" (10.5cm)



VS5-2807

Castroviejo Suturing Forcep

Castroviejo Suturing Forcep, Wide Serrated Handle With Polished Finish, Straight Shafts, 0.12mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 3/8" (11 cm). Overall Weight Of 0.5 oz



VS5-2807-S

Sanders-Castroviejo Suturing Forcep

Wide Serrated Handle With Polished Finish, Straight And Heavy Shafts, 0.12mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 1/4" (11 cm). Overall Weight Of 0.4 oz



VS5-2807-T

Castroviejo Suturing Forcep

Titanium, Wide Serrated Handle With Polished Finish, Straight Shafts, 0.12mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 1/4" (11 cm). Overall Weight Of 0.5 oz



Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS11-6-3611

Castroviejo Suturing Forcep


0.12mm 1 x 2 Teeth On A Three-Holed Handle, 20mm Straight Shafts, And An Overall Length Of 4 1/4" (11cm)



VS11-6-3611

Castroviejo Titanium Forcep

Three-holed handle, Straight and Sharp 17mm Shafts, 0.12mm Teeth, And Overall Length Of 4 1/4" (11cm).



VS11-6-3661

Snider-Castroviejo Suture Forcep

On Three Hole 0.12mm 1 x 2 Teeth With Tying Platform And Guide Pin



VS5-2808

Castroviejo Suturing Forcep

Wide Serrated Handle With Polished Finish, Straight Shafts, 0.3mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 3/8" (11cm). Overall Weight Of 0.4 oz



VS5-2808-T

Castroviejo Suturing Forcep

Wide Serrated Handle With Polished Finish, Straight Shafts, 0.3mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 1/4" (11cm)



VS11-6-3621

Castroviejo Suturing Forcep

0.3mm with Three-Holed Handle, 20mm Straight Shafts, 1x2 Teeth, And Overall Length Of 4 1/4" (11cm).



VS11-6-3621-T

Castroviejo Suturing Forcep

Titanium, 0.3mm with Three-Holed Handle, 20mm Straight Shafts, 1x2 Teeth, And Overall Length Of 4 1/4" (11cm).

Suture Forceps

Ophthalmology



VS5-2809

Castroviejo Suturing Forcep

Wide Serrated Handle With Polished Finish, Straight Shafts, 0.5mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 3/8" (11 cm). Overall Weight Of 0.5 oz

VS5-2809-T

Castroviejo Suturing Forcep

Titanium, Wide Serrated Handle With Polished Finish, Straight Shafts, 0.5mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 3/8" (11 cm)



VS11-6-3631

Castroviejo Suturing Forcep

0.5mm with Three-Holed Handle, 20mm Straight Shafts, 1x2 Teeth, Tying Platform, And Overall Length Of 4 1/4" (11 cm).

VS11-6-3631-T

Castroviejo Suturing Forcep

Titanium, 0.5mm with Three-Holed Handle, 20mm Straight Shafts, 1x2 Teeth, Tying Platform, And Overall Length Of 4 1/4" (11 cm).



VS5-2809-S

Stern-Castroviejo Suturing Forcep

Wide Serrated Handle With Thumb Catch Lock And Polished Finish, Straight Shafts, 1 x 2 (0.5mm) Teeth, 6mm Tying Platform, And Overall Length Of 4 1/2" (11.5cm). Overall Weight Of 0.5 oz



VS11-6-3663

Stern-Castroviejo Suturing Forcep

With Three-Holed Handle, 1 x 2 (0.5mm) Teeth, Tying Platform With Lock, Straight Shafts, And An Overall Length Of 4 1/4" (11 cm)



VS11-6-3641

Castroviejo Suturing Forcep

0.9mm with Three-Holed Handle, 20mm Straight Shafts, 1x2 Teeth, And Overall Length Of 4 1/4" (11cm)



VS11-6-3211

Thorpe Corneal Fixation Forcep

On A 3 Hold Handle With 0.12mm, 2 x 3 Teeth With 6mm Tying Platform



VS5-2800

Castroviejo Suturing Forcep

Wide Serrated Handle With Polished Finish, Straight Shafts, 1.5mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 1/2" (11.5cm). Overall Weight Of 0.4 oz



VS5-2911

Castroviejo Suturing Forcep

Serrated Handle With Polished Finish, Straight Shafts, 0.5mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 1/4" (10.9cm). Overall Weight Of 0.5 oz



VS5-2913

Castroviejo Suturing Forcep

Serrated Handle With Polished Finish, Straight Shafts, 0.9mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 1/4" (10.9cm)

VS5-2916

Fine-Castroviejo Suturing Forcep

Serrated Handle With Polished Finish, Straight Shafts, 0.22mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 1/4" (10.9cm)



VS5-2916-XD

Fine-Castroviejo Suturing Forceps

Serrated Handle With Polished Finish, Straight Shafts, 0.12mm Teeth, 6mm Tying Platform, And Overall Length Of 3 3/4" (9.5cm). Overall Weight Of 0.3 oz



VS5-2916-XDWH

Fine-Castroviejo Suturing Forcep

Wide Serrated Handle With Polished Finish, Straight Shafts, 0.12mm Teeth, 6mm Tying Platform, And Overall Length Of 4 1/4" (11cm)

Suture Forceps

Ophthalmology



VS5-2917

Fine-Castroviejo Suturing Forcep

Polished Finish, 45 Degree Angled Shafts, 0.5mm Teeth, 6mm Tying Platform, And Overall Length Of 4 1/4" (11cm)



VS5-2910

Bonn Suturing Forceps

Serrated Handle With Polished Finish, Straight Shafts, 0.12mm Teeth, Tying Platform, And Overall Length Of 3" (8cm)



VS11-6-3111

Bonn Suture Forcep

Short Fat 3 Hole Handle Design Delicate 1 x 2 Teeth Without Tying Platforms



VS5-2922

Precision Bonn Suturing Forcep

Serrated Handle With Polished Finish, Straight Shafts, 0.12mm Teeth, 5mm Tying Platform, And Overall Length Of 3 3/4" (9.5cm). Overall Weight Of 0.3 oz



VS11-6-3131

Bonn Suture Forcep

Long Design, 3 Hole Handle Design, Straight 20mm Shafts, Delicate with 1x2 Teeth, Tying Platforms, And Overall Length Of 4 1/4" (11cm)

VS5-2924-S

McPherson Straight Suturing-Corneal Forceps

Serrated Handle with Dull Finish, Straight Shafts, 1x2 Teeth, 5mm Tying Platform, Overall Length Of 3 3/8" (87mm)

VS11-6-3411

McPherson Straight Suture Forcep

.44 mm On a 3-Hole Handle Design With 5mm Tying Platforms 1 x 2 Teeth



VS5-2924-A

McPherson Angled Suturing–Corneal Forcep

Serrated Handle With Dull Finish, 45 Degree Angled Shafts, 1 x 2 Teeth, 5mm Tying Platform, And Overall Length Of 3 1/4" (8.3cm). Overall Weight Of 0.2 oz

VS11-6-3421

McPherson Angled Suture Forcep

On a 3–Hole Handle Design With Tying Platforms 7mm Bend To Tip Angled 45 Degree 0.4mm 1 x 2 Teeth

VS5-2924-AK

Kelman–McPherson Angled Suturing–Corneal Forcep

Serrated Handle With Dull Finish, 45 degree Angled Shafts 8mm Long, 1 x 2 Teeth, 5mm Tying Platform, And Overall Length Of 3 1/4" (85mm)

VS11-6-3441

Mcpherson Angled Suture Forcep

On a 3–Hole Handle Design With 6mm Tying Platform 8mm Bend To Tip Angled 45 Degree 0.3mm 1 x 2 Teeth

VS11-6-3461

Hoeffler Angled Suture Forcep

On a 3–Hole Handle Design With 6mm Tying Platform 11mm Bend To Tip Angled 45 Degree 0.3mm 1 x 2 Teeth



VS5-2927

Berens Curved Suturing Forcep

Serrated Handle With Polished Finish, Heavy Curved Shafts, 2 X 2 Teeth, Crisscross Serrated 6mm Tying Platforms, And Overall Length Of 4 1/4" (11cm). Overall Weight Of 0.5 oz

VS5-2935

McCullough Suturing Forcep

Serrated Handle With Polished Finish, Straight Shafts, 6mm Crisscross Serrated Tying Surfaces, And Overall Length Of 3 5/8" (92mm)

VS5-2939

Precision Fixation / Suture Forcep

1 by 2 Teeth On A Thin Serrated Handle Overall Length 4 3/8"



VS5-2941

Sauer Suturing Forcep

Wide Serrated Handle With Polished Finish, Straight Shafts, 1 x 2 Teeth, And Overall Length Of 3 3/8" (8.7cm). Overall Weight Of 0.3 oz

Suture Forceps

Ophthalmology



VS5-2942

Paufique Suturing Forceps


Serrated And Wide Handle With Polished Finish, Straight Shafts, 1 x 2 Teeth, 6mm Tying Platform, And Overall Length Of 3 3/8" (8.7cm). Overall Weight Of 0.3 oz



VS5-2943

ManHattan Eye And Ear Suturing Forcep

Smooth And Wide Handle With Polished Finish, Straight Shafts, 1 x 2 Teeth, 6mm Tying Platform, And Overall Length Of 3 1/2" (9cm). Overall Weight Of 0.3 oz



VS5-2992-12C

O'Gawa Curved Suture-Fixation Forcep-0.12mm

Round Handle In Closed Position With Polished Finish, Curved Shafts, 5.5mm Tying Platform, 0.3mm Gap At Platform Heel, Tips With 1 x 2 Teeth (0.12mm), And Overall Length Of 4 1/4" (108mm)



VS5-2992-12S

O'Gawa Straight Suture-Fixation Forcep-0.12mm

Round Handle In The Closed Position, Straight Shafts, 5.5mm Tying Platform, 0.3mm Gap At Platform Heel, Tips With 0.12mm Teeth, And Overall Length Of 4 1/4" (10.8cm). Overall Weight Of 0.5 oz



VS5-2992-3S

O'Gawa Suture And Fixation Forcep-0.3mm

Round Handle With Polished Finish, Curved Shafts, 5.5mm Tying Platform, 0.3mm Gap At Platform Heel, Tips With 0.3mm Teeth, And Overall Length Of 4 1/4" (108mm)

VS5-2992-S

O'Gawa Straight Suture-Fixation Forcep

Round And Knurled Handle In The Closed Position, Straight Shafts, 5.5mm Tying Platform, 0.3mm Gap At Platform Heel, And Overall Length Of 4 1/4" (108mm)

VS5-2992-PS

O'Gawa Straight Suture-Fixation Forcep-0.12mm

Round And Knurled Handle In The Closed Position With Polished Finish, Straight Shafts, 5.5mm Tying Platform, 0.3mm Gap At Platform Heel, And Overall Length of 4 1/4" (108mm)



VS5-2990-S

Paton Suture Forcep

Extra Delicate Straight Sharp Pointed Tips, 4.5mm Tying Jaw, With A Castroviejo Grip Handle, With An Overall Length Of 3 1/2" (9cm). Overall Weight Of 0.4 oz



Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020

VSTDK-3-211

Pierce Titanium Suture Forcep

Flat Handled Square Grip Straight 0.25mm Notched Tip With 6mm Tying Platforms, Overall Length 85mm

VSTDK-3-212

Pierce Titanium Suture Forcep

Flat Handled Square Grip Curved 0.25mm Notched Tip With 6mm Tying Platforms, Overall Length 85mm

VSTDK-3-214

Pierce Titanium Suture Forcep

Flat Handled Square Grip Curved 0.25mm Notched Tip With 6mm Tying Platforms, Overall Length 114mm

VSTDK-3-215

Pierce Titanium Suture Forcep

Flat Handled Square Grip Straight 0.25mm Notched Tip With 6mm Tying Platforms, Overall Length 115mm

VSTDK-3-217

Precision Rabkin Style Titanium Suture Forcep

Flat Handled Square Grip Straight 0.7mm Notched Tip With 6mm Tying Platforms, Overall Length 115mm



VS5-2804

Castroviejo Fixation Forcep

Wide Serrated Handle With Polished Finish, Straight Shafts, 1 x 2 (0.3mm) Teeth, And Overall Length Of 4 1/4" (11cm). Overall Weight Of 0.5 oz



VS5-2806

Castroviejo Fixation Forcep

Wide Serrated Handle With Polished Finish, Straight Shafts, 1 x 2 (0.5mm) Teeth, And Overall Length Of 4 1/2" (11.5cm). Overall Weight Of 0.5 oz



VS5-2664

Moody Curved Fixation Forcep

Serrated Handle With Slide Lock And Polished Finish, Curved Shafts, 1 x 2 0.5mm Teeth, And Overall Length Of 4 1/4" (108mm)

VS11-6-3664

Moody Curved Left Fixation Forcep

On A 3 Hole Handle With Lock And Polished Finish, Curved Shafts, 0.5mm, 1 x 2 Teeth, And Overall Length Of 4 1/4" (108mm)

VS11-6-3665

Moody Curved Right Fixation Forcep

On A 3 Hole Handle With Lock And Polished Finish, Curved Shafts, 0.5mm, 1 x 2 Teeth, And Overall Length Of 4 1/4" (108mm)

Fixation Forceps

Ophthalmology



VS5-2843

Thorpe Corneal Fixation Forcep

Serrated Handle With Polished Finish, Straight Shafts, 45 Degree 2 X 3 (0.12mm) Teeth, Tying Platform, And Overall Length Of 3 3/4" (94mm)

VS5-2841

Thorpe Conjunctival Fixation Forcep

Serrated Handle With Polished Finish, Straight Shafts, 90 Degree 2 X 3 Teeth, Tying Platform, And Overall Length Of 3 3/4" (94mm)

VS5-2828

Osher Conjunctival Fixation Forcep

Serrated And Wide Handle With Polished Finish, Straight Shafts, 45 Degree, 1 x 2 (0.3mm) Teeth, 5.5mm Tying Platform, And Overall Length Of 4 1/4" (108mm)

VS5-2675

Osher 0.12mm Forcep–Elschnig Teeth

Delicate Tying Platform, 1 x 2 (0.12mm) Teeth, And Overall Length Of 4 1/8" (104mm)

VS5-2912

Thorton Angled Fixation Forcep

With Wide Serrated Handle With Polished Finish, Straight Shafts With 25 Degrees Angled 1.5mm Tips, 4mm Platform With 0.12mm 1 x 2 Teeth, And An Overall Length Of 4.3" (109mm)



VS5-2835

Iowa State Fixation Forcep

Flat Handle With Polished Finish, Straight Shafts, 1 x 2 Teeth, And Overall Length Of 3 3/8" (8.6cm). Overall Weight Of 0.2 oz



VS5-2795

Elschnig Fixation Forceps

Serrated And Narrow Handle With Polished Finish, Straight Shafts, 1 x 2 (2.0mm) Elschnig Style Teeth, And Overall Length Of 4 3/8" (11 cm). Overall Weight Of 0.3 oz



VS5-2796

Elschnig Fixation Forcep

Wide Serrated Handle With Polished Finish, Straight Shafts, 1 x 2 (2mm) Elschnig Style Teeth, And Overall Length Of 4 1/4" (11 cm). Overall Weight Of 0.4 oz

VS5-3367

O'Connor–Elschnig Fixation Forcep

Serrated And Flat Handle With Slide Style Lock And Polished Finish, Straight Shafts, Heavy Elschnig 1 x 2 Teeth, And Overall Length Of 4 1/4" (109mm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS5-2765

Lester Fixation Forcep

Serrated Handle With Polished Finish, Straight Shafts, Fine 1 x 2 Teeth, And Overall Length Of 3 1/2" (9cm). Overall Weight Of 0.2 oz

VS5-2767

Lester Fixation Forcep

Serrated Handle With Polished Finish, Straight Shafts, Fine 2 X 3 teeth, And Overall Length Of 3 3/4" (96mm)

VS5-2731

Kuhnt Fixation Forcep

Serrated Handle With Polished Finish, Straight Shafts, Curved 1 x 2 Teeth, And Overall Length Of 4 1/4" (107mm)



VS5-2781

Bracken Fixation Forcep

Serrated Handle With Polished Finish, Straight Shafts, Delicate 1 x 2 (0.7mm) Teeth, And Overall Length Of 4 1/4" (11cm). Overall Weight Of 0.3 oz



VS5-2781-FN

Bracken-Farkas Corneal Forcep

With Neumann Modification, Serrated Handle, Straight Shafts, Delicate 1 x 2 (0.12mm) Teeth, Stabilization Pin, And Overall length Of 4 1/2" (11cm) . Overall Weight Of 0.3 oz

VS5-2782

Bores Corneal Fixation Forcep

Wide Serrated Handle With Dull Finish, Straight Shafts With Two Point Fixation, Tips Separated By 3mm, Each Tip Has 1 x 2 (0.12mm) Teeth, And Overall Length Of 4 1/4" (110mm)

VS5-2785

O'Brien Fixation Forcep

Serrated And Wide Handle With Polished Finish, Straight Shafts, 1 x 2 Elschmig Style Teeth, And Overall Length Of 4 1/4" (109mm)

Fixation Forceps

Ophthalmology



VS5-2671

Graefe Fixation Forcep

Serrated Handle Without Slide Lock But With Polished Finish, Narrow 3.5mm Jaws, Fine Teeth, And Overall Length Of 4 1/4" (10.9cm). Overall Weight Of 0.3 oz



VS5-2673

Graefe Fixation Forcep

Serrated Handle Without Slide Lock But With Polished Finish, Wide 4.5mm Jaws, Fine Teeth, And Overall Length Of 4 1/4" (108mm)



VS5-2683

Graefe Fixation Forcep

Serrated Handle With Slide Lock But With Polished Finish, Wide 4.5mm Jaws, Fine Teeth, And Overall Length Of 4 1/4" (108mm)



VS5-2711

Green Fixation Forcep

Serrated Handle With Polished Finish, Narrow 12mm T-Shaped Jaws, Fine Teeth, And Overall Length Of 4" (10cm). Overall Weight Of 0.5 oz



VS5-2715

Green Fixation Forcep

Serrated Handle With Polished Finish, Narrow 5mm T-Shaped Jaws, Fine Teeth, And Overall Length Of 4" (10cm). Overall Weight Of 0.5 oz

VS5-2731

Kuhnt Fixation Forcep

Serrated Handle With Polished Finish, Straight Shafts, Curved 1 x 2 Teeth, And Overall Length Of 4 1/4" (107mm)

VS5-2782

Bores Corneal Fixation Forcep

Wide Serrated Handle With Dull Finish, Straight Shafts With Two Point Fixation, Tips Separated By 3mm, Each Tip Has 1 x 2 (0.12mm) Teeth, And Overall Length Of 4 1/4" (110mm)

VS11-6-4361

Bores Corneal Fixation Forcep

Wide Serrated Handle With Dull Finish, Straight Shafts With Two Point Fixation, Tips Separated By 3mm, Each Tip Has (0.3mm) 1 x 2 Teeth, And Overall Length Of 4 1/4" (110mm)

VS5-2770

Kremer Fixation Forcep

Serrated And Wide Handle With Slide Lock And Dull Finish, Two Point Fixation, Tips Separated By 13mm, Each Tip Has 1 x 2 0.12mm Teeth, Straight Shafts, And Overall Length Of 4 1/8" (105mm)

VS5-2770-A

Kremer Angled Fixation Forcep

Serrated And Wide Handle With Slide Lock And Dull Finish, Two Point Fixation, Tips Separated By 13mm, Each Tip Has 1 x 2 (0.12mm) Teeth, Angled Shafts, And Overall Length Of 4 1/8" (105mm)



VSPIU-2

Original Kraff Capsulorhexis Utrata

With An Overall Length Of 4" (10.1cm)



VS5-3113

Utrata Capsulorhexis Forcep

Flat Handle With Dull Finish, Fine Angled Tips, Angled Shaft, 12mm From Bend To Tip, With Iris Stop, And Overall Length Of 3 1/4" (8.3cm)



VS5-3113-S

Utrata Capsulorhexis Forcep

Flat Handle With Dull Finish, With Fine Angled Sharp Cystotome Tips, 12mm Angled Shafts With Iris Stop, And Overall Length Of 3 1/4" (80mm)

Capsulorhexis Forceps

Ophthalmology



VS5-3113-C

Lehner-Utrata Capsulorhexis Forcep

Flat Handle With Dull Finish, Utrata-Like Fine Angled Tips, Gently Curved Shaft, 12mm From Bend To Tip, And Overall Length Of 3 1/4 (8.3cm). Overall Weight Of 0.2 oz



VS5-3113-G

Precision Gimbel Style Capsulorhexis Forcep

Flat Handle With Dull Finish, Similar In Construction To Utrata Capsulorhexis Forcep, And Overall Length Of 3 1/4" (80mm)



VS5-3113-R

O'Gawa-Utrata Capsulorhexis Forcep

Round Handle With Dull Finish, Fine Angled Tips, Angled Shaft, 12mm From Bend To Tip, And Iris Stop, And Overall Length Of 4" (10cm)



VS5-3113-RC

Precision Lehner Style Utrata Capsulorhexis Forcep

Round Handle With Dull Finish, Utrata-Like Fine Angled Tips, Gently Curved Shaft, 12mm From Bend To Tip, And Overall length Of 4 3/8" (112mm)



VS5-3115

Kellan Capsulorhexis Forcep

Round Handle With Dull Finish, Delicate 45 Degree Angled Shaft, 12mm From Bend To Tip Arruga Style Tips, And Overall Length Of 3 7/8" (10cm). Overall Weight Of 0.4 oz

VS5-3116

Precision Lehner Style Capsulorhexis Forcep

Flat Handle With Dull Finish, Gently Curved Shaft, 2mm Platform, And Overall Length Of 3 1/4" (80mm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS28274M

Lerner Croy Utrata

11.5mm Vaulted Shanks With A Cystotome Tip, Ergonomic Fenestrated Handle



VS11-6-6181

Precision Style Inamura Capsulorhexis Forcep

On A 3 Hole Spring Handle Design Very Delicate Short Tapered Shanks

VS11-6-6183

Precision Style Inamura Capsulorhexis Forcep

On A 3 Hole Spring Handle Design Very Delicate Long Tapered Shanks



VS11-6-6191

Kraff Capsule Tag Forcep

With 3 Hole Handle, Long Extremely Thin Shafts Angled 11mm From Bend To Tip, And An Overall Length Of 4 1/4" (11cm)



VS11-6-6192

Utrata Capsulorhexis Forcep

With 3 Hole Handle, Long Extremely Thin Shafts Angled 11mm From Bend To Tip, Small Bent Tips, And An Overall Length Of 4" (10cm)



VS11-6-6192-T


Titanium Utrata Capsulorhexis Forcep

With 3 Hole Handle, Long Extremely Thin Shafts Angled 11mm From Bend To Tip, Small Bent Tips, And An Overall Length Of 4" (10cm)



Capsulorhexis Forceps

Ophthalmology



VS11-6-6193

Utrata Capsulorhexis Forcep

With Round Handle With Longitudinal Grooves, Long Extremely Thin Shafts Angled 11 mm From Bend To Tip, Small Bent Tips, And An Overall Length Of 4" (10cm)

VS11-6-6193-V

Utrata Capsulorhexis Forcep

With Round Handle With Longitudinal Grooves, Long Extremely Thin Shafts Angled 11 mm From Bend To Tip Vaulted, Small Bent Tips, And An Overall Length Of 4" (10cm)

VS11-6-6194

Precision Nevyas Style Utrata Capsulorhexis Forcep

With 3 Hole Handle, Long Extremely Thin Shafts Angled 11mm From Bend To Tip, Small Sharp Cystotome Tips, And An Overall Length Of 4" (10cm)

VS11-6-6194-T

Titanium Nevyas Ca[pularhexis Forcep

On A 3 Hole Handle With Extremely Thing Shafts, 11 mm From Tip To Bend, Cystotome Shaped Tips, And An Overall Length Of 4" (10cm)



VS11-6-6195

Precision Masket Style Utrata Capsulorhexis Forcep

With 3 Hole Handle, Long Extremely Very Delicate Thin Shafts Angled 11mm From Bend To Tip Vaulted, Small Bent Triangular Tips, And An Overall Length Of 4" (10cm)

VS11-6-6197

Precision Masket Style Utrata Capsulorhexis Forcep

With Round Handle, Long Extremely Very Delicate Thin Shafts Angled 11 mm From Bend To Tip Vaulted, Small Bent Triangular Tips, And An Overall Length Of 4" (10cm)

VS11-6-6199

Precision Jaffe Style Utrata Capsulorhexis Forcep

With 3 Hole Handle, Long Extremely Very Delicate Thin Shafts Angled 11 mm From Bend To Tip Small Bent Round Tips, And An Overall Length Of 4" (10cm)

VS11-6-8722

Titanium Micro Capsulorhexis Forcep

20 Gauge Curved Shaft With Micro Triangular Tips

VS11-6-8723

Titanium Micro Capsulorhexis Forcep

20 Gauge Straight Shaft With Micro Triangular Tips

VS11-6-8762

Titanium Micro Capsulorhexis Forcep

23 Gauge Curved Shaft With Micro Grasping Tips

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020

VS15-5406-ST

Titanium Capsulorhexis Forcep

On A Round Handle With Round Grooves For Grip With Extremely Thin Shanks, 11mm From Tip To Bend, Cystotome Shaped Tips, And An Overall Length Of 4" (10cm)

VS16-3433-R

Titanium Precision Kershner Style Capsulorhexis Forcep

On A Diamond Knurling Round Handle With Extremely Thin Shanks, 11.5mm From Tip To Bend, Cystotome Shaped Tips, And An Overall Length Of 4" (10cm)

VS16-3434-R

Titanium Precision Kershner Style Capsulorhexis Forcep

On A Diamond Knurling Round Handle With Extremely Thin Vaulted Shanks, 11.5mm From Tip To Bend, Cystotome Shaped Tips, And An Overall Length Of 4" (10cm)

VSTDK-3-820

Titanium Utrata Flat Handled Square Grip

With 8.5mm Platforms, 45° angled shafts, 13mm Bend To Tip Overall Length 83mm

VSTDK-3-820-1

Titanium Utrata Flat Handled Square Grip

With 8.5mm Platforms, 45° angled shafts, 13mm Bend To Tip Overall Length 114mm

VSTDK-3-814

Titanium Precision Bloomberg Style

Utrata Flat Handled Square Grip Without Platforms, 45° angled shafts, 14mm Bend To Tip Angled 45 Degree Tip (0.12) 1 x 2 Teeth Overall Length 83mm

VSTDK-3-825

Titanium Utrata

Round Handled Square Grip With 8.5mm Platforms, 45° angled shafts, 13mm Bend To Tip Overall Length 115mm

VSTDK-3-826

Titanium Precision Bloomberg Pierce Style Utrata

Flat Handled Square Grip Without Platforms, 45° angled shafts, 14mm Bend To Tip Angled 45 Degree Tip (0.25) NotchTeeth Overall Length 84mm



VSTDK-3-827

Titanium Precision Inamura Style Utrata

Flat Spring Handled Square Grip With Cross Action Design Without Platforms, 45° Angled Slightly Curved Shafts, 14mm Bend To Tip Overall Length 105mm

VSTDK-3-827-R

Titanium Precision Inamura Style Utrata

Round Spring Handled With Cross Action Design Without Platforms, 45° Angled Slightly Curved Shafts, 14mm Bend To Tip Overall Length 126mm

VSTDK-3-829

Titanium Utrata

Flat Handled Square Grip With 8.5mm Platforms, 45° Angled Curved Shafts, 13mm Bend To Tip Overall Length 83mm



VSTDK-3-243

Titanium Troutman Barraquer Colibri Forcep

Flat Handled Square Grip (.12mm) 1 x 2 Teeth 6mm Platforms, 2.0mm Tip Length, Overall Length 84mm

VSTDK-3-246

Titanium Troutman Barraquer Colibri Forcep

Flat Handled Square Grip (.12mm) 1 x 2 Teeth 6mm Platforms, 2.0mm Tip Length, Overall Length 115mm

VSTDK-3-554

Titanium Precision Nordan Style Colibri Forcep

Round Handled With Small Flat Grip (.12mm) 1 x 2 Teeth 6mm Platforms, 2.0mm Tip Length, Overall Length 115mm

VSTDK-3-241

Titanium Precision Pierse Style Colibri

Flat Handled Square Grip, With 6.0mm Tying Platforms, 45° angled shafts Angled 45 Degree Tip (0.25) Notch Teeth 2.0mm Tip Length Overall Length 84mm

VSTDK-3-247

Titanium Precision Pierse Style Colibri

Flat Handled Square Grip, With 6.0mm Tying Platforms, 45° angled shafts Angled 45 Degree Tip (0.25) Notch Teeth 2.0mm Tip Length Overall Length 115mm



VS5-4200

Troutman-Barraquer Corneal Utility Forcep

Colibri Style, Serrated Handle With Polished Finish, 1 x 2 (0.12mm) Teeth, And Overall Length Of 3" (7.5cm). Overall Weight Of 0.2 oz



VS5-4201

Katzin-Barraquer Corneal Utility Forcep

Colibri Style, Serrated Handle With Polished Finish, 1 x 2 (0.4mm) Teeth, 7mm Tying Platform, And Overall Length Of 2 7/8" (73mm)



VS5-4202

Troutman-Barraquer Corneal Utility Forcep

Colibri Style, Serrated Handle With Polished Finish, 1 x 2 (0.3mm) Teeth, 6mm Tying Platform, And Overall Length Of 2 7/8" (73mm)

VS5-4315

Harms Corneal Utility Forceps

Colibri Style, Serrated Handle With Polished Finish, 1 x 2 (0.12mm) Teeth Set At 45 Degrees, 4.5mm Tying Platform, And Overall Length Of 2 3/4" (70mm)

VS11-6-2561

Gills Colibri Forcep

Short Handle Design With 17mm x 5.5mm Finger Grip Cut-Out Very Delicate Fine Sharp Pointed Tips With Tying Platforms Overall Length Of 2.9" (75.6mm)



VS11-6-2611

Precision Barraquer Style Colibri Forcep

Short Handle Design With 17mm x 5.5mm Finger Grip Cut-Out Very Delicate 1 x 2 (0.12mm) Teeth Without Tying Platforms Overall Length Of 2.9" (75.6mm)



VS11-6-2631

Precision Harms Style Colibri Forcep

Short Handle Design With 17mm x 5.5mm Finger Grip Cut-Out Very Delicate 1 x 2 (0.12mm) Teeth With Tying Platforms Overall Length Of 2.9" (75.6mm)

11-6-2633

Precision Dyson Style Colibri Forcep

Short Handle Design With 17mm x 5.5mm Finger Grip Cut-Out Very Delicate Short Tip Very Thin Shaft 1 x 2 (0.12mm) Teeth With Tying Platforms Overall Length Of 2.9" (75.6mm)

11-6-2661

Precision Barraquer Style Colibri Forcep

Short Handle Design With 17mm x 5.5mm Finger Grip Cut-Out (0.4mm) 1 x 2 Teeth With Tying Platforms Overall Length Of 2.9" (75.6mm)



VS5-4316

Precision Corneal Forcep

Colibri Style, Serrated And Wide Handle With Polished Finish, 1 x 2 (0.12mm) Teeth Set At 45 Degrees, 6mm Tying Platform, And Overall Length Of 4" (10.9cm). Overall Weight Of 0.4 oz

Colibri Forceps

Ophthalmology



VS11-6-2011

Castroviejo Colibri Forcep

0.12mm With 1 x 2 Teeth, On A Three Hole Handle Design, 1 x 2 (0.12mm) Teeth Set At 45 Degrees, Without Tying Platform, Overall Length Of 4 1/8" (105mm)

VS11-6-2013

Precision Harms Style Colibri Forcep

0.12mm With 1 x 2 Teeth, On A Three Hole Handle Design, 1 x 2 (0.12mm) Teeth Set At 45 Degrees, With Tying Platform, Overall Length Of 4 1/8" (105mm)



VS5-4314-12

University Of Kansas Corneal Forcep

Colibri Style, Serrated Handle With Polished Finish, 7mm Tying Platform, 2 x 2 (0.12mm) Teeth Set At 45 Degrees, And Overall Length Of 4 3/8" (11cm). Overall Weight Of 0.4 oz



VS5-4314-25

University Of Kansas Corneal Forcep

Colibri Style, Serrated Handle With Polished Finish, 7mm Tying Platform, 2 x 2 (0.25mm) Teeth Set At 45 Degrees, And Overall Length Of 4 1/4" (10.8cm). Overall Weight Of 0.4 oz

VS11-6-2031

Hyde Corneal Forcep

2x2 Teeth .12

VS5-4324

Precision Piarse Double Corneal Forcep

Colibri Style, Serrated Handle With Dull Finish, Two 3mm Long Tips Separated 1mm, With 0.2mm Piarse Tips, And Overall Length Of 2 3/4" (71mm)

VS5-4326

Polack Double Corneal Forcep

Colibri Style, Serrated Handle With Dull Finish, Two 2.75mm Long Tips Separated 1mm, 1 x 2 Teeth, And Overall Length Of 2 3/4" (71mm)

VS11-6-2671

Precision Polack Double Corneal Forcep

Colibri Style Short Handle Design With 17mm x 5.5mm Finger Grip Cut-Out With Two Tips Spaced 1mm Apart, 1 x 2 (0.12mm) Teeth, Without Tying Platforms Overall Length Of 2.9" (75.6mm)

VS11-6-2677

Precision Hoffman Polack Double Corneal Forcep

Colibri Style Long 3 Hole Handle Design With Two Tips Spaced 1mm Apart, (0.2mm) 1 x 1 Teeth, Without Tying Platforms Overall Length Of 4 1/8" (105mm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS5-4204

Girard Corneoscleral Forcep

Colibri Style, Round And Knurled Handle With Polished Finish, For Left-Handed Use, 6mm Tying Platform, 1 x 2 (0.12mm) Teeth, And Overall Length Of 4 1/4" (11cm). Overall Weight Of 0.7 oz

VS5-4205

Girard Corneoscleral Forcep

Colibri Style, Round And Knurled Handle With Polished Finish, For Left-Handed Use, 6mm Tying Platform, 1 x 2 (0.3mm) Teeth, And Overall Length Of 4 1/4" (107mm)



VS5-4206

Girard Corneoscleral Forcep

Colibri Style, For Right-Handed Use, With Round And Knurled Handle, Polished Finish, 6mm Tying Platform, 1 x 2 (0.12mm) Teeth, And An Overall Length Of 4 1/4" (11cm)

VS5-4207

O'Gawa Colibri Forcep

Round And Knurled Handle With Polished Finish, Alignment Pin, Colibri Style Tips, 1 x 2 (0.12mm) Teeth, And Overall Length Of 4 1/4" (106mm)



VS5-4208

Troutman-Barrquer Corneal Utility Forcep

Colibri Style, Serrated Handle With Polished Finish, 6mm Tying Platform, 1 x 2 (0.5mm) Teeth Set At 45 Degrees, And Overall Length Of 2 3/4" (7cm). Overall Weight Of 0.2 oz

VS5-2429-12

Catalano Corneal Scleral Forcep

Round Handle With Polished Finish, 1 x 2 Teeth, Cross Serrated 5mm 60 Degree Angled With Platforms, And Overall Length Of 4 1/4" (107mm)

VS5-2429-3

Catalano Corneal Scleral Forcep

Round And Knurled Handle With Dull Finish, 1 x 2 0.3mm Teeth, 5.5mm Tying Platform, And Overall Length Of 4 3/8" (113mm)

Colibri Forceps

Ophthalmology



VS5-4318-45

Maumenee Corneal Forcep

Colibri Style, Serrated Cross Action Handle With Polished Finish, 1 x 2 (0.1mm) Teeth Set At 45 Degrees, 5mm Tying Platform, And Overall Length Of 4 3/4" (121mm)



VS5-4318-90

Maumenee Corneal Forcep

Colibri Style, Serrated Cross Action Handle With Polished Finish, 1 x 2 (0.1mm) Teeth Set At 90 Degrees, 5mm Tying Platform, 21mm Mid Screw To Tip Length, And Overall Length Of 4 1/2" (11.5cm). Overall Weight Of 0.4 oz

VS11-6-4111

Maumenee Corneal Forcep

Colibri Style, Cross Action Squeeze Spring 3 Hole Handle Design Very Delicate (0.12mm) 1 x 2 Teeth Set At 90 Degrees, Without Tying Platforms Overall Length Of 4 1/8" (105mm)

VS5-2427

Kirby Corneoscleral Forcep

Round And Knurled Handle With Polished Finish, 1 mm x 2mm Teeth, Cross Serrated 5mm 60 Degree Angled Platform, And Overall Length Of 4 1/4" (107mm)

VS5-2577

Gill-Hess Iris Forcep

Serrated Handle With Polished Finish, Bayonet Tips, 1 x 2 Teeth, And Overall Length Of 2 3/4" (71mm)



VS5-2527

Bracken Curved Iris Forcep

Serrated Handle With Polished Finish, Delicate Cross-Serrated Lightly Curved Tips, And Overall Length Of 7/8" (98mm)



VS5-2535

Curved Iris Forcep

Serrated Handle With Polished Finish, Delicate Lightly Curved Tips, 1 x 2 Mouse Style Teeth, And Overall Length Of 3 7/8" (10cm). Overall Weight Of 0.2 oz

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020

VS5-2585

Gill Curved Iris Forcep

Serrated Handle With Polished Finish, Curved Shafts, Delicate 1 x 2 Teeth, And Overall Length Of 2 3/4" (72mm)

VS5-2586

Gill Curved Iris Forcep With Serrations

Serrated Handle With Polished Finish, Delicate Curved Shafts With Light Crisscross Serrations, And Overall Length Of 2 3/4" (72mm)

VS5-2588

Botvin Iris Forcep

Serrated Handle With Polished Finish, Bayonet Tips, Delicate Crisscross Serrations, And Overall Length Of 2 3/4" (71mm)

VS11-6-2151

Botvin Iris Forcep Bayonet Design

Short Colibri Style Handle With 17mm x 5.5mm Finger Grip Cut-Out, Bayonet Tips, Delicate Crisscross Serrations, And Overall Length Of 2.9" (75.6mm)

VS5-2587

Botvin Iris Forcep

Serrated Handle With Polished Finish, Bayonet Tips, Delicate 1 x 2 Teeth, And Overall Length Of 2 3/4" (71mm)

VS11-6-2161

Botvin Iris Forcep Bayonet Design

Short Colibri Style Handle With 17mm x 5.5mm Finger Grip Cut-Out, Bayonet Tips, Delicate 1 x 2 Teeth, Overall Length Of 2.9" (75.6mm)

VS5-2594

Graefe Straight Iris Forcep With Serrations

Serrated Handle With Polished Finish, Straight Tips With Delicate Serrations, And Overall Length Of 2 5/8" (67mm)



VS5-2595

Graefe Straight Iris Forcep With Teeth

Serrated Handle With Polished Finish, Straight Tips, Fine 1 x 2 Teeth, And Overall Length Of 2 7/8" (7cm)

VS5-2598

Graefe Curved Iris Forcep

Serrated Handle With Polished Finish, Curved Shafts, Delicate Crisscross Serrations, And Overall Length Of 2 5/8" (66mm)

Tying Forceps

Ophthalmology



VSTDK-3-611

Titanium Flat Handled Square Grip Tying Forcep
Straight 6mm Tying Platforms, Overall Length 3 1/3" (85mm)

VSTDK-3-612

Titanium Flat Handled Square Grip Tying Forcep
Curved 6mm Tying Platforms, Overall Length 3 1/4" (83mm)



VSTDK-3-615

Titanium Flat Handled Square Grip Tying Forcep
Straight 6mm Tying Platforms, Overall Length 4 1/2" (115mm)

VSTDK-3-616

Titanium Flat Handled Square Grip Tying Forcep
Curved 6mm Tying Platforms, Overall Length 4 1/2" (114mm)

VSTDK-3-616-1

Titanium Flat Handled Square Grip Tying Forcep
Curved Extra Delicate Design 6mm Tying Platforms, Overall Length 4 1/2" (114mm)

VSTDK-3-631

Titanium Flat Handled Square Grip Tying Forcep
Angled 45 Degees 5mm From Bend To Tip With A 4mm Tying Platforms, Overall Length 3 1/4" (82mm)

VSTDK-3-633

Titanium Flat Handled Square Grip Tying Forcep
Angled 45 Degees 7mm From Bend To Tip With A 6mm Tying Platforms, Overall Length 3 1/4" (82mm)



VSTDK-3-635

Titanium Flat Handled Square Grip Tying Forcep
Angled 45 Degees 10mm From Bend To Tip With A 6mm Tying Platforms, Overall Length 3 1/4" (82mm)

VSTDK-3-635-1

Titanium Flat Handled Square Grip Tying Forcep
Angled 45 Degees 10mm From Bend To Tip With A 10mm Tying Platforms, Overall Length 3 1/4" (82mm)

VSTDK-3-636

Titanium Flat Handled Square Grip Tying Forcep
Angled 45 Degees 5mm From Bend To Tip With A 4mm Tying Platforms, Overall Length 4 3/8" (112mm)

VSTDK-3-637

Titanium Flat Handled Square Grip Tying Forcep
Angled 45 Degees 12mm From Bend To Tip With A 6mm Tying Platforms, Overall Length 3 1/4" (82mm)

VSTDK-3-638

Titanium Flat Handled Square Grip Tying Forcep
Angled 45 Degees 7mm From Bend To Tip With A 6mm Tying Platforms, Overall Length 4 3/8" (112mm)



VSTDK-3-630

Titanium Flat Handled Square Grip Tying Forcep

Angled 45 Degees 10mm From Bend To Tip With A 6mm Tying Platforms, Overall Length 4 3/8" (112mm)



VS5-2807-WO

O'Gawa Castroviejo Tying Forcep

Serrated And Wide Handle With Polished Finish, Straight Shafts, 6mm Tying Platform, And Overall Length Of 4 1/4" (11cm). Overall Weight Of 0.4 oz



VS5-2926-S

McPherson Straight Tying Forcep

Serrated Handle with Dull Finish, Straight Shafts, 4mm Tying Platform, And Overall Length Of 3 1/2" (9cm). Overall Weight Of 0.3 oz



VS11-6-6111

McPherson Tying Forcep

On A Three Hole Handle Design With 18mm Straight Shafts With 5mm Tying Platforms, And An Overall Length Of 4 1/8" (10.5cm)



VS11-6-6111-T

Titanium McPherson Tying Forcep

On A Three Hole Handle Design With 18mm Straight Shafts With 5mm Tying Platforms, And An Overall Length Of 4 1/8" (10.5cm)



VS5-2926-A

McPherson Angled Tying Forcep

Serrated Handle With Dull Finish, Angled Shafts, 4mm Tying Platform, And Overall Length Of 3 1/4" (8.3cm). Overall Weight Of 0.2 oz



Tying Forceps

Ophthalmology



VS11-6-6121

McPherson Angled Tying Forcep

On A Three Hole Handle Design With 5mm Tying Platforms, And An Overall Length Of 4 1/8" (10.5cm)

VS11-6-6121-T

Titanium McPherson Angled Tying Forcep

On A Three Hole Handle Design With 5mm Tying Platforms, And An Overall Length Of 4 1/8" (10.5cm)



VS5-2926-AK

Kelman-McPherson Angled Tying Forcep

Serrated Handle With Dull Finish, 45 Degree Angled Shafts, 8mm From Tip To Bend, 5mm Tying Platform, And Overall Length Of 3 1/4" (8.3cm). Overall Weight Of 0.3 oz



VS11-6-6141


Kelman McPherson Angled Tying Forcep

On A Three Hole Handle Design With 7.5mm Tying Platforms, And An Overall Length Of 4 1/8" (10.5cm)

VS11-6-6141-T

Titanium Kelman McPherson Angled Tying Forcep

On A Three Hole Handle Design With 7.5mm Tying Platforms, And An Overall Length Of 4 1/8" (10.5cm)



VS11-6-6141-T-UF

Kelman McPherson Titanium Angled Tying Forcep

Ultra Fine With Three Holed Handle, 8mm From Tip To Bend, And An Overall Length Of 4" (10cm)



VS5-2926-AB

Bechert-McPherson Angled Tying Forcep

Serrated Handle With Dull Finish, 45 Degree Angled Shafts, 10mm Tying Platforms, And Overall Length Of 3 1/4" (8.3cm). Overall Weight Of 0.2 oz

VS11-6-6150

Bechert McPherson Angled Tying Forcep

On A Three Hole Handle Design With 10mm Tying Platforms, And An Overall Length Of 4 1/8" (10.5cm)

VS11-6-6150-T

Titanium Bechert McPherson Angled Tying Forcep

On A Three Hole Handle Design With 10mm Tying Platforms, And An Overall Length Of 4 1/8" (10.5cm)

VS5-2920

Clayman-McPherson Angled Tying Forcep

Smooth Handle With Dull Finish, 45 Degree Angled Shafts, 8mm From Bend To Tip, 4mm Tying Platforms, And Overall Length Of 3 1/4" (80mm)



VS5-2918

Sheets-McPherson Angled Tying Forcep

Serrated Handle With Dull Finish, 45 Degree Angled Shafts, 10.5mm Tying Platform, 1mm Of Tip Is Crisscross Serrated, And Overall Length Of 3 1/2" (9cm). Overall Weight Of 0.2 oz



VS5-2846

MacGregor Conjunctiva Forcep

Serrated Handle With Polished Finish, Straight Shafts, 5.5mm Tying Platform, 1mm Crisscross Serrated Tips, And Overall Length Of 4" (10cm). Overall Weight Of 0.4 oz



VS5-4203-TS

Girard Tying Forcep

Round And Knurled Handle With Dull Finish, Straight Shafts, 5mm Tying Platform, And Overall Length Of 4" (10.2cm). Overall Weight Of 0.4 oz



VS5-4203-T

Girard Tying Forcep

Round And Knurled Handle With Dull Finish, Curved Shafts, 5mm Tying Platform, And Overall Length Of 3 7/8" (99mm). Overall Weight Of 0.4 oz

Tying Forceps

Ophthalmology



VS5-4203-FT

Girard Fine Tying Forcep

Round And Knurled Handle, Curved Shafts, Ultra-Fine Tying Platform, And Overall Length Of 3 7/8" (100mm)



VS5-4203-FTS

Girard Fine Tying Forcep


Round And Knurled Handle, Straight Shafts, Ultra-Fine Tying Platform, And Overall Length Of 4" (10.2cm). Overall Weight Of 0.4 oz



VS11-6-6331

Tennant Tying Forcep

With Straight 19mm Shafts, With Longitudinal Serrated Handle Design, And An Overall Length Of 4" (10cm)



VS11-6-6341

Tennant Tying Forcep

With Curved 19mm Shafts, With Longitudinal Serrated Handle Design, And An Overall Length Of 4" (10cm)



VS5-2992

O'Gawa Curved Tying Forcep

Round And Knurled Handle , With Polished Finish, Curved Shafts, Tips, 5.5mm Tying Platform, 0.3mm Gap At Platform Heel, And Overall Length Of 4 1/8" (10.5cm). Overall Weight Of 0.4 oz



VS5-2992-S

O'Gawa Straight Tying Forcep

Round And Knurled Handle , Straight Shafts, Tips, 5.5mm Tying Platform, 0.3mm Gap At Platform Heel, And Overall Length Of 4 1/4" (108mm)

VS5-2967

Catalano Straight Micro Tying Forcep

Round And Knurled Handle With Polished Finish, Straight Shafts, 6mm Tying Platform, Crisscross Serrated Tips, And Overall Length Of 4 1/2" (114mm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020

VS5-2967-C

Catalano Angled Micro Tying Forcep

Round And Knurled Handle With Dull Finish, 45 Degree Angled Shafts, 6mm Tying Platform, Crisscross Serrated Tips, And Overall Length Of 4 3/8" (111mm)



VS5-2999-S

Harms-Tubingen Straight Tying Forcep

Serrated Handle With Polished Finish, Straight Shafts, 6mm Tying Platform, And Overall Length Of 4" (10cm). Overall Weight Of 0.4 oz



VS5-2999-C

Harms-Tubingen Curved Tying Forcep

Serrated Handle With Polished Finish, Curved Shafts And Tips, 7mm Tying Platform, And Overall Length Of 4" (10.1cm). Overall Weight Of 0.4 oz



VS5-2919

Osher Tying Forcep

Serrated And Flat Handle With Dull Finish, Curved Shafts, Tips, 6 1/2 Tying Platform, And Overall Length Of 4" (10cm). Overall Weight Of 0.4 oz



VS5-2921

Osher Ultra Fine Tying Forcep

Serrated And Flat Handle With Dull Finish, Curved Shafts, Tips, 6.5mm Tying Platform, And Overall Length Of 4" (10cm). Overall Weight Of 0.4 oz



VS11-6-6511

Harms Straight Tying Forcep

On A Three Hole Handle Design, Straight 23mm Shafts, Overall Length Of 4 1/4" (11 cm)

Tying Forceps


Ophthalmology



VS11-6-6521

Harms Curved Tying Forcep

On A Three Hole Handle Design, Curved 20mm Shafts, An Overall Length Of 4" (10cm)



VS11-6-6311

Jaffe Straight Tying Forcep

On A Three Hole Handle Design, Straight 23mm Shafts, Overall Length Of 4" (10cm)



VS11-6-6321

Jaffe Curved Tying Forcep

On A Three Hole Handle Design Straight 23mm Shafts, Overall Length Of 4" (10cm)



VS11-6-6201

Freeman Jaffe Curved Lens Holding Forcep

On A Three Hole Handle Design With Longitudinal Groove In Jaw, Overall Length Of 4" (10cm)



VS5-2966

Shepard Straight Tying Forcep

Serrated And Wide Handle With Dull Finish, Straight Shafts, 4mm Tying Platform, And Overall Length Of 4 1/4" (11mm). Overall Weight Of 0.4 oz





VS5-2966-A

Shepard Angled Tying Forcep

Serrated And Wide Handle With Dull Finish, 45 Degree Angled Shafts, Tying Platform, And Overall Length Of 4 1/4" (11 cm). Overall Weight Of 0.4 oz



VS5-2966-C

Shepard Curved Tying Forcep

Serrated And Wide Handle With Dull Finish, Curved Shafts, Tips, 6mm Tying Platform, And Overall Length Of 4 1/4" (11 cm). Overall Weight Of 0.4 oz



VS5-2908

Precision Straight Tying Forcep

Wide Handle With Dull Finish, Straight Shafts, 5mm Tying Platform, And Overall Length Of 4 3/8" (11 cm). Overall Weight Of 0.4 oz



VS5-2936

Sinskey Straight Tying Forcep

Smooth Flat Cross Action Handle With Dull Finish, Straight Fine 7mm Tying Platform, And Overall Length Of 4 3/8" (111 mm)



VS5-2909-CM

Anis Curved Tying Forcep

Round And Knurled Handle With Dull Finish, Curved Delicate 5mm Tying Platforms, And Overall Length Of 4 5/8" (118 mm)

VS11-6-6223

Anis Tying Forcep

Straight, On A Round Spring Handle Straight Smooth Jaw

VS11-6-6224

Anis Tying Forcep

Curved, On A Round Spring Handle Straight Smooth Jaw

Tying Forceps

Ophthalmology



VS5-2996

Fine Tying Forcep (also known as Max Fine)

Serrated And Wide Handle With Dull Finish, Curved Shafts, Tips, Tying Surfaces, And Overall Length Of 4 1/4" (11cm). Overall Weight Of 0.5 oz

VS11-6-6661

Precision Fine Style Curved Tying Forcep

On A Three Hole Handle Design Overall Length Of 4 1/4" (11cm)

VS5-2907

Precision Tying Forceps

Round And Knurled Handle With Dull Finish, Straight Shafts, 4mm V-Groove Tying Platform, And Overall Length Of 4 1/4" (108mm)

VS11-6-6345

Precision Maumenee Style Straight Tying Forcep

On A Round Longitudinal Serrated Handle Design, Inside Jaw Has V Shaped Groove Overall Length Of 4" (10cm)

VS11-6-6346

Precision Maumenee Style Angled Tying Forcep

On A Round Longitudinal Serrated Handle Design, Inside Jaw Has V Shaped Groove Overall Length Of 4" (10cm)

VS5-2998

Precision Tying Forcep

Serrated And Wide Handle With Polished Finish, Straight And Heavy Shaft, 7mm Tying Surfaces, And Overall Length Of 4 1/4" (106mm)



VS5-2995

Puntteney Curved Tying Forcep

Serrated Handle With Polished Finish, Curved Shafts, 12mm Tying Platform, And Overall Length Of 4" (10.2cm)

VS5-2993

Troutman Tying Forcep

Serrated Cross Action Handle With Polished Finish, Slightly Curved 10mm Jaws, And Overall Length Of 4 5/8" (118mm)

VS11-6-6167

Kurstin Flap Stretching Forcep

On A 3 Hole Handle



VS11-6-6169

Buratto Lasik Flap Forcep

On A 3 hole Handle With Serrated Disk Shaped Tips With Curved Shaft

VS11-6-6171

Mendez Lasik Forcep

On A 3 hole Handle With Spatulated Jaw, Smooth Tips Vaulted Shaft

5-2053-15

Pierse Type Micro Forcep

Smooth And Flat Handle With Dull Finish, Straight Shafts, Fine Pierse Type Tips, And Overall Length Of 3 3/4" (94mm)

5-2053-16

Pierse Type Micro Forcep

Serrated And Flat Open Handle With Dull Finish, Curved Shafts, Micro Pierse Type Tips, And Overall Length Of 3 3/8" (86mm)

5-2053-18

Pierse Type Fine Skeleton Forcep

Serrated And Flat Open Handle With Dull Finish, Straight Shafts, Fine Pierse Type Tips, And Overall Length Of 3 1/4" (82mm)

5-2053-19

Pierse Type Micro Skeleton Forcep

Serrated And Flat Open Handle With Dull Finish, Straight Shafts, Micro Pierse Type Tips, And Overall Length Of 3 1/4" (82mm)

5-2053-22

Pierse Type Colibri Forcep

Serrated And Flat Open Handle With Dull Finish, Micro Pierse Type Tips, 2mm Angled Colibri Style Jaws, And Overall Length Of 3 1/8" (80mm)

5-2053-30

Pierse Type Colibri Forcep

Serrated And Flat Open Handle With Dull Finish, Mini-Micro Pierse Type Tips, 2mm Angled Colibri Style Jaws, And Overall Length Of 3 1/8" (80mm)

5-2054-15

Pierse Micro Forcep

Straight 0.1mm Pierse Tips Overall Length 3 1/4" (82mm)

5-2054-16

Pierse Colibri Style Micro Forcep

Straight 0.1mm Pierse Tips Overall Length 3 1/4" (82mm)

5-2054-21

Pierse Type Colibri Forcep

Flat And Open Handle With Dull Finish, Colibri Style 0.3mm Pierse Tips, And Overall Length Of 3 1/8" (78mm)

5-2054-22

Pierse Type Colibri Micro Forcep

Flat And Open Handle With Dull Finish, Colibri Style 0.1mm Pierse Tips, And Overall Length Of 3 1/8" (78mm)

Nucleus / Capsule Forceps

Ophthalmology



VS16-3446-R

Titanium Salz Nucleus Splitter

Short Colibri Style Handle Squeeze Cross Action Design For Paddles To Be Spread Apart



VSTDK-3-911

Titanium B Nucleus Cracker

Cross Action Spring Squeeze Square Grip Handle With Two Paddle Tips

VS5-3188

Osher Capsule-Foreign Body Forcep

Serrated Handle With Dull Finish, Fine 3mm Platform Angled 60 Degrees, 10mm From Angled To Tip, And Overall Length Of 3 1/8" (80mm)



VS5-4200-GW

Gills-Welsh Capsule Forcep

Colibri Style, Serrated Handle With Polished Finish, Smooth Colibri Tip, And Overall Length Of 2 7/8" (73mm)



5-3139

Castroviejo Curved Capsule Forcep

Serrated Cross Action Handle With Polished Finish, 50 Degree Angled 2.5mm Cups, 10mm From Angle To Tip, And Overall Length Of 4 1/8" (104mm)



VS5-3134

Arruga Curved Capsule Forcep

Serrated Handle With Polished Finish, 50 Degree Angled 2.5mm Cups, 10mm From Angled To Tip, And Overall Length Of 4 1/4" (106mm)





VS5-2511

Straight Dressing Forcep With Serrations

Serrated Handle With Polished Finish, Delicate 11mm Straight Tips, And Overall Length Of 4" (10cm). Overall Length Of 0.3 oz



VS5-2515

Straight Dressing Forcep With Serrations

Serrated Handle With Polished Finish, Heavy 12mm Straight Tips, And Overall Length Of 3 7/8" (9.9cm). Overall Length Of 0.2 oz



VS5-2517

Curved Dressing Forcep With Serrations

Serrated Handle With Polished Finish, Gently Curved Delicate Tips, And Overall Length Of 3 7/8" (9.8cm) . Overall Length Of 0.2 oz



VS5-2519

Curved Dressing Forcep With Serrations

Serrated Handle With Polished Finish, Delicate Strong Curved Serrated Tips, And Overall Length Of 3 3/4" (97mm)



VS5-2521

Curved Dressing Forcep With Serrations

Serrated Handle With Polished Finish, Gently Curved Heavy Tips, And Overall Length Of 3 7/8" (10cm). Overall Weight Of 0.3 oz



VS5-2617

Bishop-Harmon Dressing Forcep

Flat Handle With Holes And Polished Finish, Straight Crisscross Serrated Tips, And Overall Length Of 3 3/8" (9cm). Overall Weight Of 0.2 oz



Dressing Forceps

Ophthalmology



VS15-11120

Adson Dressing Forcep
4 3/4" Delicate



VS15-11150

Adson Dressing Forcep
6" Delicate



VS15-12120

Adson Dressing Forcep
4 3/4" Standard



VS15-12150

Adson Dressing Forcep
6" Standard



VS15-13120

Adson Dressing Forcep
4 3/4" Standard Cross Serrated Jaw

VS15-14120

Adson Ewald Dressing Forcep
4 3/4" Heavy

VS15-15120

Micro Adson Dressing Forcep
4 3/4" Very Delicate



VS15-15150

Micro Adson Dressing Forcep
6" Very Delicate

VS15-16120

Adson Dressing Forcep
4 3/4" Narrow

VS86-12120

Adson TC Dressing Forcep
4 3/4" Standard With Tungsten Carbide Inserts

VS86-12150

Adson TC Dressing Forcep
6" Standard With Tungsten Carbide Inserts



VS5-2611

Bishop-Harmon Straight Tissue Forcep
Flat Handle With Holes And Polished Finish, Straight Shafts, Fine 1 x 2 Teeth, 0.8mm Tip, With An Overall Length Of 3 3/8" (8.5cm). Overall Weight Of 0.2 oz



VS5-2613

Bishop-Harmon Straight Tissue Forcep
Flat Handle With Holes And Polished Finish, Straight Shafts, Delicate 1 x 2 Teeth, 0.3mm Tip, And Overall Length Of 3 1/2" (9cm). Overall Weight Of 0.2 oz



VS5-2611-C

Bishop-Harmon Curved Tissue Forcep
Flat Handle With Holes And Polished Finish, Curved Shafts, Fine 1 x 2 Teeth, And Overall Length of 3 3/8" (8.5cm) 0.8mm Tip. Overall Weight Of 0.2 oz

Tissue Forceps

Ophthalmology



VS5-2617

Bishop-Harmon Dressing Forcep

Flat Handle With Holes And Polished Finish, Straight Crisscross Serrated Tips, And Overall Length Of 3 3/8" (9cm). Overall Weight Of 0.2 oz



VS5-2531

Straight Tissue Forcep

Serrated Handle With Polished Finish, Straight Shafts, 1 x 2 Delicate Teeth, And Overall Length Of 3 7/8" (10cm). Overall Weight Of 0.2 oz



VS5-2533

Tissue Forcep

Heavy 1 x 2 Teeth With An Overall Length Of 4 1/8" (10.4cm). Overall Weight Of 0.4 oz

VS21-14120

Adson Tissue Forcep

1 x 2 Teeth 4 3/4" Delicate

VS21-14150

Adson Tissue Forcep

1 x 2 Teeth 6" Delicate



VS21-15120

Adson Tissue Forcep

1x2 Teeth Standard 4 3/4" Overall Length

VS21-15150

Adson Tissue Forcep

1x2 Teeth Standard 6" Overall Length



VS21-16120

Adson Tissue Forcep

Serrated With 1 x 2 Teeth And Overall Length Of 4 3/4" (12.1cm)



VS21-17120

Adson Ewald Tissue Forcep

1 x 2 Teeth 4 3/4" Heavy



VS21-18120

Adson Ewald Tissue Forcep

2 x 3 Teeth 4 3/4" Heavy



VS21-18150

Adson Ewald Tissue Forcep

2 x 3 Teeth 6" Heavy



VS22-10120

Micro Adson Tissue Forcep

1 x 2 Teeth 4 3/4" Delicate



VS22-10150

Micro Adson Tissue Forcep

1 x 2 Teeth 6" Delicate

22-11120

Adson Tissue Forcep

1 x 2 Teeth 4 3/4" Narrow



VS22-13120

Adson Brown Tissue Forcep

With 7 x 7 Teeth Overall Length Of 4 3/4" (12.1cm)



Tissue/Jewelers Forceps

Ophthalmology



22-14120

Adson Brown Tissue Forcep

With 7 x 7 Teeth Angled Overall Length Of 4 3/4" (12.1cm)



VS22-15120

Adson Brown Tissue Forcep

With 9 x 9 Teeth Overall Length Of 4 3/4"



VS22-16120

Adson Brown Tissue Forcep

With 9 x 9 Teeth Angled Overall Length Of 4 3/4" (12.1cm)

VS5-2058-0

Swiss Fine Pointed Jewelers Forcep #0

With An Overall Length Of 4 5/8" (11.7cm)



VS5-2058-1

Swiss Jeweler Forcep #1

Smooth Handle With Dull Finish, Straight Shafts, Fine Pointed Tips, And Overall Length Of 4 3/4" (12cm). Overall Weight Of 0.4 oz



VS5-2058-2

Swiss Jeweler Forceps #2

Smooth Handle With Dull Finish, Straight Shafts, Fine Pointed Tips (Thickness 0.14mm (0.0055 Inch) By (Width 0.34mm (0.0134) Overall Length Of 4.728" (120mm)



VS5-2058-3

Swiss Jeweler Forcep #3

Smooth Handle With Dull Finish, Straight Shafts, Extra Fine Pointed Tips (Thickness 0.10mm (0.0039 Inch) By (Width 0.17mm (0.0067) Overall Length Of 4.728" (120mm)



VS5-2058-3C

Swiss Jeweler Forcep #3C

Smooth Handle With Dull Finish, Straight Shafts, Extra Fine Pointed Tips (Thickness 0.10mm (0.0039 Inch) By (Width 0.17mm (0.0067) Overall Length Of 4.334" (110mm)



VS5-2058-4

Swiss Jeweler Forcep #4

Smooth Handle With Dull Finish Straight Shafts Super Fine Pointed Tips (Thickness 0.08mm (0.0032 Inch) By (Width 0.13mm (0.0051) Overall Length Of 4.334" (110mm)



VS5-2058-5

Swiss Jeweler Forcep #5

Smooth Handle With Dull Finish Straight Narrowed Shafts Super Fine Pointed Tips (Thickness 0.06mm (0.0024 Inch) By (Width 0.10mm (0.0039) Overall Length Of 4.334" (110mm)



VS5-2058-5H

Swiss Hirschman Jeweler Forcep #5H

Smooth Handle With Dull Finish, Straight And Narrowed Shafts, Extra Fine Blunted Tips, And Overall Length Of 4 1/4" (110mm)



VS5-2058-7

Swiss Jeweler Forcep #7

Curved Smooth Handle With Dull Finish Curved Shafts Extra Fine Pointed Tips (Thickness 0.10mm (0.0039 Inch) By (Width 0.17mm (0.0067) Overall Length Of 4.531" (115mm)



Jewelers Forceps

Ophthalmology



VS11-6-7611

Utility And Splinter Forcep

With Sharp Tips And An Overall Length Of 3 1/2" (9cm)



VS11-6-7621

German Jeweler Fine Pointed Forcep #1

Smooth And Flat Handle, And An Overall Length Of 4 3/4" (12cm)

VS11-6-7631

German Jeweler Fine Pointed Forcep #3C

Smooth And Flat Handle, And An Overall Length Of 4 3/4" (12cm)



VS11-6-7641

German Jeweler Forcep #3

With Extra-Fine Pointed Tips, Smooth And Flat Handle, And An Overall Length Of 4 3/4" (12cm)



VS11-6-7651

German Jeweler Forcep #4

With Smooth Handle, Fine Points, Straight Shafts, And An Overall Length Of 4 1/4" (11cm)



VS11-6-7661

German Jewelers Forcep #5

With Extra-Fine Pointed Tips, Straight 28mm Shafts, Smooth And Flat Handle, And An Overall Length Of 4 1/4" (11cm)



11-6-7681

German Jewelers Forcep #7

Curved Fine Pointed Tips Smooth And Flat Handle, And An Overall Length Of 4 1/4" (11cm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS11-6-7686

Suture Removal Forcep

With Very Fine Pointed Tips, Carbide Impregnated Jaws, And An Overall Length Of 4 1/2" (11.5cm)



VS11-6-7711

Jaffe Utility Forcep

With Curved Shafts, Three-Holed Handle, And An Overall Length Of 4 1/4" (11 cm)



VS5-3691

Wies Chalazion Forcep

Serrated Handle With Polished Finish, Solid Lower Plate With Outside Width Of 15mm, Open Serrated Upper Plate With Inside Dimensions Of 9mm x 12mm, Locking Thumb Screw, And Overall Length Of 3 3/4" (9.5cm). Overall Weight Of 0.5 oz



VS5-3711

Desmarres Chalazion Forcep

Serrated Cross Action Handle With Polished Finish, Oval Solid 16mm Wide Lower Plate, Open Upper Plate With Inside Dimensions Of 12mm x 14mm, And Overall Length Of 3 3/4" (9.5cm). Overall Weight Of 0.4 oz



VS5-3715

Desmarres Chalazion Forcep

Serrated Handle With Polished Finish, Locking Thumb Screw, Oval Solid 24mm Wide Lower Plate, Open Upper Plate With Inside Dimensions Of 12mm x 20mm, And Overall Length Of 3 1/2" (9cm). Overall Weight Of 0.5 oz



Chalazion Forceps

Ophthalmology



VS5-3721

Francis Chalazion Forcep

Serrated Handle With Polished Finish, Locking Thumb Screw, Oval Solid 17mm Wide Lower Plate, Open Upper Plate With Inside Dimensions Of 12mm x 14mm, And Overall Length Of 3 3/4" (97mm)



VS5-3725

Heath Chalazion Forcep

Serrated Handle With Polished Finish, Locking Thumb Screw, Oval Solid 16mm Wide Lower Plate, Open Upper Plate With Inside Dimensions Of 12mm x 14mm, And Overall Length Of 3 3/4" (96mm)



VS5-3729

Hunt Chalazion Forcep

Serrated Handle With Polished Finish, Locking Thumb Screw, Round Solid 15mm Wide Lower Plate, Open Upper Plate With Inside Diameter Of 12mm, And Overall Length Of 3 3/4" (97mm)



VS5-3741

Lambert Chalazion Forcep

Small, Serrated Handle With Polished Finish, Locking Thumb Screw, Round Solid 10mm Wide Lower Plate, Open Upper Plate With Inside Diameter Of 8mm, And Overall Length Of 3 5/8" (92mm)





VS5-3743

Lambert Chalazion Forcep

Large, Serrated Handle With Polished Finish, Locking Thumb Screw, Round Solid Lower Plate, Open Upper Plate With Inside Diameter Of 15mm, And Overall Length Of 3 3/4" (97mm)



VS5-3751

Baird Chalazion Forcep

Serrated Handle With Polished Finish, Locking Thumb Screw, Oval Solid 15mm Wide Lower Plate, Open Upper Plate With Inside Dimensions Of 11mm x 8mm, And Overall Length Of 3 3/4" (97mm)



VS5-3625

Ehrhardt Lid Forcep

Serrated Handle With Polished Finish, 14mm Oval Solid Lower Plate, Straight 9mm Serrated Upper Jaw, Locking Thumb Screw, And Overall Length Of 3 1/2" (9cm). Overall Weight Of 0.4 oz



VS5-3633

Snellen Entropion Forcep

Left, Serrated Handle With Polished Finish, Solid Lower Plate With Outside Width Of 28mm, Open Upper Plate, Locking Thumb Screw, And Overall Length Of 3 1/2" (88mm)

VS5-3634

Snellen Entropion Forcep

Left, Serrated Handle With Polished Finish, Solid Lower Plate With Outside Width Of 28mm, Open Upper Plate, Locking Thumb Screw, And Overall Length Of 3 1/2" (88mm)

Muscle Forceps

Ophthalmology



VS5-2797

Troutman Superior Rectus Forcep

Serrated Handle With Polished Finish, 45 Degree Angled Shafts, 10mm From Tip To Bend, 1 x 2 Teeth, And Overall Length Of 4 1/4" (106mm)



VS5-2798

Osher Superior Rectus Forcep

Serrated Handle With Polished Finish, 45 Degrees Angled Shafts, 10mm From Tip To Bend, Tips Curved With 1 x 2 (0.5mm) Teeth, And Overall Length Of 4 1/8" (10.5cm). Overall Weight Of 0.4 oz



VS5-2703

Rizzuti Superior Rectus Forcep

Serrated Cross Action Handle With Polished Finish, 45 Degree Angled Shafts, 7mm From Tip To Bend, Tips Curved With 1 x 2 Teeth, And Overall Length Of 3 3/4" (9.6cm). Overall Weight Of 0.3 oz



VS5-3415

Bishop Tendon Tucker

Round And Knurled Handle With Polished Finish, Triangular Shaped Tips, Marked For 5mm And 10mm Tucks, And Overall Length Of 4 1/4" (10.8cm)



VS5-3441

Jameson Muscle Forcep—Adult Left

Left, Serrated Handle With Slide Lock And Polished Finish, Angled Shafts, 12mm Jaws From Bend To Tip, 1mm Teeth On One Jaw Fit Into Holes On Other Jaw, And Overall Length Of 4" (10cm). Overall Weight Of 0.4 oz



VS5-3442

Jameson Muscle Forcep—Adult Right

Right, Serrated Handle With Slide Lock And Polished Finish, Angled Shafts, 12mm Scaled Jaws From Bend To Tip, 1mm Teeth On One Jaw Fit Into Holes On Other Jaw, And Overall Length Of 4" (10cm). Overall Weight Of 0.4 oz

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS5-3445

Jameson Muscle Forcep-Child Left

Child Size, Left, Serrated Handle With Slide Lock And Polished Finish, Angled Shafts, 9mm Jaw From Bend To Tip, 1mm Teeth On One Jaw Fit Into Holes On Other Jaw, And Overall Length Of 3 7/8" (99mm)



VS5-3446

Jameson Muscle Forcep-Child Right

Child Size, Right, Serrated Handle With Slide Lock And Polished Finish, Angled Shafts, 9mm Jaw From Bend To Tip, 1mm Teeth On One Jaw Fit Into Holes On Other Jaw, And Overall Length Of 3 7/8" (99mm)

VS5-3451

Gunderson Muscle Forcep

Serrated Handle With Slide Lock And Polished Finish, Angled Shafts, 11mm Long Tips, Jaws With Longitudinal Meeting Serrations, And Overall Length Of 3 7/8" (100mm)

VS5-3461

Berens Muscle Clamp

Round And Knurled Handle With Polished Finish, 13.5mm Angled Adjustable Blades, Teeth In Lower Jaw Fit Into Holes In Upper Jaw, Grooves In Lower Jaw, And Overall Length f 3 1/2" (89mm)



VS5-3611

Berke Ptosis Forcep

Serrated Handle With Slide Lock And Polished Finish, Multi-Curved Shafts, Straight 20mm Blades, Longitudinal Serrations, And Overall Length Of 4 1/4" (11cm). Overall Weight Of 0.6oz



VS5-3615

Berke Ptosis Forcep

Serrated Handle With Slide Lock And Polished Finish, Multi-Curved Shafts, Straight 27mm Blades, Longitudinal Serrations, And Overall Length Of 4" (10cm). Overall Weight Of 0.6 oz

VS5-3618

Iliff Blepharochalasis Forcep

Serrated Handle With Polished Finish, Slide Lock, Multi-Curved Shafts, Curved 28mm Serrated Jaws, And Overall Length Of 4 1/4" (107mm)

Cilia Forceps

Ophthalmology



VS5-3267

Precision Cilia Forcep

Serrated Handle With Polished Finish, Wide 1 mm x 3mm Smooth Platform, And Overall Length Of 3" (76mm)



VS5-3271

Precision Cilia Forcep

Serrated Handle With Polished Finish, Narrow Smooth Jaws, And Overall Length Of 3 1/4" (8.2cm). Overall Weight Of 0.2 oz



VS5-3273

Barraquer Cilia Forcep

Serrated Handle With Polished Finish, Slightly Curved 6mm Long Platform, And Overall Length Of 4" (10.2cm) . Overall Weight Of 0.3 oz



VS5-3275

Ziegler Cilia Forcep

Serrated Handle With Polished Finish, Slender 4.5mm Long Rounded Platform, And Overall Length Of 3 1/2" (9cm). Overall Weight Of 0.3 oz

VS5-3277

Drews Cilia Forcep

Serrated Handle With Polished Finish, Narrow 4mm Long Smooth Platform, And Overall Length Of 4 1/8" (104mm)

VS5-3281

Gradle Cilia Forcep

Serrated Cross Action Handle With Metal Tip Guard And Polished Finish, Fine 2mm Wide Platform, And Overall Length Of 3 3/4" (97mm)



VS50-10090

Douglas Cilia Forcep

Narrow Longer Serrated Tip Overall Length 3 1/2" (89mm)

VS50-11090

Douglas Cilia Forcep

Wide Short Serrated Tip Overall Length 3 1/2" (89mm)



VS51-10090

Ziegler Cilia Forcep

Rounded Pointed Tips With Smooth Platforms Overall Length 3 1/2" (89mm)



VS52-10080

Beer Cilia Forcep

Rounded Blunt Tips With Smooth Platforms Overall Length 3 1/6" (80mm)

VS52-10090

Beer Cilia Forcep

Rounded Blunt Tips With Smooth Platforms Overall Length 3 1/2" (89mm)



VS54-10090

Bergh Cilia Forcep

Angled Pointed Tip Serrated Jaw Beveled On The Outside Of Jaw Overall Length 3 1/2" (89mm)



VS54-12080

Littauer Cilia Forcep

Bowed Out Upside Down Triangular Tip Serrated Jaw Beveled On The Outside Of Jaw Overall Length 3 1/2" (89mm)



VS54-15105

Barraquer Cilia Forcep,

Serrated Handle Slightly Angled 6mm Smooth Platforms, Overall Length Of 4 1/2" (11.5cm)

VS55-10105

Bergh Cilia Forcep

Angled Pointed Tip Serrated Jaw Beveled On The Outside Of Jaw Overall Length 4 1/2" (114mm)

Muscle Clamp Forceps

Ophthalmology



VSP7201I

MSI Precision "Mueller Putterman Style 3 Pins" Ptosis Muscle Clamp For Fasanella Procedure
Overall Length 3 3/4" (95mm) 3 Pins On Top Jaw With 3 Aligned Holes On Bottom Jaw, Jaw Has 4 Inside Interlocking Grooves, Jaw Head Is Concaved And 22mm Wide



VSP7201T

MSI Precision Titanium "Mueller Putterman Style 3 Pins" Ptosis Muscle Clamp For Fasanella Procedure
Overall Length 3 3/4" (95mm) 3 Pins On Top Jaw With 3 Aligned Holes On Bottom Jaw, Jaw Has 4 Inside Interlocking Grooves, Jaw Head Is Concaved And 22mm Wide



VSP7202I

MSI Precision "Mueller Putterman Style 6 Pins" Muscle Conjunctiva Resection Blepharoptosis Clamp
Overall Length 3 3/4" (95mm) 3 Pins On Top And 3 Pins On Bottom Jaw With 3 Aligning Holes On Both Jaws, Jaw Has 4 Inside Interlocking Grooves, Jaw Head Is Concaved And 24mm Wide



VSP7202T

MSI Precision Titanium "Mueller Putterman Style 6 Pins" Muscle Conjunctiva Resection Blepharoptosis Clamp
Overall Length 3 3/4" (95mm) 3 Pins On Top And 3 Pins On Bottom Jaw With 3 Aligning Holes On Both Jaws, Jaw Has 4 Inside Interlocking Grooves, Jaw Head Is Concaved And 24mm Wide

VS5-2986

Collis Micro Utility Forcep

Serrated And Wide Handle With Polished Finish, Straight Shafts, 3mm Platform, And Overall Length Of 4 1/4" (109mm)



VS5-2901

Nugent Utility Forcep

Serrated Handle With Polished Finish, 45 Degree Angled Shafts, 10mm Smooth Jaws, And Overall Length Of 4 1/4" (11cm). Overall Weight Of 0.3 oz



VS5-2903

Nugent Utility Forcep


Serrated Handle With Polished Finish, 45 Degree Angled Shafts, 10mm Serrated Jaws, And Overall Length Of 4 1/4" (11cm). Overall Weight Of 0.3 oz



VS5-2912

Thorton Angled Fixation Forcep

With Wide Serrated Handle With Polished Finish, Straight Shafts With 25 Degrees Angled 1.5mm Tips, 4mm Platform With (0.12mm)1 x 2 Teeth, And An Overall Length Of 4.3" (109mm)



VS5-2917

Castroviejo Suturing Forcep

Polished Finish, 45 Degree Angled Shafts, 0.5mm Teeth, 6mm Tying Platform, And Overall Length Of 4 1/4" (108mm)



VS5-2919

Osher Tying Forcep

Serrated And Flat Handle With Dull Finish, Curved Shafts, Tips, 6 1/2 Tying Platform, And Overall Length Of 4" (10cm). Overall Weight Of 0.4 oz



VS5-2925

O'Gawa Angled Tying Forcep

Round And Knurled Handle Bottom of Platform With Polished Finish, Angled Shafts, 10mm Tying Platform, And Overall Length Of 4 1/4" (106mm)



VS5-2029

Osher Tying Forcep

Serrated And Flat Handle With Dull Finish, Curved Shafts, Tips, 6 1/2 Tying Platform, And Overall Length Of 4" (10cm). Overall Weight Of 0.4 oz

Hemostat Forceps

Ophthalmology



VS5-4026

Hartman Straight Hemostatic Mosquito Forcep

Ring Handle With Dull Finish, Ratchet Lock, Straight 20mm Serrated Jaws, And Overall Length Of 3 7/8" (99mm)



VS5-4027

Hartman Curved Hemostatic Mosquito Forcep

Ring Handle With Dull Finish, Ratchet Lock, Curved 20mm Serrated Jaws, And Overall Length Of 3 7/8" (99mm)



VS5-4031

Halstead Straight Hemostatic Mosquito Forcep

Ring Handle With Dull Finish, Ratchet Lock, Straight 22mm Serrated Jaws, And Overall Length Of 4 7/8" (125mm)



VS5-4031

Halstead Straight Hemostatic Mosquito Forcep

Ring Handle With Dull Finish, Ratchet Lock, Straight 22mm Serrated Jaws, And Overall Length Of 4 7/8" (125mm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020

VS5-2585

Gill Curved Iris Forcep

Serrated Handle With Polished Finish, Curved Shafts, Delicate 1 x 2 Teeth, And Overall Length Of 2 3/4" (72mm),

VS5-2586

Gill Curved Iris Forcep With Serrations

Serrated Handle With Polished Finish, Delicate Curved Shafts With Light Crisscross Serrations, And Overall Length Of 2 3/4" (72mm)

VS5-2594

Graefe Straight Iris Forcep With Serrations

Serrated Handle With Polished Finish, Straight Tips With Delicate Serrations, And Overall Length Of 2 5/8" (67mm)



VS5-2595

Graefe Straight Iris Forcep

With Teeth, Serrated Handle With Polished Finish, Straight Tips, Fine 1 x 2 Teeth, And Overall Length Of 2 7/8" (7cm). Overall Weight Of 0.2 oz

VS5-3006

Kellan Intraocular Lens Forcep

Round Handle With Delicate Angled 8mm Jaws, Small groove On Interior Portion Of Platform Is Aligned With Handle, And Overall Length Of 3 7/8" (100mm)

VS5-3077

Knolle-Shepard Lens Forcep

Serrated Cross Action Handle With Dull Finish, Curved Shafts, 2mm Wide Serrated Tips, 0.4mm Pin In Upper Jaw, And Overall Length Of 4 3/4" (121mm)

VS5-3081

Sinskey Ridged Intraocular Lens Forcep

Flat Cross Action Handle With Holes And Dull Finish, 45 Degree Angled Tips, 9mm Long x 2mm Wide Lower Jaw, Notched Upper Jaw, And Overall Length Of 4 1/2" (115mm)

VS5-3084

Hall Lens Insertor With Lock

VS5-3086

Shepard Curved Intraocular Lens Forcep,

Serrated Cross Action Handle With Slide Lock And Polished Finish, Curved Shafts, Serrated 1.5mm Upper And 2mm Lower Jaws, And Overall Length Of 3 5/8" (91mm)

VS5-3086-2

Shepard Curved Intraocular Lens Forcep

Serrated And Wide Cross Action Handle With Lock And Polished Finish, Curved Shafts, 0.5mm Notch In Tip, And Overall Length Of 4 3/4" (120mm)

VS5-3086-WOL

Shepard Curved Intraocular Lens Forcep

Without Lock Serrated Cross Action Handle With Polished Finish, Curved Shafts, Serrated 1.5mm Upper And 2mm Lower Jaws, And Overall Length Of 3 5/8" (91mm)

VS5-3087

Osher Seibel Folding Forcep

Serrated Handle Gently Curved Forcep Overall Length Of 4 1/4" (107mm)

IOL Forceps

Ophthalmology



VS0A-224-422

Buratto II Acrylic IOL Implantation Forcep

Angled To Tip 8mm Thin Version With An Overall Length Of 4 1/4" (11cm)

VS5-3092

Universal Folding Forcep

Notched Jaws Align Grasping Platform At The Tip

VS5-3093

Clayman Lens Forcep

Flat Cross Action Handle With Dull Finish, Curved Shafts, 70 Degree Angled 3mm Tips, And Overall Length Of 4 3/4" (119mm)



VS5-3093-S

Clayman Lens Forcep

Serrated Cross Action Handle With Dull Finish, Curved Shafts, 70 Degree Angled 3mm Tips, And Overall Length Of 4 3/4" (119mm)



VS5-3095

Blaydes Angled Lens Forcep

Flat Cross Action Handle With Holes And Dull Finish, Shafts Angled 8mm From Tip To Bend, Smooth Jaws, And Overall Length Of 4 5/8" (117mm)



VS5-3095-S

Blaydes Angled Lens Forcep

Flat Serrated Handle, Cross Action Design With Holes And Dull Finish, Shafts Angled 12mm From Tip To Bend, Smooth Jaws, And Overall Length Of 4 5/8" (117mm)



VS5-3098

Osher Haptic Forcep

Delicate And Used To Lift Haptic Of Lens Out Of Implant Package, Groove in Platform, Angled 45 Degrees, 5mm From Tip, And Overall Length of 3 1/2" (8.2cm). Overall weight of 0.2 oz

VS5-3099

Thornton Intraocular Lens Forcep

Smooth Flat Cross Action Handle Angled 45 Degrees 11mm Jaws Placed In Same Vertical Plane As Handle, And Overall Length Of 3 3/8" (86mm)

VSP112255I

Steinert Universal II Lens Forcep

With Curved Tips, Three Ridges On The Handle, Overall Length Of 4 1/4" (11cm), Overall Weight Of 0.4 oz

VS11-6-8343

Akahoshi Phaco PreChopper

4.5mm Long Delicate Jaws With Sharp Tips On A Round Knurled Handle

VS11-6-8346

Precision Style Phaco PreChopper

On A 3 Hole Handle

VS11-6-8364

Alfonso Nucleus Grasping Forcep

Two Rows Of Delicate Teeth

VS11-6-8366

Precision Style D Nucleus Extraction & Soft IOL Insertion Forcep

VS11-6-9121

Rosenwasser Donor Lamella Inserting Forceps

On A 3 Hole Handle With Thin Vaulted 12mm Long Shafts

VS11-6-9291

Cartridge Loading Forcep

11mm Long Shaft Smooth Jaw

VS5-2094

Landers Vitrectomy Lens Forcep

Polished Finish, 45 Degree Angled Shafts, Curved Grooved Blades For Holding Landers Vitrectomy Lenses, And Overall Length Of 4 1/8" (103mm)



VS132219-2157

Scleral Plug Forcep Wide Serrated Handle

Cross Action Design With Small Cup Shape In The Inside Of Jaw Overall Length Is 4 1/8" (105mm)

Foreign Body Forceps

Ophthalmology



VS5-2071

Schaaf Foreign Body Forcep

Finger Grip Style Handle With Dull Finish, Fine Pointed Grooved Jaws With Smooth Curve From Tip To Handle, And Overall Length Of 3 3/4" (9.5cm). Overall Weight Of 0.3 oz



VS5-2072

Rappazzo Intraocular Foreign Body Forcep

Lightweight Squeeze Handle With Dull Finish, 2mm Jaws, 2mm Inner Jaw, Angled 45 Degree Tip 20 Gauge Diameter Shaft, Overall Length Of 6 3/8" (163mm)



VS5-2072-A

Rappazzo Intraocular Foreign Body Scissors

Lightweight Squeeze Handle With Dull Finish, Sharp Pointed Tips With Guillotine Action, 20 Gauge Diameter Shaft, 1mm Jaw, Overall Length Of 6 3/8" (163mm)



VS5-2072-D

Rappazzo Intraocular Foreign Body Forcep

Lightweight Squeeze Handle With Dull Finish, Stationary Distal Jaw, 20 Gauge Diameter Shaft, 1.5mm Diamond Dusted Jaw, Inner Jaw Space 2mm, 45 Degree Angled Tip, And Overall Length Of 5 3/4" (14.5cm). Overall Weight Of 0.6 oz



VS5-2072-ED

Rappazzo Intraocular Foreign Body Forcep

Lightweight Squeeze Handle, 1.5mm ED Jaw, Tip Angled 45 Degree 20 Gauge Diameter Shaft Jaw, And Overall Length 6 3/8" (164mm)

VS5-2072-S

Rappazzo Intraocular Foreign Body Forcep

Lightweight Squeeze Handle, 1.5mm Jaw, Tip Angled 45 Degree 20 Gauge Diameter Shaft Jaw, And Overall Length 6 3/8" (164mm)



VS5-6489

Watzke Sleeve Spreading Forcep

Serrated Cross-Action Handle With Polished Finish, Serrated 30 Degree Angled 10mm Tips, And Overall Length Of 4 3/4" (12cm). Overall Weight Of 0.9 oz

VS5-6487

Lambert Sleeve Spreader

Tip Angled 90 Degrees Smooth Jaw Cross Action Design With Post And Thumb Screw Locking Device Overall Length 4 1/8" (105mm)



VS5-6763

Carter Sphere Introducer



VS5743K65

Surgical Instrument Handheld Demagnetizing Probe

Instrument Contents

Ophthalmology



Manipulating Hooks	73-82
Needle Holders	83-95
Utility Scissors	96
Super Cut Scissors	97-98
Ring Handle Scissors	99-109
Spring Handle Scissors	110-127
Speculums	128-146
Calipers	147-149
Markers	150-156
Muscle Hooks	157-158
Phaco Chops	159-165
Flexible Iris Retractors	166
Wipes and Spears	167-171
Terms of Sale	172
Order Form	173
Repair Form	174



VS5-1656

Sinskey Iris And IOL Hook

Serrated And Flat Handle With Polished Finish, 0.2mm Hook With Straight Shaft, And Overall Length Of 4 3/4" (12.1cm)



VS5-1656-A

Sinskey Iris And IOL Hook

Serrated And Flat Handle With Polished Finish, 0.2mm Angled Hook, And Overall Length Of 4 5/8" (11.7cm)



VS5-1656-M

Sinskey Hook Modified

Serrated And Flat Handle With Polished Finish, 0.2mm Hook With Curved Shaft, And Overall Length Of 4 3/4" (12cm)



VS5-1656-T-A

Sinskey Iris And IOL Hook Titanium

Serrated And Flat Handle, 0.2mm Angled Hook, And Overall Length Of 4 1/2" (11.4cm)



VS11-4-6331

Sinskey II Hook

0.2mm Straight With A Smooth And Flat Handle, 22mm Shaft, And An Overall Length Of 4 3/4" (12cm)

Manipulating Hooks

Ophthalmology



VS11-4-6341

Sinksey Hook

0.2mm Angled With A Flat And Smooth Handle, 10mm From Tip To Bend, And An Overall Length Of 4 3/4" (12cm)



VS11-4-6343

Sinskey Hook

0.2mm What A Round And Knurled Handle, 18mm Straight Shaft, And An Overall Length Of 4 1/2" (11.5cm)



VS11-4-6291

Sinskey Hook Straight

0.25mm With Flat Handle And Overall Length Of 4 3/4" (12cm)



VS11-4-6203

Sinskey Hook

0.25mm Straight With A Round And Knurled Handle, And An Overall Length Of 4 3/8" (11cm)



VS11-4-6204

Sinskey Hook

0.25mm Angled With A Round An Knurled Handle, 10mm From Tip To Bend, And An Overall Length Of 4 3/8" (11cm)





VS11-4-6201

Harris Sinskey Hook

Angled 0.25mm, With Smooth Flat Handle, 9mm Tip To Bend, And An Overall Length Of 4 3/4" (12cm)



VS5-1653

Clayman Guide

Smooth And Flat Handle With Polished Finish, Angled 45 Degrees 0.2mm Sinskey Hook With Shaft Angled 45 Degrees, 10mm From Tip To Bend, And Overall Length Of 4 5/8" (11.9cm)



VS11-4-6281

Clayman Jaffe Lens Hook

15mm And Angled, With Smooth And Flat Handle, 10mm From Tip To Bend, And An Overall Length Of 4 3/4" (12cm)



VS11-4-6341

Sinksey Hook

0.2mm Angled With A Flat And Smooth Handle, 10mm From Tip To Bend, And An Overall Length Of 4 3/4" (12cm)



VS11-4-6271

Jaffe Hook

0.15mm Straight With A Flat And Smooth Handle, 21mm Shaft, And An Overall Length Of 4 3/4" (12cm)



Manipulating Hooks

Ophthalmology



VS5-1656

Adson Style Bipolar Forcep

Wide Serrated Handle Overall Length 4 3/4" (120.65mm) With A (High Polished 10mm Long Area) With A Tip Size Of 1.5mm. Overall Weight 0.8 oz



VS5-1678

Kuglen Iris Hook And Lens Manipulator

Smooth And Flat Handle With Dull Finish, 0.5mm "K" Hook With Straight Shaft, And Overall Length Of 4 3/4" (11.9cm)



VS5-1678-A

Kuglen Iris Hook And Lens Manipulator

Smooth And Flat Handle With Dull Finish, 0.5mm "K" Hook With Angled Shaft, And Overall Length Of 4 5/8" (117mm)



VS11-4-6621

Kuglen Iris Hook With Clover Leaf

Smooth And Flat Handle, 20mm Straight Shaft, And An Overall Length Of 4 5/8" (11.7cm)



VS11-4-6621-A

Kuglan Iris Hook With Clover Leaf

Smooth And Flat Handle, 10mm From Tip To Bend, And An Overall Length Of 4 1/2" (11.5cm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS11-4-6631

Kuglen Iris Hook Angled 1mm

With Smooth Flat Handle, 10mm Tip To Bend, And An Overall Length Of 4 3/4" (12cm)



VS11-4-6633

Kuglen Hook Push-Pull

With A Round And Knurled Handle, 20mm Straight Shaft, And An Overall Length Of 4 3/4" (12cm)



VS11-4-6634

Kuglen Iris Hook Manipulator Angled

With A Round And Knurled Handle, 10mm From Tip To Bend, And An Overall Length Of 4 3/8" (11cm)



VS5-1682

Lester Lens Manipulator

Round And Knurled Handle With Dull Finish, Straight Shaft With 0.2mm Hourglass Shape Tip, And Overall Length Of 4 1/2" (11.5cm)



VS5-1682-A

Lester Lens Manipulator

Round And Knurled Handle With Dull Finish, Angled Shaft With 0.2mm Hourglass Shape Tip, And Overall Length Of 4 3/8" (11.2cm)

Manipulating Hooks

Ophthalmology



VS11-4-3791

Lester Manipulator Straight With Guard

Knurled And Round Handle, And An Overall Length of 4 3/4" (12cm)



VS11-4-3701

Lester IOL Manipulator Angled

With Round And Knurled Handle, 9mm From Tip To Bend, And An Overall Length Of 4 1/2" (11.5cm)



VS11-4-3702

Lester Manipulator With Vaulted Angle

Round And Knurled Handle, 10mm Curved From Tip To Bend, And An Overall Length Of 4 1/2" (11.5cm)

VS5-1688

Osher "Y" Hook

Round And Knurled Handle With Dull Finish, Straight Shaft With "Y" Shaped Tip, And "Y" Hook, And Overall Length Of 4 3/4" (12cm)



VS11-4-6639

Osher "Y" Hook

With Round And Knurled Handle, 23mm Straight Shaft, And An Overall Length Of 4 1/2" (11.5cm)

VS5-1683

Fenzl Insertion Hook

Round And Knurled Handle With Dull Finish, 45 Degrees Angled 10mm Shaft With 0.2mm Tip, And Overall Length Of 4 3/4" (12cm)



VS5-1866

Fenzl Iris Retractor

Round And Knurled Handle With Polished Finish, Angled 10mm From Tip To Bend With Curved Shaft, And Overall Length Of 4 1/2" (115mm)



VS11-4-6641

Fenzel Lens Hook

0.15mm Straight With A Round And Knurled Handle, 24mm Shaft, And An Overall Length Of 4 1/2" (11.5cm)



VS11-4-6642

Fenzel Lens Hook

0.15mm Angled With A Round And Knurled Handle, 10mm From Tip To Bend, And An Overall Length Of 4 5/8" (11.7cm)



VS11-4-6637

Retsch Boat Hook With Guard

Round And Knurled Handle, 20mm Straight Shaft, And An Overall Length Of 4 1/2" (11.5cm)



VS11-4-6361

Lewicky Lens Hook Angled

With Round And Knurled Handle, 10mm From Tip To Bend, And An Overall Length Of 5" (12.7cm)

Manipulating Hooks

Ophthalmology



VS5-1646

Azar Lens Manipulation Hook

Smooth And Flat Handle With Polished Finish, Angled 90 Degrees 1/5mm Tip, And Overall Length Of 4 5/8" (11.7cm)

VS5-1650

Maumenee Iris Hook

Serrated And Flat Handle With Polished Finish, Sharp 1 1/2mm Hook, And Overall Length Of 4 3/4"(123mm)

VS5-1650-A

Maumenee Iris Hook

Serrated And Flat Handle With Polished Finish, Sharp 1 1/2mm Hook, And Overall Length Of 4 3/4"(123mm)

VS5-1650-B

Shepard Iris Hook

Serrated And Flat Handle With Polished Finish, Curved And Blunt 1mm Hook, And Overall Length Of 4 3/4" (123mm)



VS11-4-6221

Shepard Iris Hook

Blunt With Smooth And Flat Handle, 27mm Straight Shaft, And An Overall Length Of 5" (12.7cm)

VS5-1650-D

Bonn Iris Hook

Serrated And Flat Handle With Dull Finish, Curved And Blunt 3/10mm Hook, And Overall Length Of 4 3/4" (123mm)



VS11-4-6171

Bonn Iris Hook

With A Smooth And Flat Handle, 28mm Shaft, And An Overall Length Of 5" (12.7cm)



VS5-1650-H

Hirschman Iris Hook

Round And Knurled Handle With Polished Finish, Angled 45 Degrees H-Shaped 3/4mm Hook, 9mm From Tip To Bend, And Overall Length Of 122mm, 4 3/4" (12.2cm)



VS5-1657

Precision Iris Hook

Serrated And Flat Handle With Polished Finish, Blunt Hook, And Overall Length Of 4 7/8" (12.3cm)



VS5-1658

Maumenee Iris Hook

Serrated And Flat Handle With Polished Finish, Blunt 1mm Hook, And Overall Length Of 4 7/8" (12.5cm)



VS5-1659

Tyrell Iris Hook

Serrated And Flat Handle With Polished Finish, Sharp 2mm Hook, And Overall Length Of 4 3/4" (12cm)



VS5-1687-M

Minami "M" Hook

Round And Knurled Handle, Multi-Angled Shaft With Teardrop-Shaped Tip, And Overall Length Of 4 3/4" (12cm)

Manipulating Hooks

Ophthalmology



VS11-4-6211

Knolle Iris Hook

Blunt And Straight With A Round And Knurled Handle, A Small Hook, 28mm Shaft, And An Overall Length Of 4 3/4" (12cm)



VS5-1679

Kratz "K" Iris Push-Pull Hook

30 Gauge Cannula With Dull Finish, 23 Gauge Reinforced Shaft Angled 60 Degrees, Front Portion 18mm Long, 14mm From Tip To Bend, "K" Shape Tip On Right Side Of Cannula, And Overall Length Of 4 1/2" (11.5cm)



VS5-4939-CW

Barraquer Needle Holder

Round Knurled Handle With Dull Finish, With Lock, Delicate 9mm Straight Jaws, Overall Length Of 5 1/4" (13.3cm)



VS5-4939-CWO

Barraquer Needle Holder

Round Knurled Handle With Dull Finish, Without Lock, Delicate 9mm Straight Jaws, Overall Length Of 5 1/8" (13cm)



VS5-4939-SW

Barraquer Needle Holder

Round Knurled Handle With Dull Finish, With Lock, Delicate 9mm Straight Jaws, Overall Length Of 5" (12.7cm)



VS5-4939-SWO

Barraquer Needle Holder

Round Knurled Handle With Dull Finish, Without Lock, Delicate 9mm Straight Jaws, Overall Length Of 5 1/4" (13.3cm)



VS5-4953

Barraquer Needle Holder

Round Knurled Large Handle With Dull Finish, Without Lock, Fine 10mm Moderately Curved Jaws, Overall Length Of 5 3/8" (137mm)

VS5-4954

Barraquer Needle Holder

Round Knurled Large Handle With Dull Finish, With Lock, Fine 10mm Moderately Curved Jaws, Overall Length Of 5 3/8" (137mm)

Needle Holders

Ophthalmology



VS5-4955

Barraquer Needle Holder

Round Knurled Medium Handle With Dull Finish, Fine 10mm Gently Curved Jaws, Overall Length Of 5 1/4" (13.5cm)

VS5-4956

Barraquer Needle Holder

Serrated Round Large Handle With Polished Finish, Fine 10mm Gently Curved Jaws, Overall Length Of 4 7/8" (125mm)



VS11-7-3002

Drew's Needle Holder

On Round Knurled Handle Straight Without Lock, With , Mid Screw To Tip Length Of 15mm, Overall Length Of 4" (10cm)



VS11-7-3003

Drew's Needle Holder

Curved With Lock, Knurled Handle, 12mm Mid Screw To Tip Length, Overall Length Of 4 3/4" (12cm)

VS11-7-4111

Barraquer Needle Holder Delicate Short Model

Delicate Straight Jaw, On Round Knurled Handle With Lock Jaw 14mm From Mid Screw To Tip Overall Length Of 3 3/4" (9.5cm)



VS11-7-4121

Barraquer Needle Holder Delicate Short Model


Delicate Straight Jaw, On Round Knurled Handle Without Lock Jaw 14mm From Mid Screw To Tip Overall Length Of 3 3/4" (9.5cm)



VS11-7-4131

Baby Barraquer Needle-Holder

Curved With Lock, Knurled Handle, 12mm Mid Screw To Tip Length, Overall Length Of 4 1/8" (10.5cm)



VS11-7-4141


Baby Barraquer Needle-Holder

Curved Without Lock, With Knurled Handle, 15mm Mid Screw To Tip Length, Overall Length Of 4" (10cm)

5-4940

Troutman Needle Holder

Medium Round And "Horizontally" Serrated Handle With Dull Finish, No Lock, Medium 10mm Gently Curved Jaws, And Overall Length Of 4 3/4" (122mm)



VS5-4940-X

Troutman Needle Holder

Medium Round "Horizontally" Serrated Handle With Dull Finish, Lock, Medium 10mm Gently Curved Jaws, And Overall Length Of 4 3/4" (122mm)



VS11-7-4361

Troutman Needle-Holder

Curved Without Lock, With Knurled Handle, 15mm Mid Screw To Tip Length, And An Overall Length Of 4 7/8" (12.5cm)

Needle Holders

Ophthalmology



VS11-7-4371

Troutman Needle-Holder

Curved With Lock, With Knurled Handle, 15mm Mid Screw To Tip Length, And An Overall Length Of 4 7/8" (12.5cm)

VS11-7-4911

Precision Needle Holder

Extra Delicate Straight Jaw With Lock

VS11-7-4921

Precision Needle Holder

Extra Delicate Straight Jaw Without Lock

VS11-7-4931

Precision Needle Holder

Extra Delicate Curved Jaw With Lock

11-7-4941

Precision Needle Holder

Extra Delicate Curved Jaw Without Lock



11-7-4411

Barraquer Needle-Holder

Delicate And Curved With Lock, Knurled Handle, 12mm Mid Screw To Tip Length, Overall Length Of 5 1/4" (13.5cm)



11-7-4611

Barraquer Needle-Holder

Straight With Lock, Knurled Handle, 13mm Mid Screw To Tip Length, Overall Length Of 5 1/4" (13.5cm)



11-7-4621

Barraquer Needle-Holder

Straight Without Lock, With Knurled Handle, Mid Screw To Tip Length Of 14mm, Overall Length Of 5 1/4" (13.5cm)

VS11-7-4631

Barraquer Needle-Holder

Delicate Curved With Lock, Knurled Handle, 15mm Mid Screw To Tip Length, Overall Length Of 5 1/4" (13.5cm)

VS11-7-4641

Barraquer Needle-Holder

Delicate Curved Without Lock, With Knurled Handle, 15mm Mid Screw To Tip Length, Overall Length Of 5 1/4" (13.5cm)



VS11-7-4731

Barraquer Needle-Holder

Standard Curved With Lock, Knurled Handle, 15mm Mid Screw To Tip Length, Overall Length Of 5 1/4" (13.5cm)



VS11-7-4741

Barraquer Needle-Holder

Curved Without Lock, With Knurled Handle, 15mm Mid Screw To Tip Length, And An Overall Length Of 5 1/4" (13.5cm)

VS5-4906

Precision Needle Holder

Straight With Lock, Large And Round And Knurled Handle With Polished Finish, Medium 9mm Jaws, Overall Length Of 5 5/8" (14.3cm)



VS5-4807

Precision Needle Holder

Straight Without Lock, Large And Round And Knurled Handle With Polished Finish Medium 9mm Jaws, Overall Length Of 5 5/8" (14.3cm)

Needle Holders

Ophthalmology



VS5-4808

Precision Needle Holder

Curved With lock, Large and Round And Knurled Handle With Polished Finish, Lock, Medium 9mm Curved Jaws, Overall Length Of 5 3/8" (13.6cm)



VS5-4808-X

Precision Needle Holder

Curved Without lock, Large and Round And Knurled Handle With Polished Finish Medium 9mm Curved Jaws, Overall Length Of 5 3/8" (13.6cm)



VS5-4808-LC

Precision Needle Holder

Curved With Lock, Large And Round And knurled Handle With Polished Finish Fine 9mm Gently Curved Jaws, Overall Length Of 5 1/8" (130mm)



VS5-4996

Precision Needle Holder

Round Knurled Small Handle With Dull Finish, No Lock, Delicate 6mm Moderately Curved Jaws, 12mm Mid=Screw To Tip Length, And Overall Length Of 4 7/8" (12.5cm)



VS5-4915-WLM

Anis Needle Holder

Round Knurled Handle With Dull Finish, Lock, Delicate 8mm Straight Jaws, And Overall Length Of 4 7/8" (12.2cm)





VS5-4916-M

Anis Needle Holder

Round Knurled Handle Without Lock, Delicate 8mm Gently Curved Jaws, And Overall Length Of 4 5/8" (11.7cm)



VS5-4916-WLM

Anis Needle Holder

Round Knurled Handle With Lock, Delicate 8mm Gently Curved Jaws, And Overall Length Of 4 7/8" (12.3cm)



VS5-4918-W

Osher Needle Holder

Round And Knurled Handle With Lock, Delicate Curved And Tapered 9mm Jaw, Overall Length Of 4 1/4" (110mm)

VS5-4918-WO

Osher Needle Holder

Round And Knurled Handle Without Lock, Delicate Curved And Tapered 9mm Jaw Overall Length Of 4 1/4" (110mm)



VS5-4938

Bechert-Sinskey Needle Holder

Straight With Lock Round And Knurled Large Handle With Dull Finish Delicate 9mm Jaws, Overall Length Of 5 1/8" (12.9cm)



VS5-4938-CW

Bechert-Sinskey Needle Holder


Curved With Lock, Round And Knurled Large Handle With Dull Finish Delicate 9mm Jaws, Overall Length Of 5 1/8" (12.9cm)



VS5-4947-C

Cohan Needle Holder

Round And Knurled Medium Handle With Dull Finish, No Lock, Fine 8mm Gently Mini-Curved Jaws, And Overall Length Of 4" (10.2cm)



Needle Holders

Ophthalmology



VS5-4947-CW

Cohan Needle Holder

Round And Knurled Medium Handle With Dull Finish, Lock, Fine 8mm Gently Mini-Curved Jaws, 18mm Mid Screw To Tip Length, And Overall Length Of 4" (10cm)



VS5-4958-C

Catalano Needle Holder

Round And Knurled Large Handle With Dull Finish, Lock, Fine 9mm Moderately Curved Jaws With Flattened Bottom, And Overall Length Of 5 1/4" (135mm)

VS5-4801

Precision Castroviejo Needle Holder

Straight With Lock, Round And Knurled Handle, Straight 9mm Medium Jaws, And Overall Length Of 5 1/4" (134mm)



VS5-4972-A

Castroviejo Needle Holder

Straight With Lock, Serrated And Flat Handle, 12mm Jaws, 18mm Mid Screw To Tip Length, And Overall Length Of 5 3/4" (14.5cm)



VS5-4972

Castroviejo Needle Holder

Serrated And Flat Handle With Polished Finish, Lock, Medium 10mm Straight Jaws, And Overall Length Of 5 5/8" (14.4cm)



VS5-4972-WH

Castroviejo Needle Holder

Serrated And Wide Handle With Polished Finish, Lock, Fine 9mm Straight Jaws, 14mm Mid Screw To Tip Length, And An Overall Length Of 5 1/2" (140mm)



Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS5-4972-WHMC

Castroviejo Needle Holder

Curved With Lock, Serrated And Wide Handle With Polished Finish, Lock, Fine 9mm Curved Jaws, And Overall Length Of 5 1/8" (13cm)



VS5-4972-WHMS

Castroviejo Needle Holder


Serrated And Wide Handle With Polished Finish, Lock, Fine 9mm Curved Jaws, And Overall Length Of 5 1/8" (13cm)



VS5-4940-XWH

Troutman Needle Holder

Serrated And Wide Handle With Polished Finish, Lock, Fine 9mm Moderately Curved Jaws, And Overall Length Of 5 1/8" (130mm)



VS5-4949-WH

Collis Needle Holder

Serrated And Handle With Lock And Dull Finish, Fine 7mm Jaws Angled 30 Degrees, And Overall Length Of 5" (128mm)



VS5-4963

Castroviejo Needle Holder

Serrated And Flat Handle With Dull Finish, Lock, Medium 9mm Tungsten Carbide Inlaid Straight Jaws, And Overall Length Of 6" (15cm)



Needle Holders

Ophthalmology



VS5-4963-WO

Castroviejo Needle Holder

Tungsten-Carbide Straight With Lock, 9mm Medium Smooth Jaws, 16mm Mid Screw To Tip Length, And Overall Length Of 5 3/4" (14.5cm)



VS5-4973

Castroviejo Needle Holder

Serrated And Flat Handle With Polished Finish, Fine 10mm Straight Jaws, And Overall Length Of 5 3/8" (136mm)



VS5-4959

Castroviejo Needle Holder

Serrated And Flat Handle With Polished Finish, Lock, Medium 11mm Gently Curved Jaws, And Overall Length Of 5 3/8" (137mm)



VS5-4950

Castroviejo Needle Holder

Serrated And Flat Handle With Polished Finish, No Lock, Medium 11mm Gently Curved Jaws, And Overall Length Of 5 3/4" (14.5cm)



VS5-4961

Castroviejo Needle Holder

Serrated And Flat Handle With Polished Finish, Lock, Heavy 12mm Straight Jaws, 18mm Mid Screw To Tip Length, And Overall Length Of 5 3/4" (14.5cm)





VS5-4975

Castroviejo Needle Holder

Serrated And Flat Handle With Polished Finish, Lock, Heavy 16mm Straight Jaws, And Overall Length Of 5 7/8" (14.8cm)



VS5-4962

Azar Needle Holder

Flat And Scalloped Handle With Dull Finish, Lock, Fine 11m Gently Curved Jaws, And Overall Length Of 5 1/8" (13cm)



VS5-4964-S

Knolle Needle Holder

Serrated And Flat Handle With Lock And Polished Finish, Fine 12mm Straight Jaws, 18mm Mid Screw To Tip Length, And Overall Length Of 5 3/4" (14.5cm)



VS5-4945

McPherson Needle Holder

Serrated And Flat Handle With Polished Finish, Lock, Fine 11 mm Straight Jaws, And Overall Length Of 4" (10.2cm)



VS5-4946

McPherson Needle Holder

Serrated And Flat Handle With Polished Finish, Lock, Fine 11 mm Strongly Curved Jaws, And Overall Length Of 4 1/4" (107mm)



Needle Holders

Ophthalmology



VS5-4957

Paton Needle Holder

Serrated And Flat Handle With Satin Finish, Heavy 11mm Serrated Jaws With Length Wise Groove, And Overall Length Of 4 5/8" (11.8cm)



VS5-4809

Kalt Needle Holder

Serrated And Flat Locking Handle With Dull Finish, Heavy 12mm Straight Serrated Jaws, Tungsten Carbide Inserts, And Overall Length Of 5 5/8" (141mm)



VS5-4911

Kalt Needle Holder

Serrated And Flat Locking Handle With Polished Finish, Heavy 12mm Straight Jaws, Serrated Top Jaw Edges, And Overall Length Of 5 5/8" (141mm)



VS5-4913

Kalt Needle Holder

Serrated And Flat Locking Handle With Polished Finish, Medium 10mm Straight Serrated Jaws, And Overall Length Of 5 5/8" (141mm)





VS5-4914

Kalt Needle Holder

Serrated And Flat Locking Handle With Polished Finish, Medium 10mm Serrated Straight 20% Finer Than Original Model Jaws, And Overall Length Of 5 3/8" (13.8cm)



VS11-7-2611

Castroviejo Needle-Holder

Straight With Lock, Serrated Handle, 15mm Mid Screw To Tip Length, And An Overall Length of 5 1/2" (14cm)



VS11-7-2631

Castroviejo Needle-Holder

Curved With Serrated Handle, 13mm Mid Screw To Tip Length, And An Overall Length Of 5 1/2" (14cm)



VS11-7-2911

Kalt Needle Holders

With Smith Jaws, Serrated Handle, 15mm Mid Screw To Tip Length, And An Overall Length Of 5 3/8" (13.7cm)



VS11-7-4862

Anis Needle-Holder

Straight Without Lock, With Knurled Handle, 14mm Mid Screw To Tip Length, And An Overall Length Of 5" (12.7cm)



Utility Scissors

Ophthalmology



7 1/4" Utility Scissors

VS320B	Medicut Utility Scissors 7 1/4"	Blue
VS320BD	Medicut Utility Scissors 7 1/4"	Burgundy
VS320BK	Medicut Utility Scissors 7 1/4"	Black
VS320GR	Medicut Utility Scissors 7 1/4"	Green
VS320HG	Medicut Utility Scissors 7 1/4"	Hunter Green
VS320M	Medicut Utility Scissors 7 1/4"	Magenta
VS320N	Medicut Utility Scissors 7 1/4"	Navy
VS320NO	Medicut Utility Scissors 7 1/4"	Neon Orange
VS320NP	Medicut Utility Scissors 7 1/4"	Neon Pink
VS320NY	Medicut Utility Scissors 7 1/4"	Neon Yellow
VS320O	Medicut Utility Scissors 7 1/4"	Orange
VS320R	Medicut Utility Scissors 7 1/4"	Red
VS320RB	Medicut Utility Scissors 7 1/4"	Royal Blue
VS320TL	Medicut Utility Scissors 7 1/4"	Teal
VS320V	Medicut Utility Scissors 7 1/4"	Purple
VS320W	Medicut Utility Scissors 7 1/4"	White



7 1/4" Utility Scissors

VS321B	Medicut Utility Scissors 5 1/2"	Blue
VS321BD	Medicut Utility Scissors 5 1/2"	Burgundy
VS321BK	Medicut Utility Scissors 5 1/2"	Black
VS321GR	Medicut Utility Scissors 5 1/2"	Green
VS321HG	Medicut Utility Scissors 5 1/2"	Hunter Green
VS321M	Medicut Utility Scissors 5 1/2"	Magenta
VS321N	Medicut Utility Scissors 5 1/2"	Navy
VS321NO	Medicut Utility Scissors 5 1/2"	Neon Orange
VS321NP	Medicut Utility Scissors 5 1/2"	Neon Pink
VS321NY	Medicut Utility Scissors 5 1/2"	Neon Yellow
VS321O	Medicut Utility Scissors 5 1/2"	Orange
VS321R	Medicut Utility Scissors 5 1/2"	Red
VS321RB	Medicut Utility Scissors 5 1/2"	Royal Blue
VS321TL	Medicut Utility Scissors 5 1/2"	Teal
VS321V	Medicut Utility Scissors 5 1/2"	Purple
VS321W	Medicut Utility Scissors 5 1/2"	White

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS19-133110

Precision Iris Scissors Straight

4 1/2" (14.5cm) Pure Super Cut With Black Colored Rings For Identification



VS19-133111

Precision Iris Scissors Curved

4 1/2" (14.5cm) Pure Super Cut With Black Colored Rings For Identification



VS19-133142

Precision Plastic Scissors Straight

4 7/8" (12.5cm) Pure Super Cut With Black Colored Rings For Identification.

Super Cut Scissors

Ophthalmology



VS19-133143

Precision Plastic Scissors Curved

4 7/8" (12.5cm) Pure Super Cut With Black Colored Rings For Identification



VS19-133166

Precision Stevens Tenotomy Scissors Straight

4 1/2" (11.5cm) Pure Super Cut With Black Colored Rings For Identification



VS19-133168

Precision Stevens Tenotomy Scissors Curved

4 1/2" (11.5cm) Pure Super Cut With Black Colored Rings For Identification



VS5-4407

Precision Straight Iris Scissors

Ribbon Style Ring Handle With Polished Finish On Blades, Pointed Tips, 18mm From Mid-Screw To Tip, And Overall Length Of 3 7/8" (98mm)



VS5-4408

Precision Curved Iris Scissors

Ribbon Style Ring Handle With Polished Finish On Blades, Pointed Tips, 18mm From Mid-Screw To Tip, And Overall Length Of 3 7/8" (98mm)

VS5-4501

Aebli Straight Corneal Scissors

Ring Handle, Universal Standard Pattern, Straight Blades With Polished Finish, Blunt Tips, 19mm Mid Screw To Tip, And Overall Length Of 3 5/8" (92mm)



VS5-4511

Precision Straight Iris Scissors

Ribbon Style Ring Handle, Polished Finish On Blades, Pointed Tips, 25mm Mid Screw To Tip, And Overall Length Of 4 1/8" (105mm)

Ring Handle Scissors

Ophthalmology



VS5-4513

Precision Curved Iris Scissors

Ribbon Style Ring Handle, Polished Finish On Blades, Pointed Tips, 25mm Mid Screw To Tip, And Overall Length Of 4 1/8" (105mm)



VS5-4515

Precision Straight Iris Scissors

Ring Handle, Polished Finish On Blades, Pointed Tips, 28mm Mid Screw To Tip, And Overall Length Of 3 7/8" (100mm)



VS5-4516

Precision Curved Iris Scissors

Ring Handle, Polished Finish On Blades, Pointed Tips, 28mm Mid Screw To Tip, And Overall Length Of 3 7/8" (100mm)



VS5-4517

Precision Straight Iris Scissors

Ring Handle, Polished Finish On Blades, Light Standard Pattern, Pointed Tips, 25mm Mid Screw To Tip, And Overall Length Of 3 7/8" (98mm)



VS5-4519

Precision Curved Iris Scissors

Ring Handle, Polished Finish On Blades, Light Standard Pattern, Pointed Tips, 25mm Mid Screw To Tip, And Overall Length Of 3 7/8" (98mm)



VS5-4521

Precision Straight Eye Scissors

Ribbon Style Ring Handle, Polished Finish On Blades, Pointed Tips, 29mm Mid Screw To Tip, And Overall Length Of 4 1/4" (109mm)

Ring Handle Scissors

Ophthalmology



VS5-4523

Precision Curved Eye Scissors

Ribbon Style Ring Handle, Polished Finish On Blades, Pointed Tips, 29mm Mid Screw To Tip, And Overall Length Of 4 1/4" (109mm)



VS5-4527

Precision Straight Eye Scissors

Ribbon Style Ring Handle With Polished Finish, Blunt Tips, 29mm Mid Screw To Tip, And Overall Length Of 4 1/4" (107mm)



VS5-4529

Precision Curved Eye Scissors

Ribbon Style Ring Handle With Polished Finish, Blunt Tips, 30mm Mid Screw To Tip, And Overall Length Of 4 1/4" (107mm)



VS5-4535

Precision Straight Eye Scissors

Ring Handle, Polished Finish On Blades, Sharp Tips, 32mm Mid Screw To Tip, And Overall Length Of 4 1/4" (109mm)



VS5-4537

Precision Curved Eye Scissors

Ring Handle, Polished Finish On Blades, Sharp Tips, 32mm Mid Screw To Tip, And Overall Length Of 4 1/4" (109mm)



VS5-4547

Precision Eye Scissors

Ring Handle, Polished Finish On Blades, Blunt Probe Points, Angled Shafts, 26mm Mid Screw To Tip, And Overall Length Of 4 5/8" (116mm)

Ring Handle Scissors

Ophthalmology



VS5-4561

Wilmer Angled Conjunctival And Utility Scissors

Ring Handle, Serrated Blades With Polished Finish, Sharp Pointed Tips, 18mm Mid Screw To Tip, And Overall Length of 3 3/4" (94mm)



VS5-4571

McReynolds Pterygium Scissors

Ring Handle, Polished Finish On Blades, Blunt Angled Tips, 18mm Mid Screw To Tip, And Overall Length Of 4 1/8" (105mm)

VS5-4591

Aebli Straight Corneal Scissors

Ribbon Style Ring Handle, Polished Finish On Blades, Blunt Tips, 21mm Mid Screw To Tip, And Overall Length Of 3 7/8" (100mm)

VS5-4593

Aebli Angled Corneal Scissors

Ribbon Style Ring Handle, Polished Finish On Blades, Blunt Tips, 21mm Mid Screw To Tip, And Overall Length Of 3 7/8" (100mm)



VS5-4601

Precision Curved Stitch Scissors

Precision Curved Stitch Scissors, Ribbon Style Ring Handle With Polished Finish, Sharp Pointed Tips, 18mm Mid Screw To Tip, And Overall Length Of 3 3/4" (97mm)

VS5-4603

Precision Original Stitch Scissors

Ribbon Style Ring Handle With Polished Finish, Sharp Pointed Tips, 22mm Mid Screw To Tip, And Overall Length Of 4" (102mm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS5-4606

Reeh Stitch Scissors

Ribbon Style Ring Handle With Polished Finish, Sharp Pointed Tips, Very Small Hook At End Of One Blade, 19mm Mid Screw To Tip, And Overall Length Of 3 7/8" (98mm)

VS5-4617

McGuire Corneal Scissors

Ring Handle, Polished Finish On Angled And Curved Blades, Right Orientation, 21mm Mid Screw To Tip, And Overall Length Of 4 1/8" (104mm)

VS5-4621

Anis Micro Scissors

Squeeze Action Handle, 20mm Long 23 Gauge Very Fine Shaft, And Overall Length Of 6 1/4" (160mm)

VS5-4661

Stevens Straight Tenotomy Scissors

Ribbon Style Ring Handle, Polished Finish On Straight Blades, Heavy Model, Blunt Tips, 19mm Mid Screw To Tip, And Overall Length Of 3 3/4" (98mm)



VS5-4663

Stevens Curved Tenotomy Scissors

Ribbon Style Ring Handle, Polished Finish On Curved Blades, Heavy Model, Blunt Tips, 19mm Mid Screw To Tip, And Overall Length Of 3 3/4" (98mm)

Ring Handle Scissors

Ophthalmology



VS5-4665

Stevens Straight Tenotomy Scissors

Ribbon Style Ring Handle, Polished Finish on Straight Blades, Light Model, Blunt Tips, And Overall Length Of 3 3/4" (98mm)



VS5-4667

Stevens Curved Tenotomy Scissors

Ribbon Style Ring Handle, Polished Finish On Curved Blades, Light Model, Blunt Tips, 19mm Mid Screw To Tip, And Overall Length Of 3 3/4" (98mm)



VS5-4671

Stevens Straight Tenotomy Scissors

Ring Handle, Polished Finish On Straight Blades, Standard Model, Blunt Tips, 30mm Mid Screw To Tip, Overall Length Of 4 3/8" (112mm)



VS5-4673

Stevens Curved Tenotomy Scissors

Ring Handle, Polished Finish On Curved Blades, Standard model, Blunt Tips, 30mm Mid Screw To Tip, And Overall Length Of 4 3/8" (112mm)

VS5-4681

Precision Straight Strabismus Scissors

Ribbon Style Ring Handle, Polished Finish On Rounded Blades, Blunt Tips, 29mm Mid Screw To Tip, And Overall Length Of 4 1/4" (109mm)



VS5-4682

Precision Curved Strabismus Scissors

Ribbon Style Ring Handle, Polished Finish On Rounded Blades, Blunt Tips, 29mm Mid Screw To Tip, And Overall Length Of 4 1/4" (109mm)



VS5-4683

Precision Curved Strabismus Scissors

Ribbon Style Ring Handle, Polished Finish On Rounded Blades, Strong Curve, Blunt Tips, 27mm Mid Screw To Tip, And Overall Length Of 4 1/4" (106mm)

VS5-4687

Knapp Straight Strabismus Scissors

Ring Handle, Polished Finish On Rounded Blades, Blunt Or Lightly Curved Blunt Tips, 25mm Mid Screw To Tip, And Overall Length Of 4 5/8" (118mm)

Ring Handle Scissors

Ophthalmology



VS5-4688

Knapp Curved Strabismus Scissors

Ring Handle, Polished Finish On Rounded Blades, Blunt Or Lightly Curved Blunt Tips, 25mm Mid Screw To Tip, And Overall Length Of 4 5/8" (118mm)



VS5-4689

Knapp Curved Strabismus Scissors

Ring Handle, Polished Finish On Rounded Blades, Strong Curve, Blunt Tips, 25mm Mid Screw To Tip, And Overall Length Of 4 5/8" (118mm)



VS5-4691

Littler Curved Dissecting Scissors

Ring Handle With Polished Finish, Lightly Curved, Blunt Tips, Rounded Blades, Eye In Blades In Place Of Eye Probe, 27mm Mid Screw To Tip, And Overall Length Of 4 3/4" (119mm)



VS5-4695

Precision Straight Utility Scissors,

Serrated And Flat Handle With Polished Finish, Blunt Tips, 25mm Mid Screw To Tip, And Overall Length Of 5 1/8" (131mm)
4 1/2" (14.5cm) Pure Super Cut With Black Colored Rings For Identification

VS5-4696

Precision Curved Utility Scissors

Ribbon Style Ring Handle With Polished Finish, Blunt Tips, 24mm Mid Screw To Tip, And Overall Length Of 4 1/4" (109mm)



VS5-4699

Precision Curved Light Scissors

Ribbon Style Ring Handle With Polished Finish, Blunt Tips, 18mm Mid Screw To Tip, And Overall Length Of 7/8" (99mm)

Ribbon Style Ring Handle With Polished Finish, Blunt Tips, 24mm Mid Screw To Tip, And Overall Length Of 4 1/4" (109mm)

Spring Handle Scissors

Ophthalmology



VS5-4497

Vannas Straight Capsulotomy Scissors

Serrated Handle With Polished Finish, 6mm Blades, Sharp Pointed Tips, 8mm Mid Screw To Tip, And Overall Length Of 3 1/4" (84mm)



VS11-5-6111

Vannas Scissor

Straight Serrated Handle With Sharp 5mm Blades

VS5-4497-B

Vannas Straight Capsulotomy Scissors

Serrated Handle With Polished Finish, 6mm Blades, Blunt Tips, 8mm Mid Screw To Tip, And Overall Length Of 3 1/4" (84mm)



VS5-4498

Vannas Curved Capsulotomy Scissors

Serrated Handle With Polished Finish, 6mm Blades, Sharp Pointed Tips, 8mm Mid Screw To Tip, And Overall Length Of 3 1/8" (8cm)



VS11-5-6121

Vannas Scissors

With Curved Blades And Sharp Tips, Serrated Handle, Mid Screw To Tip length Of 10mm, And An Overall Length Of 3 1/4" (8.3cm)

VS5-4498-B

Vannas Curved Capsulotomy Scissors

Serrated Handle With Polished Finish, 6mm Blades, Blunt Tips, 8mm Mid Screw To Tip, And Overall Length Of 3 1/4" (84mm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS5-4494-S

Clayman-Vannas Straight Scissors

Smooth Handle With Polished Finish, Sharp Pointed Tips, 7mm Mid Screw To Tip, And Overall Length Of 3 1/4" (83mm)

VS5-4494-C

Clayman-Vannas Curved Scissors

Smooth Handle With Polished Finish, Sharp Pointed Tips, 7mm Mid Screw To Tip, And Overall Length Of 3 1/4" (83mm)



VS11-5-6161

Gills-Vannas Extra Long Thin Straight Scissors

Blade Length Of 7mm, Overall Length Of 3 1/8" (8.1cm)

VS11-5-6171

Gills-Vannas Extra Long Thin Curved Scissors

Blade Length Of 7mm, Overall Length Of 3 1/8" (8.1cm)



VS5-4499

Vannas Angled Capsulotomy Scissors

Serrated Handle With Polished Finish, 6mm Blades, Sharp Pointed Tips, 9mm Mid Screw To Tip, And Overall Length Of 3 3/8" (8.5cm)



VS11-5-6131

Vannas Scissors Angled Forward

Scissors Blade Length Of 7mm, Overall Length Of 3" (7.8cm)

Spring Handle Scissors

Ophthalmology



VS5-4499-GW

Gills-Welsh Vannas Angled Micro Scissors

Serrated Handle With Polished Finish, 5mm Blades, Sharp Pointed Tips, 11mm Mid Screw To Tip, And Overall Length Of 3 3/8" (86mm)



VS11-5-6181

Gills-Vannas Extra Long Thin Angled Forward Scissors

Blade Length Of 7mm, Overall Length Of 3 1/8" (8.1cm)



VS5-4499-GWL

Gills-Welsh Angled Vannas Micro Scissors

Serrated Handle With Polished Finish, 11mm Blades, Sharp Pointed Tips, 13mm Mid Screw To Tip, And Overall Length Of 3 3/8" (86mm)



VS5-4490

Vannas Angled Capsulotomy Scissors

Serrated Handle With Polished Finish, 6mm Blades, Sharp Pointed Tips, 9mm Mid Screw To Tip, And Overall Length Of 3 1/4" (8cm)



VS11-5-6141

Vannas Scissors

Angled To The Side, Scissors Blade Length Of 7mm, Overall Length Of 3" (7.8cm)

VS5-4450

Cohan-Vannas Curved Iris Scissors

Serrated And Flat Handle With Dull Finish, Sharp Pointed Tips, 12mm Mid Screw To Tip, And Overall Length Of 2 1/4" (57mm)



VS5-4454

Straight Iris Scissors

Serrated And Flat Handle With Polished Finish, Sharp Pointed Tips, 21mm Mid Screw To Tip, And Overall Length Of 4 1/2" (11.5cm)



VS5-4455

Curved Iris Scissors

Serrated And Flat Handle With Polished Finish, Sharp Pointed Tips, 21mm Mid Screw To Tip, And Overall Length Of 4 1/2" (11.5cm)



VS5-4458

McPherson-Vannas Curved Iris Scissors

Serrated And Flat Handle With Polished Finish, Sharp Pointed Tips, 11mm Mid Screw To Tip, And Overall Length Of 3 3/8" (8.5cm)



VS11-5-5511

Iris Scissors With Small Blades And Sharp Tips

On A 3 Hole Handle With Small Straight Blades Mid Screw To Tip Length Of 10mm Overall Length Of 4 1/4" (11cm)



VS11-5-5521

Iris Scissors

On A 3 Hole Handle With Small Curved Blades Mid Screw To Tip Length Of 10mm Overall Length Of 4 1/4" (11cm)



Spring Handle Scissors

Ophthalmology



VS5-4451

Precision Fine Straight Stitch Scissors

Serrated And Flat Handle With Polished Finish, Sharp Pointed Tips, 17mm Mid Screw To Tip, And Overall Length Of 4 1/4" (10.7cm)



VS11-5-5611

Fine Stitch Scissors Medium

With Straight And Sharp Blades, Three-Holed Handle, Mid Screw To Tip Length Of 20mm, And An Overall Length Of 4 3/4" (12cm)



VS5-4451-WH

Precision Fine Straight Stitch Scissors

Serrated And Wide Handle With Polished Finish, Sharp Pointed Tips, 16mm Mid Screw To Tip, And Overall Length Of 5 1/4" (13cm)



VS11-5-5621

Fine Stitch Scissors Medium With Curved And Sharp Tips

Three-Holed Handle, Mid Screw To Tip Length Of 20mm, And An Overall Length Of 4 3/4" (12cm)



VS11-5-6221

Precision-Vannas Scissors With Straight And Extra-Thin Tips,

Three-Holed Handle, Mid Screw To Tip Length Of 10mm, Blade Length 7mm And An Overall Length Of 4 1/2" (11cm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS11-5-6222

Precision-Vannas Scissors With Curved And Extra-Thin Tips,
Three-Holed Handle, Mid Screw To Tip Length Of 10mm, Blade Length 7mm And An Overall Length Of 4 1/2" (11cm)



VS11-5-6223

Precision-Vannas Scissors With Angled Forward And Extra-Thin Tips
Three-Holed Handle, Mid Screw To Tip Length Of 10mm, Blade Length 7mm And An Overall Length Of 4 1/2" (11cm)

VS5-4401-A

Capsule Scissors 10mm Angled Blades
Delicate 10mm Blades Angled 45 Degrees, And Overall Length of 4 1/4" (110mm)

VS5-4401-C

Capsule Scissors 10mm Curved Blades
Delicate 10mm Blade Angled 45 Degrees, Gently Curved Blades, And Overall Length Of 4 1/4" (110mm)



VS11-5-6231

Precision-Vannas Scissors With Straight And Extra-Thin Tips
Three-Holed Handle, Mid Screw To Tip Length Of 13mm, Blade Length 10mm And An Overall Length Of 4 1/2" (11cm)



VS11-5-6232

Precision-Vannas Scissors With Curved And Extra-Thin Tips
Three-Holed Handle, Mid Screw To Tip Length Of 13mm, Blade Length 10mm And An Overall Length Of 4 1/2" (11cm)

Spring Handle Scissors

Ophthalmology



VS11-5-6233

Precision Stern-Gills Style Vannas Scissor With Angled Forward And Straight Extra-Thin Tips
Three-Holed Handle, Mid Screw To Tip Length Of 13mm, Blade Length 10mm And An Overall Length Of 4 1/2" (11cm)



VS11-5-6234

Precision Uribe-Stern Style Vannas Scissor
Angled Forward Blades With Curved Extra-Thin Tips On A 3 Hole Handle, Mid Screw To Tip Length Of 13mm, Blade Length 10mm And An Overall Length Of 4 1/2" (11cm)

VS11-5-6234-XC NEW Scissors Design

Precision Uribe-Stern (XC) Style Vannas Scissor
Angled Forward With Extra Curved Extra-Thin Blades Tips On A 3 Hole Handle, Mid Screw To Tip Length Of 13mm, Blade Length 10mm And An Overall Length Of 4 1/2" (11cm)



VS11-5-6235

Precision Hoffer-Stern Style Vannas Scissor Angled Forward Vaulted Curved Extra-Thin Sharp Tips
Three-Holed Handle, Mid Screw To Tip Length Of 13mm, Blade Length 10mm And An Overall Length Of 4 1/2" (11cm)

VS11-5-6236

Precision Gills Capsulotomy Scissor
Angled Forward Blades With Curved Extra-Thin Sharp Tips On A 3 Hole Handle Mid Screw To Tip Length Of 13mm, Blade Length 10mm And An Overall Length Of 4 1/2" (11cm)

VS11-5-6237

Precision Uthoff-Gills Style Vannas Scissor
Angled Forward Blades With Curved Extra-Thin Blunt Tips, On A 3 Hole Handle Mid Screw To Tip Length Of 13mm, Blade Length 10mm And An Overall Length Of 4 1/2" (11cm)

VS11-5-6239

Precision Dahan Style Vannas Scissor Angled To The Side And Curved Extra-Thin Blunt Tips
Three-Holed Handle, Mid Screw To Tip Length Of 12mm, Blade Length 9mm And An Overall Length Of 4 1/2" (11cm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS5-4431-R

Westcott Curved Tenotomy Scissors

Right, Serrated And Flat Handle With Polished Finish, Rounded Westcott Blades, Blunt Tips, 19mm Mid Screw To Tip, And Overall Length Of 4 1/2" (11.4cm)



VS5-4431-L

Westcott Curved Tenotomy Scissors

Left, Serrated And Flat Handle With Polished Finish, Rounded Westcott Blades, Blunt Tips, 19mm Mid Screw To Tip, And Overall Length Of 4 1/2" (11.4cm)



VS5-4431-RWH

Westcott Curved Tenotomy Scissors

Serrated And Wide Handle With Polished Finish, Right Model, Rounded Westcott Blades, Blunt Tips, 23mm Mid Screw To Tip, And Overall Length Of 5 1/2" (14cm)



VS11-5-4115

Westcott Tenotomy Scissors,

On A 3 Hole Handle, Curved Blunt Tips Medium Blades Length Of 15mm, Overall Length Of 4 1/4" (10.9cm)



VS5-4431-RS

Shepard-Westcott Curved Tenotomy Scissors

Serrated And Flat Handle With Polished Finish, Right Model, Westcott Blades, Blunt Tips, Serrated Upper Blade, 21mm Mid Screw To Tip, And Overall Length Of 4 5/8" (11.7cm)

Spring Handle Scissors

Ophthalmology



VS11-5-4111

McPherson–Westcott Conjunctival Scissors

On A 3 Hole Handle Curved, Blunt Tips, Mid Screw To Tip Length Of 3/8" (10mm) Overall Length Of 4 3/8" (11cm)

VS5-4431-RX

Westcott Curved Tenotomy Scissors

Serrated Handle With Polished Finish, Right Sub–Miniature Model, Round Westcott Blades, Blunt Tips, 9mm Mid Screw To Tip, And Overall Length Of 3 3/8" (88mm)

VS5-4431-RXC

Cohan–Westcott Curved Tenotomy scissors

Serrated Handle With Dull Finish, Sub–Miniature Model, Westcott Blades, 8mm Mid Screw To Tip, And Overall Length Of 2 1/8" (54mm)



VS5-4354

Clayman–Westcott Curved Scissors

Smooth Handle With Dull Finish, Blunt Tips, 11mm From Mid–Screw To Tip, And Overall Length Of 4 1/4" (107mm)



VS5-4432

Westcott Type Stitch Scissors

Serrated And Flat Handle With Polished Finish, Sharp Pointed Tips, 21mm Mid Screw To Tip, And Overall Length Of 4 3/4" (12cm)



VS5-4432-WH

Westcott Type Stitch Scissors

Serrated And Wide Handle With Polished Finish, Sharp Pointed Tips, 22mm MidScrew To Tip, And Overall Length Of 5 1/2" (140mm)



VS11-5-5111

McPherson Westcott Stitch Scissor

With Small And Very Sharp Blade, Three-Holed Handle, Mid Screw To Tip Length Of 3/8" (10mm), And An Overall Length Of 4 3/8" (11cm)



VS11-5-5261

Jaffe Stitch Scissors

On A 3 Hole Handle, Medium Curved Blades Very Sharp tips, 1/2" (1.3cm) Mid Screw To Tip Length, And Overall Length Of 4 1/2" (11.5cm)



VS5-4432-MWH

Castroviejo Mini Keratoplasty-Utility-Stitch Scissors

Serrated And Wide Handle With Polished Finish, Very Small Pointed Tips, 13mm Mid Screw To Tip, And Overall Length Of 4 7/8" (12.4cm)

VS5-4432-C

Collis Micro Utility Scissors

Serrated And Wide Handle With Polished Finish, Curved, Sharp Pointed Tips, 12mm From Mid Screw To Tip, And Overall Length Of 5" (128mm)



VS5-4433

Westcott Type Utility Scissors

Serrated And Flat Handle With Polished Finish, Rounded Sharp Tips, 25mm Mid Screw To Tip, And Overall Length Of 4 1/2" (11.5cm)

Spring Handle Scissors

Ophthalmology



VS5-4427

Precision Westcott Conjunctival Scissors

Serrated And Wide Handle With Polished Finish, Lightly Curved Blunt Tips, 16mm Mid Screw To Tip, And Overall Length Of 4 7/8" (12.4cm)

VS5-4426

Clayman–Westcott Curved Scissors

Smooth Handle With Dull Finish, Blunt Tips, 11mm From Mid–Screw To Tip, And Overall Length Of 4 1/4" (107mm)



VS5-4477

Barraquer Iris Scissors

Serrated Grip With Small Small Hole In Center, Squeeze Action Handle With Polished Finish, Blunt Tips, 8mm From Bend To Tip, And Overall Length Of 2 1/8" (54mm)

VS5-4478

Barraquer Vitreous Scissors

Serrated Grip With Small Small Hole In Center, Squeeze Handle With Polished Finish, Delicate Angled 10mm Blades, Blunt Tips, And Overall Length Of 2 1/8" (54mm)



VS11-5-5011

Barraquer Iris Scissor Round Looped Fenestrated Squeeze Action Handles

Blades Are Angled With A 7mm Blunt Tips

VS5-4470

Nadler Superior Radial Scissor

Serrated Grip With Small Small Hole In Center, Squeeze Handle With Polished Finish, Delicate Angled 7mm Blade, Blunt Tips, And Overall Length Of 2 1/8" (54mm)

5-4471

DeWecker Iris Scissors

Serrated And Flat Squeeze Action Handle With Dull Finish, Blunt And Sharp Points, 12mm From Bend To Tip, And Overall Length Of 4 1/2" (114mm)

5-4472

DeWecker Iris Scissors

Serrated And Flat Squeeze Action Handle With Dull Finish, Blunt And Blunt Points, 11mm From Tip To Bend, And Overall Length Of 4 1/2" (114mm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020

VS5-4415

Castroviejo Corneal Section Scissors

Right, Serrated And Flat Handle With Polished Finish, Lower Blade 0.5mm Longer Than Upper Blade, 14mm Mid Screw To Tip, And Overall Length Of 4" (10.3cm)



VS5-4415-M

Castroviejo Corneal Section Scissors

Right, Serrated And Flat Handle With Polished Finish, Mid-Miniature Model With Curved Blades, Blunt Tips, Lower Blade 0.5mm Longer Than Upper Blade, 13mm Mid Screw To Tip, And Overall Length Of 4" (10cm)

VS5-4415-MC

McPherson-Castroviejo Miniature Corneal Section Scissors

Right, Serrated And Flat Handle With Polished Finish, Miniature Model Has Curved Blades, Blunt Tips, Lower Blade 0.5mm Longer Than Upper Blades, 13mm From Mid Screw To Tip, And Overall Length Of 3 3/4" (93mm)

VS5-4415-WH

Castroviejo Corneal Section Scissors

Right, Serrated And Wide Handle With Polished Finish, Miniature Model Had Curved Blades, Blunt Tips, 10mm Mid Screw To Tip, And Overall Length Of 4 7/8" (125mm)



VS5-4415-X

Castroviejo Corneal Section Scissors

Right, Serrated Handle With Polished Finish, Sub-Miniature Model Has Curved Blades, Blunt Tips, Lower Blade 0.5mm Longer Than Upper Blades, 9mm Mid Screw To Tip, And Overall Length Of 4 1/8" (103mm)

VS5-4415-XWH

Castroviejo Corneal Section Scissors

Right, Serrated And Wide Handle With Polished Finish, Sub-Miniature Model Had Curved Blades, Blunt Tips, Lower Blade 0.5mm Longer Than Upper Blade, 7mm Mid Screw To Tip, And Overall Length Of 4 3/4" (120mm)

Spring Handle Scissors

Ophthalmology



VS5-4416

Castroviejo Corneal Section Scissors

Left, Serrated And Flat Handle With Polished Finish, Lower Blade 0.5mm Longer Than Upper Blade, 14mm Mid Screw To Tip, And Overall Length Of 4 1/8" (10.5cm)



VS5-4416-M

Castroviejo Corneal Section Scissors

Left, Serrated And Flat Handle With Polished Finish, Mid-Miniature Model With Curved Blades, Blunt Tips, Lower Blade 0.5mm Longer Than Upper Blade, 13mm Mid Screw To Tip, And Overall Length Of 4" (10.2cm)

VS5-4416-MC

McPherson-Castroviejo Miniature Corneal Section scissors

Left, Serrated And Flat Handle With Polished Finish, Miniature Model Has Curved Blades, Blunt Tips, Lower Blade 0.5mm Long Than Upper Blade, 13mm Mid Screw To Tip, And Overall Length Of 3 3/4" (93mm)

VS5-4416-WH

Castroviejo Corneal Section Scissors

Left, Serrated And Wide Handle With Polished Finish, Miniature Model Has Curved Blades With Blunt Tips, 10mm Mid Screw To Tip, And Overall Length of 4 7/8" (125mm)



VS5-4416-X

Castroviejo Corneal Section Scissors

Left, Serrated Handle With Polished Finish, Sub-Miniature Model Has Curved Blades, Blunt Tips, Lower Blade 0.5mm Longer Than Upper Blade, 9mm Mid Screw To Tip, And Overall Length Of 4 1/8" (103mm)

VS5-4416-XWH

Castroviejo Corneal Section Scissors

Left, Serrated And Wide Handle With Polished Finish, Sub-Miniature Model Has Curved Blades, Blunt Tips, Lower Blade 0.5mm Longer Than Upper Blade, 7mm Mid Screw To Tip, And Overall Length Of 4 3/4" (120mm)

VS5-4411

Castroviejo Corneal Section Scissors

Castroviejo Corneal Section Scissors, Right Serrated Flat Handle, Polished Finish, Angled Blade And Curved, Lower Blade 1mm Long Than Upper Blade And Rounded Like Spatula, 18mm MidScrew To Tip, And Overall Length Of 4 1/4" (110mm)



VS5-4412

Castroviejo Corneal Section Scissors

Left Serrated Flat Handle, Polished Finish, Angulated Blade And Curved, Lower Blade 1mm Longer Than Upper Blade And Rounded Like Spatula, 18mm MidScrew To Tip, And Overall Length Of 4 1/4" (110mm)



VS5-4417

Troutman-Castroviejo Corneal Section Scissors

Right, Serrated And Flat Handle With Polished Finish, Barraquer Stop/Lock, Miniature Model Has Curved Blades, Blunt Tips, Both Blades Of Equal Length, And Overall Length Of 4 5/8" (116mm)



VS5-4418

Troutman-Castroviejo Corneal Section Scissors

Left, Serrated And Flat Handle With Polished Finish, Barraquer Stop/Lock, Miniature Model Has Curved Blades, Blunt Tips, Equal Length 10mm Blades, 13mm Mid Screw To Tip Length, And Overall Length of 4" (10cm)



VS5-4329-R

Castroviejo Curved Corneal Scissors

Right, Serrated And Flat Handle With Polished Finish, Miniature Model With Blunt Tips, 11mm From Mid-Screw To Tip, And Overall length Of 4 1/8" (10.5cm)



VS5-4329-L

Castroviejo Curved Corneal Scissors

Left, Serrated And Flat Handle With Polished Finish, Miniature Model With Blunt Tips, 11mm From Mid-Screw To Tip, And Overall Length Of 4 1/8" (10.5cm)

Spring Handle Scissors

Ophthalmology



VS5-4329-MM

Castroviejo Mini Corneal Scissors

Serrated And Flat Handle With Polished Finish, Universal Miniature model With Blunt Tips, 11mm From Mid-Screw To Tip, And Overall Length Of 4 1/4" (106mm)



VS5-4342-R

Precision Miniature Corneal Scissors

Right, Delicate Corneal Scissors, 9mm Blades With Lower Blade Extended 0.5mm, And Overall Length Of 3 7/8" (100mm)

VS5-4342-L

Precision Miniature Corneal Scissors

Left, Delicate Corneal Scissors, 9mm Blades With Lower Blade Extended 0.5mm, And Overall Length Of 3 7/8" (100mm)

VS11-5-2111

Castroviejo Corneal Scleral Scissor

On 3 Hole Handle Small Blades Curved Right

VS11-5-2121

Castroviejo Corneal Scleral Scissor

On 3 Hole Handle Small Blades Curved Left

VS11-5-2131

Castroviejo Corneoscleral Scissors Medium

Right, On 3 Hole Handle, Mid Screw To Tip Length Of 12mm, And An Overall Length Of 4" (10cm)

VS11-5-2141

Castroviejo Corneoscleral Scissors Medium

Left, On 3 Hole Handle, Mid Screw To Tip Length Of 12mm, And An Overall Length Of 4" (10cm)

VS11-5-2311

Troutman Castroviejo Corneoscleral Scissor

On 3 Hole Handle With Stop Small Blades Curved Right

VS11-5-2321

Troutman Castroviejo Corneoscleral Scissor

On 3 Hole Handle With Stop Small Blades Curved Left

VS11-5-2331

Troutman Castroviejo Corneoscleral Scissor

On 3 Hole Handle With Stop Medium Blades Curved Right

VS11-5-2341

Troutman Castroviejo Corneoscleral Scissor

On 3 Hole Handle With Stop Medium Blades Curved Left

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020

VS11-5-3121

Troutman Katzin Corneal Transplant Scissor
On 3 Hole Handle With Stop Small Blades Curved Left

VS11-5-3125

Rich Corneal Transplant Scissor
On 3 Hole Handle Extra Delicate Blades Curved Right

VS11-5-3126

Rich Corneal Transplant Scissor
On 3 Hole Handle Extra Delicate Blades Curved Left



VS5-4344

Katzin Curved Corneal Transplant Scissors

Right, Serrated And Flat Handle With Polished Finish, Miniature Model With Strongly Curved Blunt Blades Of Uneven Lengths, Lower Blade 12mm, Upper Blade 11mm (From Mid-Screw To Tip), And Overall Length Of 4" (10cm)

VS11-5-3131

Katzin Corneal Transplant Scissor
On 3 Hole Handle Medium Blades Curved Right



VS5-4343

Katzin Curved Corneal Transplant Scissors

Left, Serrated And Flat Handle With Polished Finish, Miniature Model With Strongly Curved Blunt Blades Of Uneven Lengths, Lower Blade 12mm, Upper Blade 11mm (From Mid-Screw To Tip), And Overall Length Of 4" (10cm)

VS11-5-3141

Katzin Corneal Transplant Scissor
On 3 Hole Handle Medium Blades Curved Left

VS11-5-3152

Rosenwasser Scissor
On 3 Hole Handle Medium Curved Blades

VS11-5-3153

Rosenwasser Scissor
On 3 Hole Handle Strong Curved Blades

VS5-4338

Azar Curved Calibrated Corneal Scissors

Flat Scalloped Handle With Dull Finish, Universal Medium Model With Blunt Blades, Calibration Markings Of 1mm To 6mm, 16mm From Mid-Screw To Tip, And Overall Length Of 4 1/2" (114mm)

Spring Handle Scissors

Ophthalmology



VS5-4331

Castroviejo Curved Corneal Scissors

Serrated Flat Handle With Polished Finish, Universal Medium Model With Blunt Tips, 16mm From Mid-Screw To Tip, And Overall Length Of 4 1/2" (115mm)



VS11-5-3311

Castroviejo Corneal Scissors Curved With Small Blunt Jaws

Three-Holed Handle, 10mm Mid Screw To Tip Length, And An Overall Length Of 4 3/8" (111mm)

VS5-4331-WH

Castroviejo Keratoplasty And Keratectomy Scissors

Serrated And Wide Handle With Polished Finish, Universal Miniature Model With Medium Curved Blades, 17mm From Mid-Screw To Tip, And Overall Length Of 5 1/4" (133mm)



VS11-5-3331

Castroviejo Corneal Scissors Curved

Blunt, Medium, With Mid Screw To Tip Length Of 10mm, And An Overall Length Of 4 3/8" (11cm)



VS11-5-3351

Castroviejo Corneal Scissors Large Blades With Curved And Blunt Tips

Three-Holed Handle, Mid Screw To Tip Length Of 15mm, And An Overall Length Of 4 1/2" (11.5cm)

VS11-5-3371

Castroviejo Corneal Scissors Large Blades With Strong Curve And Blunt Tips

Three-Holed Handle, Mid Screw To Tip Length Of 15mm, And An Overall Length Of 4 1/2" (11.5cm)



VS5-4334**Castroviejo Angled Corneal Scissors**

Serrated And Flat Handle With Polished Finish, Universal Miniature Model With Blunt Blades, 19mm From Midscrew To Tip, And An Overall Length Of 4 1/4" (106mm)

VS11-5-3961**Castro Keratoplasty Scissor**

Angled Blade To Side On A Three Hole Handle With Medium Blunt Blades

VS5-4332**Castroviejo Kertoplasty Scissors**

Serrated And Flat Handle With Polished Finish, Universal With Blunt Angled Blades, Lower Blade 1mm Longer Than Upper Blade, 16mm From Mid-Screw To Tip, And Overall Length Of 4 1/4" (110mm)

VS5-4390**Aebli Corneal Section Scissors**

Right, Serrated And Flat Handle With Polished Finish, Angled Blades With Blunt Tips, 18mm From Mid-Screw To Tip, And Overall Length Of 4 3/4" (120mm)

VS5-4301**Aebli Corneal Section Scissors**

Left, Serrated And Flat Handle With Polished Finish, Angled Blade With Blunt Tips, 18mm From Mid-Screw To Tip, And Overall Length Of 4 3/4" (120mm)

Speculums

Ophthalmology



VS5-5217

Barraquer Eye Speculum


Wire 14mm Blades With Dull Finish, 20mm Spread, And Overall Length Of 1 1/2" (3.8cm). Overall Weight Of 0.1 oz



VS11-2-6121

Barraquer Wire Speculum

Large, V-Shaped End Angled Upward, Round Blades, Blade Length Of 16mm, Spread Of 29mm, Overall Length Of 1 1/2" (3.9cm). Overall Weight Of 0.1 oz



VS5-5217-15

Barraquer Eye Speculum

Wire 15mm Blades With Dull Finish, 18.5mm Spread, And Overall Length Of 1 1/2" (3.8cm)



VS5-5217-H

Barraquer Eye Speculum

Open 14mm Blades Of Heavy 1mm Wire With Dull Finish, 16mm Spread, And Overall Length Of 1 5/8" (4cm) Overall Weight Of 0.1 oz



VS11-2-6110

Barraquer Speculum Large

Heavy Wire, U-Shaped Blades, V-Shaped End, Blade Length Of 15mm, Spread Of 24mm, Overall Length Of 1 3/4" (4.3cm) Overall Weight Of 0.1 oz



VS5-5218

Barraquer Eye Speculum

Open 11mm Blades With Dull Finish, 19mm Spread, And Overall Length Of 1 3/8" (3.6cm) Overall Weight Of 0.1 oz



VS11-2-6111

Barraquer Speculum Small

Curved, Open Wire Blades, V-Shaped End, Blade Length Of 13mm, Spread Of 14mm, Overall Length Of 1 1/4" (3.3cm) Overall Weight Of 0.1 oz



VS5-5218-P

Barraquer Pediatric Speculum

Closed Wire 5mm Blade With An Overall Length Of 1 1/4" (30mm) Overall Weight Of 0.1 oz



VS5-5217-M

Connor Barraquer Eye Speculum

Wire With 14mm Blades, Dull Finish, 25mm Spread, And An Overall Length Of 1 1/2" (40mm) Overall Weight Of 0.1 oz



VS5-5217-MT

Connor Barraquer Temporal Approach Lid Speculum

Wire With 14mm Rounded Blades For Patient Comfort, Dull Finish, 25mm Spread, Angled, And Overall Length Of 1 1/2" (40mm) Overall Weight Of 0.1 oz



VS5-5217-K

Kratz-Barraquer Eye Speculum

Open 14mm Blades With Dull Finish, 18mm Spread, And Overall Length Of 1 5/8" (41mm) Overall Weight Of 0.1 oz




VS11-2-6161

Kratz-Barr Wire Speculum Large

V-Shaped End Slanted Sharply Upward, Open Wire Blades, Blade Length Of 16mm, Spread Of 15mm, Overall Length Of 1 5/8" (4.2cm) Overall Weight Of 0.1 oz

Speculums

Ophthalmology



VS5-5217-KH

Kratz-Barraquer Eye Speculum

Open 14mm Blades With Dull Finish, Heavy 1mm Wire, 18mm Spread, And Overall Length Of 1 5/8" (41mm) Overall Weight Of 0.1 oz



VS11-2-6150

Kratz-Barraquer Speculum


Large, Heavy Wire, Straight, V-Shaped End, Blade Length Of 16mm, Spread Of 19mm, Overall Length Of 1 1/2" (3.7cm) Overall Weight Of 0.1 oz



VS5-5218-K

Kratz-Barraquer Eye Speculum

Open 11mm Blades With Dull Finish, 1gmm Spread, And Overall Length Of 1 3/8" (3.6cm), Overall Weight Of 0.1 oz



VS11-2-6162

Kratz-Barraquer Speculum

Small, Straight, V-Shaped End, Blade Length Of 12mm, Spread Of 16mm, Overall Length Of 1 1/4" (3.1 cm), Overall Weight Of 0.1 oz

VS5-5218-PK

Kratz Barraquer Pediatric Speculum

Open Wire 5mm Blade With Wing Tip Twist, Spread Of 16mm, And Overall Length Of 1" (2.6cm), Overall Weight Of 0.1 oz

VS5-5217-SQT

Temporal Barraquer Eye Speculum

Square Wire 14mm Blades With Dull Finish, 25mm Spread, And Overall Length Of 1 1/2" (3.8cm) Same As 5-5217 But Temporal



VS5-5217-T

Temporal Approach Lid Specul

Wire 14mm Rounded Blades For Patient Comfort With Dull Finish, 16mm Spread, Angled, And Overall Length Of 1 1/2" (3.8cm), Overall Weight Of 0.1 oz



VS11-2-6123

Temporal Closed Wire Speculum

U-Shaped Blades, Blade Length Of 16mm, Spread Of 18mm, Overall Length Of 1 1/2" (3.8cm), Overall Weight Of 0.1 oz



VS11-2-6123-W

Temporal Closed Wire Speculum

Wide Model, U-Shaped Blades, Blade Length Of 21mm, Spread Of 40mm, Overall Length Of 2 1/4" (5.6cm)

VS5-5216

Lehner Temporal Speculum

Open 14mm Blades With 18mm Spread, And Angled 45 Degrees Upward For Temporal Approach, Overall Weight Of 0.1 oz



VS11-2-6126

Temporal Speculum With C-Shaped End

Open Wire Blades, Blade Length Of 16mm, Spread Of 21mm, Overall Length Of 1 5/8" (4cm), Overall Weight Of 0.1 oz



VS5-5217-KT

Kratz Barraquer Temporal Approach Lid Speculum

Open 14mm Rounded Blades With Dull Finish, 18mm Spread, Angled, And Overall Length Of 1 5/8" (41 mm), Overall Weight Of 0.1 oz



VS11-2-6126-R

Temporal Open Wire Speculum

Round Blades, V-Shaped End, Blade Length Of 15mm, Spread Of 17mm, Overall Length Of 1 3/4" (4.4cm), Overall Weight Of 0.1 oz



VS5-5218-KT

Temporal Speculum

Kratz Type, Blade Length Of 12mm, Spread Of 15mm, And An Overall Length Of 1 1/2" (3.7cm), Overall Weight Of 0.1 oz



VS5-5218-KTR

Temporal Speculum

Kratz Type, Round Blade, V-Shaped End, Blade Length Of 11 mm, Spread Of 16mm, Overall Length Of 1 5/8" (4cm), Overall Weight Of 0.1 oz



VS5-5217-N

Kratz–Nichamin Speculum

Open Wire Style 16mm Blade With Wing Tip Twist, Spread Of 24mm, And Overall Length Of 1 5/8" (4cm)



VS5-5217-S

Barraquer Eye Speculum

Solid 14mm Blades With Dull Finish, 19mm Spread, And Overall Length Of 1 5/8" (41mm), Overall Weight Of 0.1 oz

VS5-5218-PS

Barraquer Pediatric Speculum

Solid 6mm Blade Speculum For Newborn Or Pediatric Use, And Overall Length Of 1 1/4" (30mm)



VS11-2-6143

Feast–Barr Speculum With Wide

Solid Blades, V–Shaped End, Blade Length Of 19mm, Spread Of 20mm, Overall Length Of 1 3/4" (4.3cm), Overall Weight Of 0.1 oz

VS11-2-6125

Temporal Approach Lid Speculum

Solid 14mm Blades With Dull Finish, 16mm Spread, Angled Up With U Shape Wire Back, And Overall Length Of 1 5/8" (4.1cm)

VS11-2-6261

Kratz–Barraquer Wire Speculum Large

U–Shaped End With Open Wire Blades, Blade Length Of 16mm, Spread Of 15mm, Overall Length Of 1 5/8" (4.2cm) With Silicone Tubing And Luer Lock Adaptor



VS11-2-6131

Small Solid Blade Speculum With V–Shaped End

Blade Length Of 12mm, Spread Of 20mm, Overall Length Of 1 3/8" (3.6cm), Overall Weight Of 0.1 oz

VS11-2-6155

Kratz–Barraquer Wire Eye Speculum Right

With One Solid 14mm Blade Right Side and Left Side With Open Wire Blades, Blade Length Of 16mm Dull Finish, 15m Spread, And Overall Length Of 1 5/8" (41mm)

VS11-2-6156

Kratz–Barraquer Wire Eye Speculum Left

With One Solid 14mm Blade Left Side and Right Side With Open Wire Blades, Blade Length Of 16mm Dull Finish, 15m Spread, And Overall Length Of 1 5/8" (41mm)



VS5-5217-ST

Temporal Approach Lid Speculum


Solid 14mm Blades With Dull Finish, 19mm Spread, Angled, And Overall Length Of 1 5/8" (4.1cm), Overall Weight Of 0.1 oz



VS11-2-6311

Tenant Solid Lid Retractor

Sold In A Pair, Blade Length Is 14mm, Overall Length Is 1" (2.5cm), Overall Weight Of 0.1 oz



VS11-2-6181

McIntyre Closed Wire Eye Speculum

V Shaped Blades, Blade Length Of 14mm, Spread Of 28mm, Overall Length Of 1 3/4" (4.5cm), Overall Weight Of 0.1 oz

VS5-5219

Floyd-Barraquer Eye Speculum

Open 14mm Blades With Dull Finish, 6mm Lifting Wings, 18mm Spread, And Overall Length Of 1 5/8" (40mm), Overall Weight Of 0.1 oz

VS5-5227

Thornton Open Wire Lid Speculum With Lifting Wings

Open 15mm Blades, Dull Finish, 16mm Lifting Wings, 24mm Spread, And Overall Length Of 1 1/2" (40mm), Overall Weight Of 0.1 oz

VS5-5225

Rappazzo Lid Speculum

Open 15mm Blades With Dull Finish, 44mm Spread, Screw Locking Mechanism, Lifting Wings, And Overall Length Of 2 1/10" (53mm), Overall Weight Of 0.1 oz

VS11-2-6171

Kratz Barraquer Open Wire Lid Speculum With Lifting Wings, Open 15mm Round Bar Blades, Dull Finish, 16mm Lifting Wings, 24mm Spread, And Overall Length Of 1 1/2" (40mm)



VS11-2-6182

McIntyre Open Wire Speculum

V Shaped Blades, Blade Length Of 16mm, Spread Of 27mm, Overall Length Of 1 3/4" (4.5cm), Overall Weight Of 0.1 oz

Speculums

Ophthalmology



VS11-2-6184

Nevyas Speculum With V-Shaped End

Round Open Wire Blades, Blade Length Of 16mm, Spread Of 31mm, Overall Length Of 1 7/8" (4.7cm), Overall Weight Of 0.1 oz



VS11-2-6184-U

Nevyas Wire Speculum

U-Shaped End, U-Shaped Blades, Blade Length Of 22mm, Spread Length Of 40mm, Overall Length Of 2 3/8" (6cm), Overall Weight Of 0.1 oz



VS11-2-6461

Infant Wire Speculum

U-Shaped End Angled Slightly Upward, Round Blades, Blade Length Of 4mm, Spread Of 20mm, Overall Length Of 1 1/2" (3.8cm), Overall Weight Of 0.1 oz



VS5-5210

Pierse Eye Speculum

Open 11mm Blades With Polished Finish, 20mm Spread, And Overall Length Of 1 1/4" (30mm), Overall Weight Of 0.1 oz



VS5-5223

Alfonso Eye Speculum

Open 5mm Blades With Dull Finish, 27mm Spread, And Overall Length Of 1 3/4" (45mm), Overall Weight Of 0.1 oz



VS5-54918S

Premature Infant Self Retaining Retractor Small

Overall Length 1 1/2" (3.8cm) With Round Wire Blades, Blade Length Of 4mm

VS5-54919S

Premature Infant Self Retaining Retractor Small
Overall Length 1 1/2" (3.8cm) With Round Wire Blades, Blade Length Of 5mm

VS5-54921S

Infant Self Retaining Retractor Small
Overall Length 1 1/2" (3.8cm) With Round Wire Blades, Blade Length Of 7mm

VS5-54923S

Infant Self Retaining Retractor Small
Overall Length 1 1/2" (3.8cm) With Round Wire Blades, Blade Length Of 9mm

VS5-54924S

Child Self Retaining Retractor Small
Overall Length 1 1/2" (3.8cm) With Round Wire Blades, Blade Length Of 11mm

VS5-54920S

Premature Infant Self Retaining Retractor Small
Overall Length 1 1/2" (3.8cm) With Round Wire Blades, Blade Length Of 6mm

VS5-54922S

Infant Self Retaining Retractor Small
Overall Length 1 1/2" (3.8cm) With Round Wire Blades, Blade Length Of 8mm

VS5-54924S

Child Self Retaining Retractor Small
Overall Length 1 1/2" (3.8cm) With Round Wire Blades, Blade Length Of 10mm

VS5-54924S

Child Self Retaining Retractor Small
Overall Length 1 1/2" (3.8cm) With Round Wire Blades, Blade Length Of 12mm



VS5-5221

Sauer Infant (A) Eye Speculum
Solid 8mm Blades With Dull Finish, 20mm Spread, Overall Length Of 1 3/8" (35mm)

VS11-2-6413

Sauer Premature Infant Lid Speculum
Tiny Small Solid Blades, Blade Length Of Blank mm, Spread Of 20mm, Overall Length Of 1 1/8" (2.8cm)



VS5-5224

Precision Maloney Style Adjustable Lid Wire Speculum
Straight, Blade Length Of 15mm, Spread Of 51 mm, Overall Length Of 3 1/4" (8.4cm), Overall Weight Of 0.1 oz



VS11-2-6411

Sauer Infant (B) Eye Speculum
Solid 11mm Blades With Dull Finish, 18mm Spread, Overall Length Of 1 3/8" (34mm)

VS11-2-6421

Frankel Infant Adjustable Twist Lid Speculum
Small Solid Blades, Blade Length Of Blank mm, Spread Of Blank mm, Overall Length Of Blank" (Blankcm)



VS5-5176

Lieberman Open Wire Speculum
Open 17mm Blades With Adjustable Wire, Adjustable Lid Tension With Thumb Screw Control, And Available With Straight Shaft Or Angled 45 Degrees, 50mm Spread And Overall Length Of 3 3/8" (8.4cm), Overall Weight Of 0.1 oz

Speculums

Ophthalmology



VS5-5176-A

Lieberman Open Wire Speculum

Open 17mm Blades With Adjustable Wire, Adjustable Lid Tension With Thumb Screw Control, And Angled 45 Degrees, Spread Of 37mm, Overall Length Of 3 1/8" (8.1cm), Overall Weight Of 0.1 oz



VS5-5176-AT

Lieberman Open Wire Temporal Speculum

Open 16mm Blades With Adjustable Wire, Adjustable Lid Tension With Thumb Screw Control, Angled 45 Degrees, In Titanium, With An Overall Length Of 3 1/4" (8.4cm), Overall Weight Of 0.1 oz



VS11-2-6782

Lieberman Adult K-Wire Speculum With Adjustable Mechanism

Kratz-Style Open Wire Blades Angled Downward Towards The Temple, Blade Length Of 15mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-R

Lieberman Adult Round-Wire Speculum With Adjustable Mechanism

Kratz Round Style Open Wire Blades Angled Downward Towards The Temple, Blade Length Of 15mm Overall Length Of Speculum Is 3" (7.6cm)



VS11-2-6782-LS

Special Order Item *

Frankel Infant Adjustable Twist Lid Speculum

Small Solid Blades, Blade Length Of Blank mm, Spread Of Blank mm, Overall Length Of Blank" (Blankcm)

VS11-2-6782-R

Lieberman Adult Round-Wire Speculum With Adjustable Mechanism

Kratz Round Style Open Wire Blades Angled Downward Towards The Temple, Blade Length Of 15mm Overall Length Of Speculum Is 3" (7.6cm)

***Special Order Items are Non-Cancelable**

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS11-2-6782-RS

Special Order Item*

Lieberman Adult K-Wire And Solid Blade Speculum With Adjustable Mechanism

Kratz-Style Open Wire Left Blade And Solid Adult Right Blade Angled Downward Towards Towards The Temple, Both Blade 15mm, Overall Length Of Speculum Is 3" (7.6cm)



VS11-2-6782-C

Lieberman Child K-Wire Speculum With Adjustable Mechanism

Kratz-Style Open Wire Blades Angled Downward Towards The Temple, Blade Length Of 11mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-LSC

Special Order Item*

Lieberman Child K-Wire And Solid Blade Speculum With Adjustable Mechanism

Kratz-Style Open Wire Right Blade And Solid Adult Left Blade Angled Downward Towards Towards The Temple, Both Blade 11mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-RSC

Special Order Item*

Lieberman Child K-Wire And Solid Blade Speculum With Adjustable Mechanism

Kratz-Style Open Wire Left Blade And Solid Adult Right Blade Angled Downward Towards Towards The Temple, Both Blade 11mm, Overall Length Of Speculum Is 3" (7.6cm)



VS11-2-6782-N

Lieberman Nasal Adult K-Wire Speculum With Adjustable Mechanism

Kratz-Style Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 15mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-NLS

Special Order Item*

Lieberman Nasal Adult K-Wire And Solid Blade Speculum With Adjustable Mechanism

Kratz-Style Open Wire Right Blade And Solid Adult Left Blade Angled Upward For The Bridge Of The Nose, Blade Length Of 15mm, Overall Length Of Speculum Is 3" (7.6cm)

***Special Order Items are Non-Cancelable**



VS11-2-6782-NRS

Special Order Item*

Lieberman Nasal Child K-Wire And Solid Blade Speculum With Adjustable Mechanism
Kratz-Style Open Wire Left Blade And Solid Adult Right Blade Angled Upward For The Bridge Of The Nose, Blade Length Of 11mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-NLSC

Special Order Item*

Lieberman Nasal Child K-Wire And Solid Blade Speculum With Adjustable Mechanism
Kratz-Style Open Wire Right Blade And Solid Adult Left Blade Angled Upward For The Bridge Of The Nose, Blade Length Of 11mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-P

Lieberman Pediatric K-Wire Speculum With Adjustable Mechanism
Kratz-Style Open Wire Blades Angled Downward Towards The Temple, Blade Length Of 5mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-RC

Lieberman Child Round-Wire Speculum With Adjustable Mechanism
Kratz-Style Open Wire Blades Angled Downward Towards The Temple, Blade Length Of 11mm Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-NRC

Lieberman Nasal Child Round-Wire Nasal Speculum With Adjustable Mechanism
Round Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 11mm Overall Length Of Speculum Is 3" (7.6cm)



VS11-2-6786

Lieberman Adult V-Wire Speculum With Adjustable Mechanism
V-Shaped Style Open Wire Blades Angled Downward Towards Towards The Temple, Blade Length Of 14mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-NC

Lieberman Nasal Child K-Wire Speculum With Adjustable Mechanism
Kratz-Style Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 11mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-NRSC

Special Order Item*

Lieberman Nasal Child K-Wire And Solid Blade Speculum With Adjustable Mechanism
Kratz-Style Open Wire Left Blade And Solid Adult Right Blade Angled Upward For The Bridge Of The Nose, Blade Length Of 11mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-R

Lieberman Adult Round-Wire Speculum With Adjustable Mechanism
Kratz-Style Open Wire Blades Angled Downward Towards Towards The Temple, Blade Length Of 15mm Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-NR

Lieberman Nasal Adult Round-Wire Nasal Speculum With Adjustable Mechanism
Round Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 15mm Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6782-T

Titanium Lieberman Adult K-Wire Speculum With Adjustable Mechanism
Kratz-Style Open Wire Blades Angled Downward Towards Towards The Temple, Blade Length Of 15mm, Overall Length Of Speculum Is 3" (7.6cm)



VS11-2-6786-LS

Special Order Item*

Lieberman Adult V-Wire And Solid Blade Speculum With Adjustable Mechanism
V-Shaped Style Open Wire Right Blade And Solid Child Left Blade Angled Downward Towards Towards The Temple, Both Blade 14mm, Overall Length Of Speculum Is 3" (7.6cm)

***Special Order Items are Non-Cancelable**



VS11-2-6786-RS

Special Order Item *

Lieberman Adult V-Wire And Solid Blade Speculum With Adjustable Mechanism
V-Shaped Style Open Wire Left Blade And Solid Child Right Blade Angled Downward Towards Towards The Temple, Both Blade 14mm, Overall Length Of Speculum Is 3" (7.6cm)



VS11-2-6786-C

Lieberman Child V-Wire Speculum With Adjustable Mechanism

V-Shaped Style Open Wire Blades Angled Downward Towards Towards The Temple, Blade Length Of 10mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6786-LSC

Special Order Item *

Lieberman Child V-Wire And Solid Blade Speculum With Adjustable Mechanism
V-Shaped Style Open Wire Right Blade And Solid Child Left Blade Angled Downward Towards Towards The Temple, Both Blade 11mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6786-RSC

Special Order Item *

Lieberman Child V-Wire And Solid Blade Speculum With Adjustable Mechanism
V-Shaped Style Open Wire Left Blade And Solid Child Right Blade Angled Downward Towards Towards The Temple, Both Blade 11mm, Overall Length Of Speculum Is 3" (7.6cm)



VS11-2-6786-N

Lieberman Nasal Wire Speculum
Open Wire V-Shaped Blades Angled Slightly Upward, Blade Length Of 15mm, Spread Of 44mm, Overall Length Of 3 1/8" (8.1cm), Overall Weight Of 0.1 oz

VS11-2-6786-NLS

Special Order Item *

Lieberman Nasal Adult V-Wire And Solid Blade Speculum With Adjustable Mechanism
V-Shaped Style Open Wire Right Blade And Solid Adult Left Blade Angled Upward For The Bridge Of The Nose, Blade Length Of 14mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6786-NRS

Special Order Item *

Lieberman Nasal Adult V-Wire And Solid Blade Speculum With Adjustable Mechanism
V-Shaped Style Open Wire Left Blade And Solid Adult Right Blade Angled Upward For The Bridge Of The Nose, Blade Length Of 14mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6786-NC

Lieberman Nasal Child V-Wire Speculum With Adjustable Mechanism
V-Shaped Style Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 10mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6786-NLSC

Lieberman Nasal Child V-Wire And Solid Blade Speculum With Adjustable Mechanism
V-Shaped Style Open Wire Left Blade And Solid Adult Right Blade Angled Upward For The Bridge Of The Nose, Blade Length Of 10mm, Overall Length Of Speculum Is 3" (7.6cm)


VS11-2-6786-NRSC

Lieberman Nasal Child V-Wire And Solid Blade Speculum With Adjustable Mechanism
V-Shaped Style Open Wire Right Blade And Solid Adult Left Blade Angled Upward For The Bridge Of The Nose, Blade Length Of 10mm, Overall Length Of Speculum Is 3" (7.6cm)

***Special Order Items are Non-Cancelable**

Speculums

Ophthalmology



VS11-2-6282

Lieberman Aspirating Adult K-Wire Speculum With Adjustable Mechanism

Kratz-Style Open Wire Blades Angled Downward Towards Towards The Temple, Blade Length Of 15mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing And Cannula (Length Of 1 7/8" 3.8cm)

VS11-2-6282-N

Lieberman Nasal Aspirating Adult K-Wire Speculum With Adjustable Mechanism

Kratz-Style Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 15mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing And Cannula (Length Of 1 7/8" 3.8cm)

VS11-2-6282-NR

Lieberman Nasal Aspirating Adult Round-Wire Nasal Speculum With Adjustable Mechanism

Round Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 15mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing

VS11-2-6286-N

Lieberman Nasal Aspirating Child V-Wire Speculum With Adjustable Mechanism

V-Shaped Style Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 14mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing And Cannula (Length Of 1 7/8" 3.8cm)



VS11-2-6282-R

Lieberman Aspirating Adult Round-Wire Speculum With Adjustable Mechanism

Kratz-Style Open Wire Blades Angled Downward Towards Towards The Temple, Blade Length Of 15mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing And Cannula (Length Of 1 7/8" 3.8cm)



VS11-2-6282-C

Lieberman Aspirating Child K-Wire Speculum With Adjustable Mechanism

Kratz-Style Open Wire Blades Angled Downward Towards The Temple, Blade Length Of 11mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing And Cannula (Length Of 1 7/8" 3.8cm)

VS11-2-6282-NC

Lieberman Nasal Aspirating Child K-Wire Speculum With Adjustable Mechanism

Kratz-Style Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 11mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing And Cannula (Length Of 1 7/8" 3.8cm)

VS11-2-6282-NRC

Lieberman Nasal Aspirating Child Round-Wire Nasal Speculum With Adjustable Mechanism

Round Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 11mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing

VS11-2-6286-NC

Lieberman Nasal Aspirating Adult V-Wire Speculum With Adjustable Mechanism


V-Shaped Style Open Wire Blades Angled Upward For The Bridge Of The Nose, Blade Length Of 11mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing And Cannula (Length Of 1 7/8" 3.8cm)

VS11-2-6282-RC

Lieberman Aspirating Child Round-Wire Speculum With Adjustable Mechanism

Kratz-Style Open Wire Blades Angled Downward Towards The Temple, Blade Length Of 11mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing And Cannula (Length Of 1 7/8" 3.8cm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS11-2-6286

Lieberman Aspirating Adult V-Wire Speculum With Adjustable Mechanism

V-Shaped Style Open Wire Blades Angled Downward Towards Towards The Temple, Blade Length Of 14mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing And Cannula (Length Of 1 7/8" 3.8cm)

VS11-2-6286-C

Lieberman Aspirating Child V-Wire Speculum With Adjustable Mechanism

V-Shaped Style Open Wire Blades Angled Downward Towards The Temple, Blade Length Of 11mm, Spread Of 31mm, Overall Length Of Speculum Is 3" (7.6cm), Includes Silicone Tubing And Cannula (Length Of 1 7/8" 3.8cm)



VS11-2-6780

Ginsberg Adjustable Wire Eye Speculum

Blades Half Solid, Half Open-Wire, Angled Upward, Blade Length Of 15mm, Spread Of 47mm, Overall Length Of 3 1/8" (8.1cm)

VS11-2-6788

Lieberman Child Solid Blade Speculum With Adjustable Mechanism

Angled Downward Towards The Temple, Blade Length Of 10mm, Overall Length Of Speculum Is 3" (7.6cm)



VS11-2-6789

Lieberman Adult Solid Blade Speculum With Adjustable Mechanism

Angled Downward Towards The Temple, Blade Length Of 14mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6789-N

Lieberman Adult Solid Blade Speculum With Adjustable Mechanism

Angled Upward Towards The Bridge Of The Nose, Blade Length Of 14mm, Overall Length Of Speculum Is 3" (7.6cm)

Speculums

Ophthalmology



VS11-2-6791

Lieberman Reversible Adult Solid Blade Speculum With Adjustable Mechanism
Can Be Used For Both Temporal Or Nasal Approach, Blade Length Of 14mm, Overall Length Of Speculum Is 3" (7.6cm)

VS11-2-6791-C

Lieberman Reversible Child Solid Blade Speculum With Adjustable Mechanism
Can Be Used For Both Temporal Or Nasal Approach, Blade Length Of 10mm, Overall Length Of Speculum Is 3" (7.6cm)



VSTDK-2036-ES

Schapira Twist Speculum Titanium And Curved With Barraquer Style
13mm Blades, And An Overall Length Of 3 3/8" (8.5cm)



VS5-5222

Saunders Eye Speculum
Solid 5mm Blades With Dull Finish, 15mm Spread, Locking Thumb Screw, Holes In Bar, And Overall Length Of 1 1/8" (3cm), Overall Weight Of 0.1 oz



VS5-5193

Cook Eye Speculum
Pediatric Model With Dull Finish, Solid 11mm Blades, 26mm Spread, Locking-Mechanism, And Overall Length Of 1 1/8" (4cm), Overall Weight Of 0.1 oz



VS11-2-6512

Cook Speculum
Small #1, Lock Screw, Round Solid Blades Curved Slightly Upward, Blade Length Of 8mm, Spread Of 22mm, Overall Length Of 1 1/2" (3.8cm), Overall Weight Of 0.1 oz



VS11-2-6513

Cook Speculum

Medium #2, Lock Screw, Round Solid Blades, Blade Length Of 11mm, Spread Of 26mm, Overall Length Of 1 3/4" (4.3cm), Overall Weight Of 0.1 oz



VS11-2-6514

Cook Speculum

Large With #3 Lock Screw, Round Solid Blades Curved Upward, Blade Length Of 16mm, Spread Of 32mm, Overall Length Of 1 3/4" (4.5cm), Overall Weight Of 0.1 oz



VS5-5197

McKinney Eye Speculum

Solid 15mm Blades With Dull Finish, 31mm Spread, Self-Locking Mechanism, And Overall Length Of 2 3/8" (60mm), Overall Weight Of 0.1 oz



VS5-5105

Murdoch Eye Speculum

Fenestrated 16mm Blades With Dull Finish, 25mm Spread, Self-Locking Mechanism, And Overall Length Of 2 3/8" (6cm), Overall Weight Of 0.1 oz



VS11-2-6517

Murdoch Eye Speculum

#2 With Lock Screw, Round Open-Wire Blades Curved Slightly Upward, Blade Length Of 15mm, Spread Of 26mm, Overall Length Of 1 3/4" (4.6cm), Overall Weight Of 0.1 oz



VS11-2-6518

Murdoch Eye Speculum

#3 With Lock Screw, Open-Wire Blades Curved Slightly Upward, Blade Length Of 17mm, Spread Of 27mm, Overall Length Of 1 3/4" (4.6cm), Overall Weight Of 0.1 oz

Speculums

Ophthalmology



VS11-2-6611

Williams Small Speculum

Open Wire Blades Curved Slightly Upward, Blade Length Of 12mm, Spread Of 30mm, Overall Length Of 2 3/8" (6.6cm), Overall Weight Of 0.1 oz



VS5-5211

Wiener Eye Speculum

Fenestrated 14mm Blades With Dull Finish, 22mm Spread, Wide Extension Of Blade For Lid Protection, Self-Locking Mechanism, And Overall Length Of 3 1/4" (8.2cm), Overall Weight Of 0.1 oz

VS5-5215

Osher Lid Speculum

Open Wire Style Blades, And Overall Length Of 5 7/8" (150mm), Overall Weight Of 0.1 oz

VS5-5183

Mellinger Eye Speculum

Fenestrated 17mm Blades With Dull Finish, 29mm Spread, Self-Locking Mechanism, And Overall Length Of 2 5/8" (6.6cm), Overall Weight Of 0.1 oz



VS5-5183-G

Gimble-Mellinger Speculum

Fenestrated 16mm Blades With Dull Finish, 29mm Spread, Self-Locking Mechanism, And Overall Length Of 2 5/8" (6.6cm), Overall Weight Of 0.1 oz



VS5-5184

Mellinger Eye Speculum

Solid 16mm Blades With Dull Finish, 30mm Spread, Self-Locking Mechanism, And Overall Length Of 2 5/8" (6.8cm), Overall Weight Of 0.1 oz



VS5-5161

Williams Eye Speculum

Adult Model With Dull Finish, 14mm Blades With Spring-Controlled 34mm Spread And Locking Mechanism, And Overall Length Of 3" (7.7cm), Overall Weight Of 0.1 oz



VS5-5162

Williams Eye Speculum

Pediatric Model With Dull Finish, 11mm Blades With Spring-Controlled 41mm Spread And Locking Mechanism, And Overall Length Of 2 3/4" (7cm), Overall Weight Of 0.1 oz



VS5-5175

Knapp Eye Speculum

Open 13mm Blades With 60mm Spread, With Thumb Screw Locking Mechanism, Dull Finish, And Overall Length Of 3 1/4" (85mm)



VS5-5167

Lancaster Eye Speculum

Solid 16mm Blades. Screw Controlled Locking Mechanism with Stabilizing Disk. Blade Spread 38mm. Dull Finish, Overall Length 3" (7.7cm), Overall Weight Of 0.1 oz



VS5-5171

Lester-Burch Eye Speculum

Solid 21mm Blades With 32mm Spread, With Thumb Screw Locking Mechanism, Dull Finish, And Overall Length Of 3 1/8" (80mm)

VS5-5169

Lancaster K Eye Speculum

Solid 15mm Blades With 8mm Groove, 30mm Spread, With Thumb Screw Dull Finish, And Overall Length Of 3 1/4" (82mm)



VS5-5151

Castroviejo Speculum

Fenestrated 15mm Blades With With Thumb Screw 28mm Spread And Dull Finish, And Overall Length Of 3" (75mm)



VS5-5152

Castroviejo Speculum

Fenestrated 16mm Blades With With Thumb Screw 31mm Spread And Dull Finish, And Overall Length Of 3 1/4" (82mm)

Speculums

Ophthalmology



VS5-5118

Maumenee-Park Speculum

Solid 14mm Blades With Dull Finish, 36mm Spread, Center Retractor Blade With Tightening Screw, Adjustable Blade Locking Thumb Screws On Outside Of Retractor Blades And Overall Length Of 3 1/8" (80mm)

VS5-5117

Maumenee-Park Speculum

Fenestrated 14mm Blades With Dull Finish, 36mm Spread, Center Retractor Blade With Tightening Screw, Adjustable Blade Locking Thumb Screws On Outside Of Retractor Blades And Overall Length Of 3 1/8" (80mm)



VS5-5112

Guyton-Park Speculum

Solid 14mm Blades With Dull Finish, 45mm Spread, Suture Posts, Canthus Bar, And Overall Length Of 3 1/4" (82mm), Overall Weight Of 0.1 oz

VS5-5111

Guyton-Park Speculum

Fenestrated 14mm Blades With Dull Finish, 53mm Spread, Suture Posts, Canthus Bar, And Overall Length Of 3 1/4" (8.2cm), Overall Weight Of 0.1 oz



VS5-5125

McPherson Solid Blade Speculum

Blades Are 14mm Wide, 30mm Spread, Stabilizing Bar With Thumb Screw, And Overall Length Of 3 1/4" (82mm)

VS5-5119

Park Speculum

Fenestrated 15mm Blades With Dull Finish, 36mm Spread, Adjustable Blade Locking Thumb Screws On Outside Of Retractor Blades, And Overall Length Of 3 1/4" (82mm)



VS11-2-6066

Schott Eye Speculum Wire Blade
With Rotatable Adjustment

VS11-2-6062

Smirmaul Eye Speculum Wire Blade
With Adjustable Blades

VS11-2-6061

Smirmaul Eye Speculum Solid Blade
With Adjustable Blades

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS5-3512

Stahl Caliper

Dull Finish, Measures: 7.0mm, 10.0mm, 10.5mm, 11.0mm, 11.5mm, 12.0mm, 12.5mm, And 13.0mm, And Overall Length Of 1.6" (40mm)



VS11-4-0411

Stahl A/C Caliper Square

For Measuring corneal Diameter From "White to White" With Size Markings on Both Sides; Measures: 10.0mm, 10.5mm, 11.0mm, 11.5mm, 12.0mm, 12.5mm, And 13.0mm

VS11-4-0412

Stahl A/C Caliper Square

For Measuring corneal Diameter From "White to White" With Size Markings on Both Sides; Measures: 10.25mm, 10.75mm, 11.25mm, 11.75mm, 12.25mm, 12.75mm, And 13.25mm



VS5-3513

Braunstein Fixed Caliper

Double Ended Fine Points Overall Length 3 1/4" 83mm With 3.5mm One End And 4.0mm On The Other End



VS5-3514

Braunstein Fixed Caliper

Double Ended Fine Points Overall Length 3 1/4" 83mm With 3.0mm One End And 3.5mm On The Other End



VS5-3515

Castroviejo Caliper

Polished Finish, Measures From 0mm To 20mm In 1mm Increments With Scale Readings On Both Sides, Thumb Screw Hold Caliper in Place for Correct Measurement, With 3.8cm Mid Screw To Tip, And Overall Length Of 3 3/8" (8.6cm). Overall Weight 0.6 oz



VS5-3515-T

Castroviejo Caliper

Titanium, Polished Finish, Measures From 0mm To 20mm In 1mm Increments With Scale Readings On Both Sides, Thumb Screw Hold Caliper in Place for Correct Measurement, With 3.8cm Mid Screw To Tip, And Overall Length Of 3 3/8" (8.6cm). Overall Weight 0.6 oz



VS5-3516

Castroviejo Marking Caliper

Polished Finish, 45 Degree Curved Shafts With Grooves, Measures From 0mm To 20mm In 1mm Increments With Scale Readings On Both Sides, Thumbscrew Holds Tips Apart At Desired Setting, And Overall Length Of 3 1/4" (8.3cm). Overall Weight 0.7 oz

VS5-3516-D

Castroviejo Angled Caliper

Polished Finish, 45 degree Curved Shafts With Grooves, Measures From 0mm To 20mm In 1mm Increments On Both Sides, Fine Tips, Thumbscrew Holds Tips Apart, And Overall Length Of 3 3/8" (8mm)

VS5-3520


Osher Internal Caliper

Dull Finish, Measures Up To 10mm In 1mm increments, Uses Fine Pointed Tips, And Overall Length Of 2 3/4" (68mm)

VS5-3521

Jameson Caliper

Dull Finish, Measures 0mm To 80mm In 1mm Increments And 0" To 3" In 1/16" Increments, Narrow Flat 3mm Tips, Thumbscrew For Locking Caliper, And Overall Length Of 3 3/4" (9.4cm). Overall Weight 1.2 oz



VS5-3525

Thorpe Caliper


Dull Finish, Measures Up To 80mm with 1mm Increments And 3" In 1/16" Increments, Fine 3mm Wide Angulated Tips And Overall Length Of 4 3/4" (122mm).



VS5-3531

Stainless Steel Flexible Ruler

Dull Finish, Deeply Etched Graduations, Measures Up To 152mm In 1mm Increments, 6" In 1/32" Increments, And Overall Length Of 6 1/4" (16cm). Overall Weight 0.3 oz



VS5-3532

Mendez Degree Gauge

Stainless, Round And Knurled Handle with Dull finish, Calibrated Every 10 Degrees From 0 Degrees To 180 Degrees, Inner Dimensions Are 3/4" (12mm), And Overall Length Of 4 5/8" (11.7cm). Overall Weight 0.5 oz



VS11-4-8005

Lindstrom Acruate Marker

With Smooth And Flat Handle, 6mm Marker, And An Overall Length Of 4" (10cm)

VS11-4-8007

Lindstrom Acruate Marker

With Smooth And Flat Handle, 7mm Marker, And An Overall Length Of 4" (10cm)

VS11-4-8011

Mendez Degree Gauge

Stainless, With Smooth And Flat Handle Dull finish, Calibrated Every 10 Degrees From 0 Degrees To 180 Degrees, Inner Dimensions Are 3/4" (12mm), And Overall Length Of 4 5/8" (11.7cm) Works With 11-4-8021



VS11-4-8021

Bores Axis Marker

With A Smooth And Flat Handle, 12mm Long Marker, And An Overall Length Of 4" (10cm) Works With 11-4-8011



VS5-3532

Mendez Degree Gauge

Stainless, Round And Knurled Handle with Dull finish, Calibrated Every 10 Degrees From 0 Degrees To 180 Degrees, Inner Dimensions Are 3/4" (12mm), And Overall Length Of 4 5/8" (11.7cm)



VS5-0111

Thornton Corneal Press On Ruler For Limbal Incisions

Curved Marking Blade With 2mm Lines 1mm Apart Measuring Up To 10mm, Round And Knurled Handle With Dull Finish, And Overall Length of 4 1/4" (11cm)



VS11-4-8053

Grandon T Incision Marker

With Round And Knurled Handle, 2.5 Incision Marker, And An Overall Length Of 4 1/4" (11cm)



VS11-4-8054

Grandon T Incision Marker

With A Smooth And Flat Handle, 3.0mm Marker Blade, And An Overall Length Of 4 1/4" (11cm)



VS11-4-8831

Mendez Corneal Marker

With Smooth And Flat Handle, 4mm Inner Diameter Marker, And An Overall Length Of 4 1/4" (11cm)



VS11-4-8841

Chayet Lasik Corneal Marker

With A Smooth And Flat Handle, With 3mm Inner Diameter And 9mm Outer Diameter Marker, And An Overall Length Of 4 1/4" (11cm)



VS11-4-8846

LU Lasik Corneal Marker

With Smooth And Flat Handle, 2.5mm Inner Diameter And 11.5mm Outer Diameter Marker, and An Overall Length Of 4 1/4" (11cm)



VS11-4-8848


Dulaney Lasik Marker

With A Smooth And Flat Handle, 5mm Long Marker, And An Overall Length Of 4" (10cm)

VS11-4-8903

Mandel Keratometer

With 3.5mm Ring Used To Project A Circular Shadow On The Cornea For Quantitative Evaluation Of Surgically Induced Astigmatism



VS11-4-8999

Maloney Keratometer

Titanium With A Bottom Diameter Of 30mm, Upper Diameter Of 10mm, And An Overall Length Of 1" (2.5cm)



VS11-4-9212

Greene Visual Axis Marker

With Smooth And Flat Handle, 1.5mm Bend At Tip, And An Overall Length Of 4 1/2" (11.5cm)

VS11-4-9228

Precision Mcdonald Style Laser Optical Zone Marker

With Cross Hairs Double Ended Size 6.50mm and 7.00mm



Bores O-Z Marker

With Smooth Handle, 2.00mm to 9.5mm Diameter Marker, And An Overall Length Of 4" (10cm)

VS11-4-9237	2.00mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9238	2.25mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9239	2.50mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9230	2.75mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9241	3.00mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9242	3.25mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9243	3.50mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9244	3.75mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9245	4.00mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9246	4.25mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9247	4.50mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9248	4.75mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9249	5.00mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9240	5.25mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9251	5.50mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9252	5.75mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9253	6.00mm Diameter Marker	Overall Length Of 4 1/4" (11 cm)
VS11-4-9255	6.50mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9257	7.00mm Diameter Marker	Overall Length Of 4 1/4" (11 cm)
VS11-4-9259	7.50mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9261	8.00mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9263	8.50mm Diameter Marker	Overall Length Of 4" (10cm)
VS11-4-9267	9.50mm Diameter Marker	Overall Length Of 4 1/4" (11 cm)



Berkeley Marker

With Smooth And Flat Handle, 3.00mm And 3.25mm to 5.0mm And 6.00mm Diameter Marker, And An Overall Length Of 4 3/8" (11 cm)

VS11-4-9281	3.0mm And 3.25mm Diameter Marker
VS11-4-9283	3.5mm And 3.75mm Diameter Marker
VS11-4-9285	4.0mm And 4.25mm Diameter Marker
VS11-4-9286	4.0mm And 4.50mm Diameter Marker
VS11-4-9287	4.5mm And 5.00mm Diameter Marker
VS11-4-9289	5.0mm And 6.00mm Diameter Marker



Thornton Low Profile Optical Zone Marker

Round And Knurled Handle With Dull Finish, Size 3.00mm to 11 mm , 1.5mm Ring Height Beveled To Outside, And Overall Length Of 3 1/8" (8cm)

VS5-0141-3.00	Size 3.00mm
VS5-0141-3.25	Size 3.25mm
VS5-0141-3.50	Size 3.50mm
VS5-0141-3.75	Size 3.75mm
VS5-0141-4.00	Size 4.00mm
VS5-0141-4.25	Size 4.25mm
VS5-0141-4.50	Size 4.50mm
VS5-0141-4.75	Size 4.75mm
VS5-0141-5.00	Size 5.00mm
VS5-0141-5.25	Size 5.25mm
VS5-0141-5.50	Size 5.50mm
VS5-0141-5.75	Size 5.75mm
VS5-0141-6.00	Size 6.00mm
VS5-0141-7.00	Size 7.00mm
VS5-0141-8.00	Size 8.00mm
VS5-0141-9.00	Size 9.00mm
VS5-0141-10.00	Size 10.00mm
VS5-0141-11.00	Size 11.00mm

VS5-0142

Neumann 6mm X 9mm Double Optical Zone Marker

Concentric 6mm And 9mm Rings, Ring Height Is 3mm, Round And Knurled Handle With Dull Finish, And Overall Length Of 3 1/4"(82mm)



VS5-1656

Adson Style Bipolar Forcep

Wide Serrated Handle Overall Length 4 3/4" (120.65mm) With A (High Polished 10mm Long Area) With A Tip Size Of 1.5mm.
Overall Weight 0.8 oz

VS11-4-9441	3.00mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9451	3.25mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9461	3.50mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9471	3.75mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9481	4.00mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9491	4.25mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9501	4.50mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9511	4.75mm Diameter	Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9521	5.00mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9523	5.50mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9525	6.00mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9527	6.50mm Diameter	Cross Hairs	Overall Length Of 4" (10cm)
VS11-4-9529	7.00mm Diameter	Extra Fine Cross Hairs	Overall Length Of 4 1/4" (11cm)
VS11-4-9531	7.50mm Diameter	Cross Hairs	Overall Length Of 4 1/4" (11cm)
VS11-4-9533	8.00mm Diameter	Cross Hairs	Overall Length Of 4 1/4" (11cm)
VS11-4-9537	9.00mm Diameter	Cross Hairs	Overall Length Of 4 1/4" (11cm)



Hoffer Optical Zone Marker

With Cross Hairs- Size 3.00mm to 9.50mm, Ring Height Is 3mm Beveled To Outside, Round And Knurled Handle With Dull Finish, And Overall Length Of 4 1/4" (11cm)

VS5-0122-3.00	With Cross Hairs- Size 3.00mm
VS5-0122-3.25	With Cross Hairs- Size 3.25mm
VS5-0122-3.50	With Cross Hairs- Size 3.50mm
VS5-0122-3.75	With Cross Hairs- Size 3.75mm
VS5-0122-4.00	With Cross Hairs- Size 4.00mm
VS5-0122-4.25	With Cross Hairs- Size 4.25mm
VS5-0122-4.50	With Cross Hairs- Size 4.50mm
VS5-0122-4.75	With Cross Hairs- Size 4.75mm
VS5-0122-5.00	With Cross Hairs- Size 5.00mm
VS5-0122-5.25	With Cross Hairs- Size 5.25mm
VS5-0122-5.50	With Cross Hairs- Size 5.50mm
VS5-0122-6.00	With Cross Hairs- Size 6.00mm
VS5-0122-7.00	With Cross Hairs- Size 7.00mm
VS5-0122-8.00	With Cross Hairs- Size 8.00mm
VS5-0122-9.00	With Cross Hairs- Size 9.00mm
VS5-0122-9.50	With Cross Hairs- Size 9.50mm

Markers

Ophthalmology



VS5-0133-4

Osher Neuman Corneal Marker

4 Lines And Marking 4mm Long Radial Lines, 4.5mm Inner Diameter Of Blades, And Overall Length Of 4 5/8" (11.6cm)

VS5-0133-6

6 Lines And Marking 4mm Long Radial Lines Overall Length Of 4 1/2" (11.5cm)

VS5-0133-8

8 Lines And Marking 4mm Long Radial Lines Overall Length Of 4 5/8" (11.6cm)

VS5-0133-10

10 Lines And Marking 4mm Long Radial Lines Overall Length Of 4 1/2" (11.5cm)

VS5-0133-12

12 Lines And Marking 4mm Long Radial Lines Overall Length Of 4 1/2" (11.5cm)

VS5-0133-16

16 Lines And Marking 4mm Long Radial Lines Overall Length Of 4 1/2" (11.5cm)



VS11-4-9636

Ellis Astigmatism Marker

Double-Ended With A Smooth And Flat Handle, 6mm Diameter And 7mm Diameter Markers, And An Overall Length Of 4 1/2" (11.5cm)



VS11-4-9961

Anis Radial Marker

With 04 Blades, Smooth And Flat Handle, And An Overall Length Of 4 1/4" (11 cm)

VS11-4-9962

With 06 Blades

VS11-4-9963

With 08 Blades

VS11-4-9964

With 10 Blades

VS11-4-9965

With 12 Blades

VS11-4-9966

With 16 Blades



VS5-0144

Pettigrove PRK Scraper And Marker

6.5mm Marker End, Highly Polished Stainless Steel Spatula Tip, And Overall Length Of 4" (10cm)

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020



VS5-1602

Graefe Strabismus Hook

Serrated Flat Handle With Polished Finish, Flattened Hook #1 size 8mm Overall Length Of 5 3/8" (135mm)

VS5-1603

Graefe Strabismus Hook

Serrated Flat Handle With Polished Finish, Flattened Hook #2 size 10mm Overall Length Of 5 1/2" (140mm)

VS5-1604

Graefe Strabismus Hook

Serrated Flat Handle With Polished Finish, Flattened Hook #3 size 10mm Overall Length Of 5 5/8" (143mm)



VS5-1675

Scobee Oblique Muscle Hook

Serrated Flat Handle With Polished Finish, Curved Shaft With Flat 4.2mm Hook, and Overall Length Of 5 1/4" (13.3cm)



VS5-1691

Culler Rectus Muscle Hook

Serrated And Flat Handle With Polished Finish, 7mm Flattened Hook, And Overall Length Of 5" (13cm)

VS5-1693

Culler Rectus Muscle Hook

Serrated And Flat Handle With Polished Finish, Flattened 10mm Hook With Angled Tip, And Overall Length Of 4 13/16" (122mm)

Muscle Hooks

Ophthalmology



VS5-1697

Jameson Muscle Hook

Serrated And Flat Handle With Polished Finish, 2mm Bulbous Tip With 6mm Flat Portion, And Overall Length Of 5 1/4" (13.4cm)



VS5-1697-S

Jameson Muscle Hook

Serrated And Flat Handle With Polished Finish, Small Bulbous Tip With Flat Sides, And Overall Length Of 5 1/4" (13.5cm)



VS5-1699

Green Strabismus Hook

Serrated And Flat Handle With Polished Finish, Flat 7mm Hook, And Overall Length Of 5 1/8" (13cm)

VS11-4-7871	Helveston Finder Hook	Large
VS11-4-7872	Helveston Finder Hook	Small
VS11-4-7889	Helveston Finder Hook	8mm
VS11-4-7891	Helveston Finder Hook	10mm
VS11-4-7893	Helveston Finder Hook	12mm



VS18228

Precision Style Haefliger Phaco Chop
Left, With An Overall Length Of 4 3/4" (12cm)



VS18227SPT

Precision Style Haefliger Phaco
Right Handed, 45 Degree Angle, With An Overall Length Of 4 5/8" (11.6cm)



VS5-1681-N

Nichamin Style Triple Chop and Castroviejo Spatula
0.5mm With Knurled Handle And Overall Length Of 5 1/2" (14cm)



VS5-1689-R

Nagahara Style Nucleus Chopper
Round And Knurled Handle, Nucleus Divider With 1.25mm Tip And Blade Offset At 60 Degrees Angle, For Right-Handed Use, And Overall Length Of 4 1/2" (11.5cm)

VS5-1689-L

Nagahara Style Nucleus Chopper-Left
Round And Knurled Handle, Nucleus Divider With 1.25mm Tip And Blade Offset At 60 Degrees Angle, For Left-Handed Use, And Overall Length Of 4 1/2" (11.5cm)



VS5-1689-LT

Nagahara Style Phaco Chop

Titanium, Round And Knurled Handle, Nucleus Divider With 1.25mm Tip And Blade Offset At 60 Degrees Angle, For Right-Handed Use, And Overall Length Of 4 1/2" (11.5cm)

VS5-1689-RT

Nagahara Style Phaco Chop

Titanium, Round And Knurled Handle, Nucleus Divider With 1.25mm Tip And Blade Offset At 60 Degrees Angle, For Left-Handed Use, And Overall Length Of 4 1/2" (11.5cm) And Blade Offset At 60 Degrees Angle, For Right-Handed Use, And Overall Length Of 4 1/2" (11.5cm)



VS11-4-3479

Chang Style Phaco Chop Right And Double-Ended

With A Round And Knurled Handle, 1mm Bent Tip And A 2mm Bent Tip, And An Overall Length Of 5 1/8" (13cm)



VS11-4-3470

Chang Style Phaco Chopper Left And Double-Ended

With Knurled Handle, 10mm Tips To Bends, And An Overall Length Of 5 1/8" (13cm)

VS11-4-3491

Precision Style Brierley Nucleus Splitter Single Ended

VS11-4-3495

Precision Style Dodick-Kammann Nucleus Chopper

Straight, 2mm Long

VS11-4-3496

Precision Style Dodick-Kammann Nucleus Chopper

Angled 2mm Long Tip

VS11-4-3401

Precision Style Nagahara Nucleus Spatula With Flat Edge



VS11-4-3403

Fukasaku Style Snapper Hook

With Serrated Handle, 10mm Tip To Bend, And An Overall Length Of 4 1/2" (11cm)



VS11-4-3404

Fukasaku Style Pupil Snapper

Small With Knurled Handle, 10mm Tip To Bend, And An Overall Length Of 4 3/8" (11cm)



VS11-4-3406

Rosen Style Phaco

Splitter 90 Degrees With Knurled Handle, 10mm Tip To Bend, And Overall Length Of 4 3/8" (11cm))



VS11-4-3407

Rosen Style Phaco

Splitter 60 Degrees With A Knurled Handle, 10mm Tip To Bend, And An Overall Length Of 4 3/8" (11cm))



VS11-4-3409

Rowen Style Phaco

Spatula Double-Ended With Knurled Handle, 10mm Tips To Bends, And An Overall Length Of 5 3/8" (13.5cm)



VS11-4-3511

Bechert Nucleus Rotator

With Round And Knurled Handle, 7mm Tip To Bend, And An Overall Length Of 4 1/4" (11cm)

VS9-18222-R

Precision Style Goldberg Single Ended Side Port Nucleus Splitter And Manipulator



VS9-18225-R

Connor Wand Ball Manipulator

Straight Shaft With An Overall Length Of 4 1/2" (11.5cm)



VS9-18225-AR

Connor Wand Ball Manipulator

Angled Shaft With An Overall Length Of 4 1/2" (11.5cm))



VS9-18225-CR

Connor Wand Ball Manipulator

Curved Shaft With An Overall Length Of 4 1/2" (11.5cm)



VS9-24246-R

Livernois Style Fixation Loop

With An Overall Length Of 4 1/8 (10.5cm)



VS9-24325-R

Shepard Style Anterior/Posterior Capsule Polisher
Right With An Overall Length Of. 4 3/4" (12cm)



VS9-24326-R

Shepard Style Anterior/Posterior Capsule Polisher
Left With An Overall Length Of. 4 3/4" (12cm)



VS9-25613-R

Nagahara Style Phaco-Chopper
Right With Serrated Handle And Overall Length Of 4 1/2" (11.5cm)



VS9-25614-L-R

Fine Nagahara Style Phaco Chop Left
With An Overall Length Of 4 1/2" (11.5cm)



VS9-25614-R-R

Fine Nagahara Style Phaco Chopper Right
With An Overall Length Of 4 7/8" (12.4cm)



VS9-25615R

Precision Nucleus Hook Phaco Chopper Single Ended

VS9-25617-R

Precision Style Nichamin Triple Chopper Double Ended

VS9-25627

Precision Style Ball Seibel Nucleus Chopper

VS9-25627-L-R

Precision Style Ball Seibel Nucleus Chopper
Left, Angled Shaft, Stainless Steel



VS9-25632-R

Precision Style Nichamin Style Vertical Nucleus Chopper
Ambidextrous, Stainless, With An Overall Length Of 4 7/8" (12.4cm)

VS9-25640-R

Precision Style Vergus Style Phaco Chopper

VS9-25641-LR

Precision Style Nichamin Style Quick Chopper
Left

VS9-25641-R

Precision Style Nichamin Style Quick Chopper
Right

VS9-25643-R

Precision Style Nichamin Style Triple Chopper + Nichamin Quick Chopper
Right Hand Dominant Stainless Steel

VS9-25646-LR

Precision Style Seibel Style Vertical Quick
Left

VS9-25647-RR

Seibel Style Nucleus Chopper & Safety Quick Chopper
Right Hand Dominant

VS9-25648-L-R

Fine Nagahara Style & Nichamin Style I Quick Double Chopper
Left Hand Dominant Stainless Steel

VS9-25648-R-R

Fine Nagahara Style & Nichamin Style I Quick Double Chopper
Right Hand Dominant Stainless Steel

VS16-5151-L

Seibel Style Nucleus Chopper
Titanium Left

VS16-5151-R

Seibel Style Nucleus Chopper
Titanium Right

VS16-5154-R

Nichamin Style Titanium Quick Chopper
Right Hand Dominant



VS16-5166-AR

Connor Wand Style Chopper
Titanium Angled Steel



VS16-5166-CR

Connor Wand Style Chopper
Titanium Curved



VS16-5166-R

Connor Wand Style Chopper
Titanium Straight

VS16-5168-R

Nichamin Style Quick Chopper
Nucleus Pick Double Ended Titanium Right Hand Dominant

VS16-5176-L-R

Seibel Style Nucleus Double Ended SC Chopper And Spatula
Left Handed

VS16-5176-RR

Seibel Style Nucleus Double Ended SC Chopper And Spatula
Right Handed

Flexible Iris Retractors

Ophthalmology



- Ideal for temporary iris enlargement
- Easy insertions through 0.5mm, self-sealing incision
- Flat, easy to adjust silicone stopper
- Colored fiber for easy identification
- User friendly packaging

Single Use Iris Retractors

VS30-134-1



1 unit 5 pcs (5 retractors), Sterile

VS30-134-5

5 unit 5 pcs (25 retractors), Sterile

VS30-134-6

6 unit 5 pcs (30 retractors), Sterile



Reusable Iris Retractors

Supplied in reusable, steam autoclavable sterilization container.

VS30-254

1 unit 6pcs (6 retractors), Sterile



VS18-400 / VS18-405

Lint Free PVA Surgical Spears

Provides rapid recovery of fluids during ophthalmic and microsurgical procedures.

- Fiber free construction will not introduce lint into the surgical site
- Can be cut into custom shapes to reach recesses during procedures
- Gentle to tissue and highly absorbent
- Ideal for microsurgical procedures




VS18-410 / VS18-415

Cellulose Surgical Spears

Used during ophthalmic and microsurgical procedures to provide rapid recovery of fluids

- Made from all natural hydrocellular material
- Highly absorbent and fast wicking
- Soft, yet with enough rigidity to withstand gentle pressure
- Can be gas sterilized (unused Spears or Points)
- Cost effective



VS18-417 / VS18-417e Lasik Pure Spear / Extended
VS18-418 / VS18-418e Lasik Pure Spear / Extended

Product Number	Description	Quantity
VS18-400	Sterile PVA Surgical Spears	5/pack, 20 pkg/box
VS18-405	Sterile PVA Surgical Spears	10/pack, 20 pkg/box
VS18-410	Sterile Cellulose Spears	5/pack, 20 pkg/box
VS18-415	Sterile Cellulose Spears	10/pack, 20 pkg/box
VS18-417	Sterile Lasik Pure Spears	3/pack, 20pkg/box (60)
VS18-417E	Sterile Lasik Pure Spears – Expanded	2/pack, 20pkg/box (40)
VS18-418	Sterile Lasik Pure Spears	10/pack, 20pkg/box (200)
VS18-418E	Sterile Lasik Pure Spears – Expanded	10/pack, 10pkg/box (100)

Wipes and Spears

Ophthalmology



VS18-419

Lasik

Pure Disk on Handle

VS18-420

Corneal Light Shields

Smaller pore size blocks excessive light that can be damaging to the retina

- Hydrophilic sponge retains moisture on cornea
- Contours to the surface of the cornea, will not slip out of position
- Fiber free construction will not introduce lint into the surgical site
- Makes an excellent platform for retrieving sutures




VS18-421 / VS18-422 / VS18-423 / VS18-424

VS18-425 / VS18-426 / VS18-427

Lasik Eye Shield

This sponge is put in place after corneal flap is created to drain irrigating fluids from the orbit; eliminates fluid buildup in corneal bed.

- PVAc construction will not introduce lint into the surgical site
- Soft, pliable sponge easily conforms to the cornea
- Fast wicking sponge reduces fluid build-up in the canal bed



Product Number	Description	Quantity
VS18-419	Lasik Pure Disk On Handle (2mm x 9mm)	1/pack, 20pkg/box
VS18-420	Sterile Corneal Light Shields (1mm thick, 7mm diameter)	1/pack, 20pkg/box
VS18-421	Lasik Pure Shield With 6mm Flat Spot on Disk (8mm Diameter By 1mm Thick)	1/pack, 20pkg/box (20 units)
VS18-422	Chayet Shield (12.5mm)	1/pack, 20pkg/box
VS18-423	Lasik Pure Shield with 7mm Flat Spot on Disk (9mm Diameter By 1mm Thick)	1/pack, 20pkg/box (20 units)
VS18-424	Lasik Pure Drain with Round Hole on Oval (10.5mm) (40mm Long By 14mm Wide)	1/pack, 20pkg/box (20 units)
VS18-425	Lasik Eye Shield (10mm)	1/pack, 20pkg/box
VS18-426	Lasik Eye Shield (12mm)	1/pack, 20pkg/box
VS18-427	Lasik Pure Drain 11.5mm Round Hole in an Oval (42mm Long By 17mm Wide)	1/pack, 20pkg/box
VS18-428	Pled-Gets Sponge	10/pack, 20pkg/box

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020

VS18-430 / 18-431 / 18-435

Precision Eye Drains 80cc / 400cc

Collects and controls excess fluids during ophthalmic procedures

- Wick material creates a self-wicking drainage system which removes irrigating and body fluids from the site
- Adhesive backed collection bag adheres to the surgical drape to eliminate any interference with the procedure
- Easily trimmed to fit all patient sizes

VS18-440

Bloomberg Ophthalmic Ring





- Was designed to retain moisture on the sclera
- Has a unique small pore construction which allows it to contour to the surface of the eye.
- Is highly absorbant, holding many times its weight in liquid
- Will not slip out of positions
- Contains no fibers there will not shed fibrous debris into the surgical site.
- Has a central opening that permits irrigation of the cornea.

VS18-450

Eye Kit

This two piece kit contains 1 Wipe/Wick VS18-470 and a Corneal Light Shield VS18-420 procedures to provide rapid recovery of fluids

- The Wipe/Wick consists of an instrument Wipe and has a 3" removable Eye Wick to clear fluids (see page 7 for full description of the Wipe/Wick)
- The Corneal Light Shield blocks excessive light that can damage the retina (see page 4 for full description of the Corneal Light Shield)
- Combination Kit saves time, money and offers easier inventory control



Product Number	Description	Quantity
VS18-430	Sterile EyeDrain, 80cc	1/pack, 20 pkg/box
VS18-435	Sterile Eye Drain, 400cc	1/pack, 10 pkg/box
VS18-431	Sterile Eye Wicks, w/o bag, 8 inches (20 cm)	2/pack, 10 pkg/box
VS18-440	Sterile Bloomberg Ophthalmic Ring	1/pkg, 20 pkg/box
VS18-450	Sterile Precision Eye Kit (2mm 3.25 x 3.25)	1/pack, 20 pkg/box

Wipes and Spears

Ophthalmology




VS18-451

All-in-One Cataract Kit

This 4-product kit contains 5 Cellulose Eye Spears VS18-410, 1 Corneal Light Shield VS18-420, 1 Eye Drain VS18-430 and 1 Instrument Wipe VS18-900

- Provides an easy-to-use one unit system when performing cataract surgery
- Saves time ordering and setting up for procedures
- Offers easier inventory control
- More economical than ordering individual items



VS18-460

PVA Diamond Knife Cleaning Block

Removes debris from the diamond which prevents the blade from dragging in the tissue

- 4.5cm long block made of hydroxylated polyvinyl acetate
- Continuous cell structure contains no fibers and, therefore, will not shed fibrous debris
- Small pore structure especially suited for absorbing bodily fluids, while applying ample pressure to wipe the diamond clean



Product Number	Description	Quantity
VS18-451	Sterile Cataract Kit	1/pack, 10 pkg/box
VS18-460	Knife Cleaning Block	1/pack, 10pkg/box


Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020

VS18-900 / VS18-910

Instrument Wipe

Ideal for procedures where gentle cleaning of a delicate surgical instrument is important

- 100% fiber-free material will not shed fibrous debris onto the instrument
- Soft chamois-like feel will not damage sensitive instruments or optics
- Will not snag on the box locks or serrations of an instrument's jaws when moistened
- Especially gentle on the sharps, fine points, and tips of delicate micro instruments
- Available in two different thicknesses – 2mm and 1mm




VS18-470

Wipe / Wick

Separates into two products; a standard 3"x 3" Instrument Wipe and a 3" long Eye Wick. The Wick is placed in the lower canthal fold to continually drain excess fluid from the surgical site.

- Contains no fibers
- Saves time and money
- Makes inventory control easier



VS18-990

Lasik Wipe

Product No.	Description	Quantity
VS18-900	Sterile Instrument Wipe (2mm thick) (3 x 3 x 2mm)	1 pack, 20 pkg/box
VS18-910	Sterile Instrument Wipe (1mm thick) (3 x 3 x 1mm)	1 pack, 20 pkg/box
VS18-470	Sterile Wipe/Wick	1 pack, 20 pkg/box
VS18-990	Lasik Instrument Wipe (2mm thick) (2 x 3 2mm)	1 pack, 20 pkg/box



Terms of Sale

- This price list supersedes all price lists.
- All prices are subject to change without notice
- Our terms are Net 30 days to approved accounts
- List prices are subject to change
- All Vision Surgical instruments are guaranteed to be free of defects in materials and workmanship. Any instrument which proves defective in workmanship or material will be repaired or replaced at our discretion with out charge. This guarantee is void for instruments serviced by any company other than Vision Surgical instrument repair service.
- All Vision Surgical instruments are carefully inspected before shipment. Any claims must be made within (2) two days of receipt of shipment. Most orders are shipped same day and, at the latest, by UPS within 72 hours.

Return Policy

- With our prior written authorization, returns may be made within 30 days of purchase for full credit. Instruments must be unused and in their original package.
- To insure proper credit, you must request our "Return Authorization Form" and return the completed to us. A restocking and inspection charge of 20% will be charged on all authorized returns more than (30) days after invoice date. No returns will be accepted after 90 days. A restocking charge may be assess on any return that in our determination was not an error of Vision Surgical.
- A restocking and inspection charge of 35% will be charged for on all authorized returns more than (60) sixty days after our invoice date.
- Goods accepted for return will be on a merchandise credit basis only and we can only accept prepaid shipments.
- Any merchandise returned without approval will be subject to an automatic restocking fee with a minimum of 35%.

Procedure

1. Notify our Customer Service Department in writing or by phone of your request for return authorization.
2. We will mail a "Return Authorization Request" form. Please complete the form with the following information:
 - (A) Our stock number
 - (B) Invoice Number
 - (C) Date of invoice
 - (D) Reason for return
3. When your form has been approved, we will issue you an authorization to return and a label.
4. The following types of goods are not eligible for credit:
 - (A) Custom designed instruments
 - (B) Instruments which have been used, marked or altered in anyway
 - (C) Engraved instruments
 - (D) Discontinued products
 - (E) Sterile products
 - (F) Color coated instruments

Prices subject to change without notice. Vision Surgical is not accountable for any errors or omissions that may appear in the catalog.

(Please fill out completely)

Your Address

**B
I
L
L
T
O**

Attention: Accounts Payable _____

Check if Same As Bill To Address

**S
H
I
P
T
O**

Contact Name _____

Phone (____) _____

Purchase _____

Order Number _____

Authorization _____

Signature _____

Qty.	Product Number	Page No.	Description of Product	Internal Use Only

*** Please call for address**

Repair Form

Ophthalmology



(Please fill out completely)

B I L L T O	Your Address

	Attention: Accounts Payable

S H I P T O	<input type="checkbox"/> Check if Same As Bill To Address

Contact Name _____
Phone (_____) _____
Purchase _____
Order Number _____

Repair Authorization _____
Signature _____
 Please call before repairs are started.

Qty.	Product Number	Description of Problem	Internal Use Only

*** Please call for address**

Toll-Free. 1-888-518-4746 | Fax. 1-214-722-2020

Vision Surgical

Toll-Free. 1-888-518-4746
Fax. 1-214-722-2020

Website. www.Vision-Surgical.com
Email. CustSer@Vision-Surgical.com